

## ALKALINE PHOSPHATASE

2 x 50 / 2 x 15 mL  
12006002

### Intended Use

This reagent is intended for *in vitro* quantitative determination of Alkaline Phosphatase in serum or plasma.

-DGKC – SCE recommended procedure.

-Linear up to 700 U/L

### Clinical Significance


Alkaline phosphatase (ALP) is widely distributed throughout the body, but clinically important one for diagnostic reasons are in bone, liver, placenta & intestine. Growing bone is associated with the release of ALP and so in childhood the level of ALP is around 3 times of that of adult. During pregnancy in 2<sup>nd</sup> & 3<sup>rd</sup> trimester the enzyme rises considerably due to placenta releasing ALP. It can be used to examine placental function.

Elevated levels are seen in bone diseases, e.g. Paget's disease, Rickets, Osteoblastic metastatic & in obstructive disease of biliary tract.

Decreased levels are rarely seen. e.g. in Vitamin A resistant rickets.

### Principle

Kinetic determination of Alkaline Phosphatase (ALP) is based upon the following reaction.



### Kit Components

| Reagent/Component       | Product Code 12006002 | Description                                                                   |
|-------------------------|-----------------------|-------------------------------------------------------------------------------|
| Alkaline Phosphatase R1 | 2 x 50 mL             | Diethanolamine Buffer (pH 10.2) 125 mmol/L<br>Magnesium Chloride 0.625 mmol/L |
| Alkaline Phosphatase R2 | 2 x 15 mL             | P-Nitrophenyl phosphate 50 mmol /L                                            |

### Risk & Safety

Material Safety data sheets (MSDS) will be provided on request

### Reagent Preparation

Alkaline Phosphatase reagents are ready to use.

### Reagent Storage and Stability

The sealed reagents are stable upto the expiry date stated on the label, when stored at 2-8°C.

### Open Vial Stability

Once opened the reagents are stable up to 30 days if contamination is avoided.

### On-board Calibration Stability

Calibration is stable for 7 days.

### Reagent Deterioration

Turbidity or precipitation in any kit component indicates deterioration and the component must be discarded. Values outside the recommended acceptable range for the Agappe Qualicheck Norm & Path control may also be an indication of reagent instability and associated results are invalid. Sample should be retested using fresh vial of reagent.

### Precaution

To avoid contamination, use clean laboratory wares. Close reagent bottles immediately after use. Avoid direct exposure of reagent to light. Do not blow into the reagent bottles.

This reagent is only for IVD use and follow the normal precautions required for handling all laboratory reagents.

### Waste Management

Reagents must be disposed off in accordance with local regulations.

### Sample

Fresh serum / plasma (Do not use lipemic or hemolysed sample)

### Interferences

No interference for

Bilirubin up to 10 mg/dL

Haemoglobin up to 1000 mg/dL

Ascorbic Acid up to 50 mg/dL

### Materials provided

Alkaline Phosphatase reagents

### Reagents required but not provided

Multicalibrator (Product Code: 11610001), Qualicheck Norm (Product Code: 11601003), Qualicheck Path (Product Code: 11601002)

### Calibration

Agappe Multicalibrator (Product Code: 11610001) is recommended for calibration of the assay.

### Quality control

It is recommended to use Qualicheck Norm (Product Code: 11601003) or Qualicheck Path (Product Code: 11601002) to verify the performance of the measurement procedure.

Each Laboratory has to establish its own internal quality control scheme and procedures for corrective action if controls do not recover within the acceptable tolerance.

### Reference Range

It is recommended that each laboratory should establish its own reference values. The following value may be used as guide line.

Women : 64 – 306 U/L

Men : 80 – 306 U/L

Children : 180 - 1200 U/L

Results obtained for patient samples are to be correlated with clinical findings of patient for interpretation and diagnosis.

### Performance

#### 1. Linearity

The reagent is linear up to 700 U/L. If the concentration is greater than linearity (700 U/L), dilute the sample with normal saline and repeat the assay. Multiply the result with dilution factor.

#### 2. Comparison

A comparison study has been performed between Agappe reagent and another internationally available reagent yielded a correlation coefficient of  $r^2 = 0.9529$  and a regression equation of  $y = 0.9697x$ .

#### 3. Precision

|            | Intra Run |         | Inter Run |         |
|------------|-----------|---------|-----------|---------|
| Control    | Level 1   | Level 2 | Level 1   | Level 2 |
| n          | 20        | 20      | 20        | 20      |
| Mean (U/L) | 174.1     | 665.6   | 170.6     | 660.0   |
| SD         | 6.02      | 20.31   | 6.45      | 25.32   |
| CV(%)      | 3.45      | 3.05    | 3.7       | 3.8     |

| Accuracy (U/L)  |                |                |
|-----------------|----------------|----------------|
| Control         | Expected Value | Measured Value |
| Control Level 1 | 160 ± 55       | 162.3          |
| Control Level 2 | 510 ± 95       | 511.2          |
| Qualicheck Norm | 170 ± 18       | 175.3          |
| Qualicheck Path | 400 ± 51       | 396.2          |

#### 4. Sensitivity

Lower detection Limit is 1.5 U/L

### Bibliography

- Schlebusch, H. *et al.*, Dtsh.med.Wschr.99,765(1974)
- Z.Klin. Chem., Klin. Biochem.8,658(1980)10,182(1972)

#### SYMBOLS USED ON THE LABELS

IVD IN VITRO DIAGNOSTIC USE  SEE PACKAGE INSERT FOR PROCEDURE  LOT LOT NUMBER  MANUFACTURER'S ADDRESS  MANUFACTURING DATE  EXPIRY DATE  TEMPERATURE LIMIT



**AGAPPE DIAGNOSTICS LTD.**

REV. NO.: ADL/IFU/ALP/40FR/R01

 ISO 9001:2015  
EN ISO 13485:2016

Accute 40 FR Assay Parameter

|                                                                |                      |
|----------------------------------------------------------------|----------------------|
| Page 1                                                         |                      |
| Test                                                           | ALKALINE PHOSPHATASE |
| Reaction Mode                                                  | Rate                 |
| Reference Test ID                                              | **                   |
| Test WL                                                        | 404                  |
| Blank WL                                                       | NA                   |
| Test Read Timing                                               | 40-52                |
| Blank Read Timing                                              | NA                   |
| Sample                                                         | 5 µL                 |
| R1                                                             | 200 µL               |
| R2                                                             | 50 µL                |
| Stirrer                                                        | ON                   |
| Cal Mode                                                       | std                  |
| Page 3                                                         |                      |
| Calibration Mode                                               | Linear               |
| ** Not applicable                                              |                      |
| # - Input Programme Number only for sample Blanking parameters |                      |

SYMBOLS USED ON THE LABELS

 IN VITRO DIAGNOSTIC USE
  SEE PACKAGE INSERT FOR PROCEDURE
  LOT NUMBER
  MANUFACTURER'S ADDRESS
  MANUFACTURING DATE
  EXPIRY DATE
  TEMPERATURE LIMIT



AGAPPE DIAGNOSTICS LTD.

REV. NO.: ADL/IFU/ALP/40FR/R01

 ISO 9001:2015  
EN ISO 13485:2016