

INTRODUCTION

1. AutoZyme Creatinine is a reagent set for determination of creatinine based on **initial rate method** using Alkaline Picrate.
2. AutoZyme Creatinine is a **single reagent system**, using one step procedure.
3. AutoZyme Creatinine has **one step reconstitution**. It involves mixing of Picrate and Diluent reagent.
4. AutoZyme Creatinine is a **High Stability Reagent**.
5. AutoZyme Creatinine is **linear** upto 30 mg%.
6. Creatinine can be determined in **180 seconds**.
7. AutoZyme Creatinine can be used on any **Spectrophotometer, Discrete semiautomated and Automated analyzer**. Programme can be designed for any specific analyzer upon request.

PRINCIPLE

Creatinine in alkaline medium reacts with picrate to produce orange colour. This colour absorbs light at 492 nm.(490 - 510 nm.). The rate of increase in absorbance is directly proportional to the concentration of creatinine in specimen.

Creatinine + Picrate $\xrightarrow{\text{Alkaline medium}}$ Orange colour

PREPARATION OF WORKING SOLUTION

Prepare working solution by mixing **equal volume of Picrate Reagent and Diluent Reagent**.

REAGENT STORAGE & STABILITY

The reagents are stable till the expiry date stated on the bottle label, when stored at R.T. (25-30°C).

The working solution is stable for 30 days at 2-8°C.

COMPONENTS & CONCENTRATION OF WORKING SOLUTION

Component	Concentration
• Sodium Picrate	7.7 mmol/l
• Sodium Hydroxide	500 mmol/l

SPECIMEN COLLECTION & PRESERVATION

Blood should be collected in a clean dry container. Avoid use of plastic or siliconized container which may prolong clotting time. Samples should not be collected during PSP/BSP clearance test. For plasma separation Heparin (200 IU/ml blood) may be used as anticoagulant.

Creatinine in serum and plasma is stable for 2 days when stored at 2-8°C.

PROCEDURE

- ☐ Reaction type Initial rate
- ☐ Reaction direction Up
- ☐ Wavelength 492 nm.(490 - 510 nm.)
- ☐ Flowcell temperature 30°C / 37°C
- ☐ Zero setting with Distilled water
- ☐ Delay time 30 seconds
- ☐ No. of readings 2
- ☐ Interval 60 seconds
- ☐ Sample volume 0.05 ml (50 µl)
- ☐ Reagent volume 1.0 ml
- ☐ Standard concentration 2 mg %
- ☐ Factor $2 \div \Delta \text{ Abs. of standard}$
- ☐ Linearity 30 mg/dl

Manual assay procedure

Prewarm the required amount of working solution to 30°C/37°C before use.

1.0 ml procedure

Standard / Sample	0.05 ml (50 µl)
Working Solution	1.0 ml

Mix and start stopwatch simultaneously. Record absorbance assay mixture at exactly 30 seconds after Standard / Specimen addition and then again at 90 seconds.

Note : It is recommended to run a creatinine standard with each batch of assay.

Calculation:

Calculate the average change in absorbance per minute ($\Delta \text{ Abs.}$) of standard & specimen(s).

$\Delta \text{ Abs.} = \text{Abs. at 90 sec.} - \text{Abs. at 30 sec.}$

$$\text{Serum Creatinine (mg\%)} = \frac{\Delta \text{ Abs. of Specimen}}{\Delta \text{ Abs. of Standard}} \times 2$$

Initial Rate