DCM113-7

Ed. 01/2024

Anti Cardiolipin IgG

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18. BIBLIOGRAFIA

1. Dusse LM, Silva FD, Freitas LG, Rios DR, *et al.* Antiphospholipid syndrome: a clinical and laboratorial challenge. Rev Assoc Med Bras (1992). 2014 Mar-Apr;60(2):181-6.
2. Roggenbuck D, Somma V, Schierack P, Borghi MO, Meroni PL. Autoantibody profiling in APS. Lupus. 2014 Oct;23(12):1262-4.
3. Banzato A, Pengo V. Clinical relevance of β_2 -glycoprotein-I plasma levels in antiphospholipid syndrome (APS). Curr Rheumatol Rep. 2014 Jun;16(6):424.
4. Amiral J, Peyrafitte M, Dunois C, Vissac AM, Seghatchian J. Anti-phospholipid syndrome: Current opinion on mechanisms involved, laboratory characterization and diagnostic aspects. Transfus Apher Sci. 2017 Aug;56(4):612-625.
5. Miyakis S, Lockshin MD, Atsumi T, Branch, *et al.* international consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). J Thromb Haemost. 2006 Feb;4(2):295-306.
6. Carreras LO, Forastiero RR, Martinuzzo ME. Which are the best biological markers of the antiphospholipid syndrome? J Autoimmun. 2000 Sep;15(2):163-72.
7. Cervera R, Piette JC, Font J, Khamashta MA, Shoenfeld Y, *et al.*; Euro-Phospholipid Project Group. Antiphospholipid syndrome: clinical and immunologic manifestations and patterns of disease expression in a cohort of 1,000 patients. Arthritis Rheum. 2002 Apr;46(4):1019-27.
8. Linnemann B. Antiphospholipid syndrome - an update. Vasa. 2018 Oct;47(6):451-464.
9. Calvo-Alén J, Toloza SM, Fernández M, Bastian HM, *et al.*; LUMINA Study Group. Systemic lupus erythematosus in a multiethnic US cohort (LUMINA). XXV. Smoking, older age, disease activity, lupus anticoagulant, and glucocorticoid dose as risk factors for the occurrence of venous thrombosis in lupus patients. Arthritis Rheum. 2005 Jul;52(7):2060-8.
10. Janardhan V, Wolf PA, Kase CS, Massaro JM, *et al.* Anticardiolipin antibodies and risk of ischemic stroke and transient ischemic attack: the Framingham cohort and offspring study. Stroke. 2004 Mar;35(3):736-41.

11. Harris EN, Pierangeli SS. Revisiting the anticardiolipin test and its standardization. *Lupus*. 2002;11(5):269-75.
12. McClain MT, Arbuckle MR, Heinlen LD, Dennis GJ, *et al.* The prevalence, onset, and clinical significance of antiphospholipid antibodies prior to diagnosis of systemic lupus erythematosus. *Arthritis Rheum*. 2004 Apr;50(4):1226-32.
13. Singh AK. Immunopathogenesis of the antiphospholipid antibody syndrome: an update. *Curr Opin Nephrol Hypertens*. 2001 May;10(3):355-8.
14. de Laat B, de Groot PG. Autoantibodies directed against domain I of beta2-glycoprotein I. *Curr Rheumatol Rep*. 2011 Feb;13(1):70-6.
15. Suh-Lailam BB, Cromar A, Davis KW, Tebo AE. APbL antibody ELISA as an alternative to anticardiolipin test for the diagnosis of antiphospholipid syndrome. *Int J Clin Exp Pathol*. 2012;5(3):210-5.
16. Ofer-Shiber S, Molad Y. Frequency of vascular and pregnancy morbidity in patients with low vs. moderate-to-high titers of antiphospholipid antibodies. *Blood Coagul Fibrinolysis*. 2015 Apr;26(3):261-6.
17. Meroni PL, Borghi MO, Raschi E, Tedesco F. Pathogenesis of antiphospholipid syndrome: understanding the antibodies. *Nat Rev Rheumatol*. 2011 Jun;7(6):330-9.
18. de Laat B, Pengo V, Pabinger I, Musial J, *et al.* The association between circulating antibodies against domain I of beta2-glycoprotein I and thrombosis: an international multicenter study. *J Thromb Haemost*. 2009 Nov;7(11):1767-73.
19. Sanmarco M, Roll P, Gayet S, Oksman F, *et al.* J. Combined search for anti-beta2-glycoprotein I and anticardiolipin antibodies in antiphospholipid syndrome: contribution to diagnosis. *J Lab Clin Med*. 2004 Sep;144(3):141-7.
20. Otomo K, Atsumi T, Amengual O, Fujieda Y, *et al.* Efficacy of the antiphospholipid score for the diagnosis of antiphospholipid syndrome and its predictive value for thrombotic events. *Arthritis Rheum*. 2012 Feb;64(2):504-12.
21. Basic QC Practices On-line Course; <http://www.Westgard.com>.
22. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.

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Anti Cardiolipin IgG

for routine analysis

Quantitative determination of IgG autoantibodies against Cardiolipin in human serum or plasma

IVD

LOT

See external label

2°C  8°C $\Sigma = 96$ tests

REF DKO113

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use
For Laboratory Professional Use

Anti Cardiolipin IgG is a manual *in vitro* diagnostic device intended for the quantitative determination of IgG class antibodies directed against cardiolipin in human serum or plasma from an adult population. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of antiphospholipid syndrome (APS).

2. CLINICAL SIGNIFICANCE

Cardiolipin is a negatively charged phospholipid which is typically located in the inner mitochondrial membrane¹. Autoantibodies directed against cardiolipin are part of a group known as antiphospholipid antibodies which includes anti- β 2 glycoprotein 1 autoantibodies. Measurement of anti-cardiolipin autoantibodies is considered to be one of the most important markers to support diagnosis of antiphospholipid syndrome (APS)^{2,3}.

APS is a systemic autoimmune disorder characterised by a combination of arterial and/or venous thromboses, pregnancy complications, such as recurrent foetal loss, together with elevated levels of antiphospholipid antibodies⁴. APS was first described in patients with systemic lupus erythematosus (SLE), though it has subsequently been established that SLE may be independent of an underlying disease⁵.

APS can occur alone – Primary APS, or in association with other conditions, such as SLE – Secondary APS⁶ (6). However, it has been demonstrated that anti-cardiolipin antibodies (aCLs) can be detected in patients with SLE that do not develop Secondary APS⁷⁻⁹. However, thromboembolic events are the most common clinical manifestation of APS

Anti-cardiolipin antibodies can recognise both cardiolipin and portions of the phospholipid-protein complex β 2 glycoprotein 1-cardiolipin^{10,11}.

Studies have indicated an association of anti-cardiolipin IgG and IgM antibodies^{2,7,11-16}, with thrombotic events whereas

others suggest these to be linked with IgG isotype, but not IgM^{6,10,17}.

aCL IgM antibodies have been shown to occur in infections such as chronic hepatitis C, leprosy, syphilis, but they are not directly involved in thrombotic events¹⁷.

It is likely that the presence of antiphospholipid antibodies, including aCL IgG and IgM constitutes the single most recognisable risk factor in cases of recurrent pregnancy loss and late placenta-mediated obstetric complications^{4-6,11,14-15,18,19}. Where patients can present with only adverse pregnancy outcomes with isolated vascular events or with both obstetric and thrombotic manifestations⁶. It has been suggested that anti- β 2-glycoprotein I – cardiolipin antibodies are able to recognise the antigen on placental tissues, inhibiting the growth and differentiation of trophoblasts which may eventually cause defective placentation²⁰.

3. PRINCIPLE OF THE METHOD

The Anti Cardiolipin IgG assay is a two-step sandwich enzyme immunometric assay (ELISA) where patient samples, calibrators or controls are incubated on microtitre plates coated with the antigenic cardiolipin- β 2 glycoprotein complex. During the incubation, antibodies present in the test sample bind to the immobilised antigen complex. After the incubation, the bound/free separation is performed by a simple solid phase washing.

A subsequent incubation occurs with anti-human IgG conjugated with horseradish peroxidase (HRP), which binds to the immobilised antibodies. A further wash step is performed to remove excess conjugate. Then, a chromogenic substrate solution containing TMB is dispensed into the wells which reacts with the conjugated HRP and a blue colour develops that changes into yellow when the Stop Solution (H₂SO₄) is added. The colour intensity is directly proportional to the Anti Cardiolipin IgG concentration in the sample.

Anti-cardiolipin antibody concentration in the sample is calculated through a calibration curve.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (5 vials, 1.2 mL each)

Phosphate buffer 0.1M, NaN_3 <0.1%, BSA 3%

CAL0 REF DCE002/11306-0

CAL1 REF DCE002/11307-0

CAL2 REF DCE002/11308-0

CAL3 REF DCE002/11309-0

CAL4 REF DCE002/11310-0

2. Controls (2 vials, 1.2 mL each, ready for use)

Phosphate buffer 0.1 M, NaN_3 <0.1%, BSA 3%

Negative Control REF DCE045/11301-0

Positive Control REF DCE045/11302-0

3. Sample Diluent (1 vial, 100 mL)

Phosphate buffer 0.1M, NaN_3 <0.1%

REF DCE053-0

4. Conjugate (1 vial, 15 mL)

Anti h-IgG conjugate with horseradish peroxidase (HRP),

BSA 0.1%, ProClin >0.0015% REF DCE002/11302-0

5. Coated Microplate (1 breakable microplate)

Microplate coated with antigenic Cardiolipin- β 2-

Glycoprotein complex REF DCE002/11303-0

6. TMB Substrate (1 vial, 15 mL)

H_2O_2 -TMB (0.26 g/L) (avoid any skin contact), ProClin

<0.0015% REF DCE004-0

7. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

8. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M pH 7.4 ProClin >0.0015%

REF DCE054-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.

- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- Sodium Azide may be toxic if ingested or absorbed Some reagents (calibrators, controls and sample diluent) contain small amounts of Sodium Azide (NaN_3) <0.1% which may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through with large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ H_2O_2 to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.

- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinisation and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analysers which are not equipped with disposable tips.**
For this purpose, DiaMetra supplies a separate decontamination reagent for cleaning needles.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin, potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	96 hours
Freeze/thaw cycles	3 cycles

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

The calibrators and controls are ready to use.

The calibrators have approximately the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄
AU/mL	0	5	10	20	80

Calibrator / control concentrations are printed on the label.

9.2. Preparation of the Wash Solution

Dilute the content of the vial “10X Conc. Wash Solution” with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.3. Preparation of Samples

All serum and plasma samples must be diluted 1:100 with sample diluent.

e.g., 10 µL of sample should be diluted with 990 µL of sample diluent.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2-8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₄), two for each control, two for each sample, one for Blank.

Reagents	Calibrator	Sample/ Controls	Blank
Calibrator C ₀ -C ₄	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	

Incubate for 60 minutes at room temperature (22-28°C). Remove the content from each well; wash the wells 3 times with 300 µL diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

Conjugate	100 µL	100 µL	
-----------	--------	--------	--

Incubate for 60 minutes at room temperature (22-28°C). Remove the content from each well; wash the wells 3 times with 300 µL diluted wash solution.

Washing: follow the same indications of the previous point.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate for 15 minutes in the dark at room temperature (22-28°C).

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake the microplate gently.
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate. This test is only valid if the optical density at 450 nm for the Positive Control as well as for the Calibrators (C₀-C₄) fall within with the respective range indicated on the Quality Control Certificate enclosed in each test kit.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give

control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories²¹.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit with lin-log coordinates, **including Calibrator 0 is required**. A smoothed spline fit including Calibrator 0 can be used. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

12. MEASURING RANGE

The assay measuring range (AMR) is 2 – 80 AU/mL. Any value that reads below 2 AU/mL should be reported as “< 2 AU/mL”. Any value that reads above 80 AU/mL should be reported as “> 80 AU/mL”.

13. METROLOGY AND TRACEABILITY

The calibrators of this kit are traceable to the Centres for Disease Control (CDC) Human IgG Anti-Cardiolipin Monoclonal Antibody HCal - Catalogue [IS2717].

14. INTERPRETATION OF RESULTS

Index / Concentration	Interpretation
< 8 AU/mL	The sample should be considered negative
8 – 10 AU/mL	The sample should be graded equivocal and repeat testing / sampling should be performed according to internal practices
>10 AU/mL	The sample should be considered positive

Determination of a range of expected values for a “normal” population of a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore, each

laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population.

Positive results should be verified concerning the entire clinical status of the patient, with the decision for therapy being taken on an individual basis. It is recommended that each laboratory establishes its own normal and pathological ranges of Anti-Cardiolipin IgG values.

15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.59 AU/mL
Limit of Detection (LoD)	1.25 AU/mL
Limit of Quantitation (LoQ)	2.00 AU/mL

15.2. Trueness

Trueness of the Anti Cardiolipin IgG assay has been demonstrated through performance of a recovery test using the CDC Human IgG Anti-Cardiolipin Monoclonal Antibody HCAL - Catalogue [IS2717].

15.3. Diagnostic Sensitivity and Specificity

The sensitivity and specificity were determined with guidance from CLSI EP-24 "Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves" using 50 negative and 51 positive samples run on two reagent lots.

		DKO113		Total
		Positive	Negative	
True state	Positive	47	4	51
	Negative	0	57	57
Total		47	61	108

Diagnostic sensitivity: 92%

Diagnostic specificity: 100%

15.4. Precision

Precision of the Anti Cardiolipin IgG was determined by performing a complex precision study.

Repeatability: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (AU/mL)	Within run (Repeatability)	
			SD	CV
1	75	6.95	0.38	5.5%
2	75	11.03	0.51	4.6%
3	75	20.12	0.94	4.7%
4	75	30.26	1.86	6.1%
5	75	50.11	2.58	5.1%
6	75	71.71	2.53	3.5%

Reproducibility: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (AU/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	6.98	0.49	7.0%
2	150	11.17	0.84	7.5%
3	150	20.16	1.87	9.3%
4	150	30.38	3.06	10.1%
5	150	50.76	5.02	9.9%
6	150	72.30	3.93	5.4%

15.5. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For anti-cardiolipin IgG concentration by Anti Cardiolipin IgG, the measurement procedure shows linearity for the interval from 0.84 to 83.68 ng/mL within the allowable deviation of linearity (ADL) of ±15 %.

15.6. Analytical Specificity

The following substances do not interfere with a bias of > ±15% in the Anti Cardiolipin IgG assay when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	10 g/dL
Triglyceride	500 mg/dL

15.7. Serum-plasma study

The Anti Cardiolipin IgG matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 22 samples (18 native, 4 spiked) to cover the range were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	0.96 [0.92 – 0.98]	0.64 [-0.38 to 1.66]	1.00
Lithium Heparin	0.92 [0.88 to 0.96]	0.87 [-0.24 to 1.99]	1.00
Sodium Heparin	0.94 [0.89 to 0.98]	0.66 [-0.75 to 2.06]	0.99
EDTA	0.95 [0.92 to 0.99]	0.54 [-0.69 to 1.77]	1.00

16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays²². Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.
- The presence of immune complexes or other immunoglobulin aggregates in the patient sample may cause an increased level of non-specific binding and produce false positives in this assay.

17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations

18. BIBLIOGRAPHY

1. Dusse LM, Silva FD, Freitas LG, Rios DR, *et al.* Antiphospholipid syndrome: a clinical and laboratorial challenge. *Rev Assoc Med Bras* (1992). 2014 Mar-Apr;60(2):181-6.
2. Roggenbuck D, Somma V, Schierack P, Borghi MO, Meroni PL. Autoantibody profiling in APS. *Lupus*. 2014 Oct;23(12):1262-4.
3. Banzato A, Pengo V. Clinical relevance of β_2 -glycoprotein-I plasma levels in antiphospholipid syndrome (APS). *Curr Rheumatol Rep*. 2014 Jun;16(6):424.
4. Amiral J, Peyrafitte M, Dunois C, Vissac AM, Seghatchian J. Anti-phospholipid syndrome: Current opinion on mechanisms involved, laboratory characterization and diagnostic aspects. *Transfus Apher Sci*. 2017 Aug;56(4):612-625.
5. Miyakis S, Lockshin MD, Atsumi T, Branch, *et al.* international consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost*. 2006 Feb;4(2):295-306.
6. Carreras LO, Forastiero RR, Martinuzzo ME. Which are the best biological markers of the antiphospholipid syndrome? *J Autoimmun*. 2000 Sep;15(2):163-72.
7. Cervera R, Piette JC, Font J, Khamashta MA, Shoenfeld Y, *et al.* Euro-Phospholipid Project Group. Antiphospholipid syndrome: clinical and immunologic manifestations and patterns of disease expression in a cohort of 1,000 patients. *Arthritis Rheum*. 2002 Apr;46(4):1019-27.
8. Linnemann B. Antiphospholipid syndrome - an update. *Vasa*. 2018 Oct;47(6):451-464.
9. Calvo-Alén J, Toloza SM, Fernández M, Bastian HM, *et al.*; LUMINA Study Group. Systemic lupus erythematosus in a multiethnic US cohort (LUMINA). XXV. Smoking, older age, disease activity, lupus anticoagulant, and glucocorticoid dose as risk factors for the occurrence of venous thrombosis in lupus patients. *Arthritis Rheum*. 2005 Jul;52(7):2060-8.
10. Janardhan V, Wolf PA, Kase CS, Massaro JM, *et al.* Anticardiolipin antibodies and risk of ischemic stroke and transient ischemic attack: the Framingham cohort and offspring study. *Stroke*. 2004 Mar;35(3):736-41.
11. Harris EN, Pierangeli SS. Revisiting the anticardiolipin test and its standardization. *Lupus*. 2002;11(5):269-75.
12. McClain MT, Arbuckle MR, Heinlen LD, Dennis GJ, *et al.* The prevalence, onset, and clinical significance of antiphospholipid antibodies prior to diagnosis of systemic lupus erythematosus. *Arthritis Rheum*. 2004 Apr;50(4):1226-32.
13. Singh AK. Immunopathogenesis of the antiphospholipid antibody syndrome: an update. *Curr Opin Nephrol Hypertens*. 2001 May;10(3):355-8.
14. de Laat B, de Groot PG. Autoantibodies directed against domain I of beta2-glycoprotein I. *Curr Rheumatol Rep*. 2011 Feb;13(1):70-6.
15. Suh-Lailam BB, Cromar A, Davis KW, Tebo AE. APHL antibody ELISA as an alternative to anticardiolipin test for the diagnosis of antiphospholipid syndrome. *Int J Clin Exp Pathol*. 2012;5(3):210-5.
16. Ofer-Shiber S, Molad Y. Frequency of vascular and pregnancy morbidity in patients with low vs. moderate-to-high titers of antiphospholipid antibodies. *Blood Coagul Fibrinolysis*. 2015 Apr;26(3):261-6
17. Meroni PL, Borghi MO, Raschi E, Tedesco F. Pathogenesis of antiphospholipid syndrome: understanding the antibodies. *Nat Rev Rheumatol*. 2011 Jun;7(6):330-9.

18. de Laat B, Pengo V, Pabinger I, Musial J, *et al.* The association between circulating antibodies against domain I of beta2-glycoprotein I and thrombosis: an international multicenter study. *J Thromb Haemost.* 2009 Nov;7(11):1767-73.
19. Sanmarco M, Roll P, Gayet S, Oksman F, *et al.* J. Combined search for anti-beta2-glycoprotein I and anticardiolipin antibodies in antiphospholipid syndrome: contribution to diagnosis. *J Lab Clin Med.* 2004 Sep;144(3):141-7.
20. Otomo K, Atsumi T, Amengual O, Fujieda Y, *et al.* Efficacy of the antiphospholipid score for the diagnosis of antiphospholipid syndrome and its predictive value for thrombotic events. *Arthritis Rheum.* 2012 Feb;64(2):504-12.
21. Basic QC Practices On-line Course; <http://www.Westgard.com>.
22. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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5. Miyakis S, Lockshin MD, Atsumi T, Branch, *et al.* international consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost.* 2006 Feb;4(2):295-306.
6. Carreras LO, Forastiero RR, Martinuzzo ME. Which are the best biological markers of the antiphospholipid syndrome? *J Autoimmun.* 2000 Sep;15(2):163-72.
7. Cervera R, Piette JC, Font J, Khamashta MA, Shoenfeld Y, *et al.*; Euro-Phospholipid Project Group. Antiphospholipid syndrome: clinical and immunologic manifestations and patterns of disease expression in a cohort of 1,000 patients. *Arthritis Rheum.* 2002 Apr;46(4):1019-27.
8. Linnemann B. Antiphospholipid syndrome - an update. *Vasa.* 2018 Oct;47(6):451-464.
9. Calvo-Alén J, Toloza SM, Fernández M, Bastian HM, *et al.*; LUMINA Study Group. Systemic lupus erythematosus in a multiethnic US cohort (LUMINA). XXV. Smoking, older age, disease activity, lupus anticoagulant, and glucocorticoid dose as risk factors for the occurrence of venous thrombosis in lupus patients. *Arthritis Rheum.* 2005 Jul;52(7):2060-8.
10. Janardhan V, Wolf PA, Kase CS, Massaro JM, *et al.* Anticardiolipin antibodies and risk of ischemic stroke and transient ischemic attack: the Framingham cohort and offspring study. *Stroke.* 2004 Mar;35(3):736-41.
11. Harris EN, Pierangeli SS. Revisiting the anticardiolipin test and its standardization. *Lupus.* 2002;11(5):269-75.
12. McClain MT, Arbuckle MR, Heinlen LD, Dennis GJ, *et al.* The prevalence, onset, and clinical significance of antiphospholipid antibodies prior to diagnosis of systemic lupus erythematosus. *Arthritis Rheum.* 2004 Apr;50(4):1226-32.
13. Singh AK. Immunopathogenesis of the antiphospholipid antibody syndrome: an update. *Curr Opin Nephrol Hypertens.* 2001 May;10(3):355-8.
14. de Laat B, de Groot PG. Autoantibodies directed against domain I of beta2-glycoprotein I. *Curr Rheumatol Rep.* 2011 Feb;13(1):70-6.
15. Suh-Lailam BB, Cromar A, Davis KW, Tebo AE. APhL antibody ELISA as an alternative to anticardiolipin test for the diagnosis of antiphospholipid syndrome. *Int J Clin Exp Pathol.* 2012;5(3):210-5.
16. Ofer-Shiber S, Molad Y. Frequency of vascular and pregnancy morbidity in patients with low vs. moderate-to-high titers of antiphospholipid antibodies. *Blood Coagul Fibrinolysis.* 2015 Apr;26(3):261-6.
17. Meroni PL, Borghi MO, Raschi E, Tedesco F. Pathogenesis of antiphospholipid syndrome: understanding the antibodies. *Nat Rev Rheumatol.* 2011 Jun;7(6):330-9.
18. de Laat B, Pengo V, Pabinger I, Musial J, *et al.* The association between circulating antibodies against domain I of beta2-glycoprotein I and thrombosis: an international multicenter study. *J Thromb Haemost.* 2009 Nov;7(11):1767-73.
19. Sanmarco M, Roll P, Gayet S, Oksman F, *et al.* J. Combined search for anti-beta2-glycoprotein I and anticardiolipin antibodies in antiphospholipid syndrome: contribution to diagnosis. *J Lab Clin Med.* 2004 Sep;144(3):141-7.
20. Otomo K, Atsumi T, Amengual O, Fujieda Y, *et al.* Efficacy of the antiphospholipid score for the diagnosis of antiphospholipid syndrome and its predictive value for thrombotic events. *Arthritis Rheum.* 2012 Feb;64(2):504-12.
21. Basic QC Practices On-line Course; <http://www.Westgard.com>.
22. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33.

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