DCM114-8
Ed. 01/2024

Anti Phospholipid Screen

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Anti Phospholipid Screen

for routine analysis

Quantitative determination of autoantibodies against phospholipids in human serum or plasma

IVD

LOT

See external label

2°C 8°C

$\Sigma = 96$ tests

REF DKO114

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

Anti Phospholipid Screen is a manual *in vitro* diagnostic device intended for the quantitative determination of IgG and IgM class autoantibodies against β 2-glycoprotein mediated anionic phospholipids (including Cardiolipin, Phosphatidyl serine, Phosphatidyl inositol, Phosphatidic acid, Phosphatidyl choline, Lyso-phosphatidyl choline, Phosphatidyl ethanolamine) in human serum or plasma from an adult population. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of increased risk of thrombosis in patients with Systemic Lupus Erythematosus (SLE) or similar disorders.

2. CLINICAL SIGNIFICANCE

The first study on the anti-phospholipid antibodies began in 1906 when Wasserman introduced a serological test for syphilis. In 1942 it was found the active component that is a phospholipid indicated by the name of Cardiolipin. In the 50's it was observed that a large number of people appeared to be positive for syphilis tests but did not show any evidence of disease. At the beginning the phenomenon was classified as a series of false positive syphilis test, then a more accurate analysis revealed, for this group of patients, a high prevalence of autoimmune disorders including Systemic Lupus Erythematosus (SLE) and Sjögrens syndrome.

The term lupus anticoagulant (LA), used for the first time in 1972, derives from experimental observations in which it was observed an increased risk of thrombosis, paradoxically, with the presence of some anticoagulant factors; the term LA is not totally correct, in fact the disease is present more frequently in patients without lupus and it is associated with thrombosis rather than to abnormal bleeding.

Some years later the role of a cofactor has been investigated, the β 2-glycoprotein I (apolipoprotein H) also said β 2GPI, and its interactions with anionic phospholipids in human serum / plasma. This cofactor is a β -globulin with a molecular weight of 50 kDa that has the concentration of 200 μ g / mL in plasma. The β 2GPI is involved in the regulation of blood coagulation, inhibiting the intrinsic way. β 2GPI *in vivo* is associated with negatively charged substances such as anionic phospholipids, heparin and

lipoproteins. The region that binds phospholipids is in its fifth domain.

The acronym "aPL" (anti-phospholipid antibodies) indicates improperly antibodies directed against negatively charged phospholipids like Cardiolipin (CL), Phosphatidyl serine (PS) Phosphatidyl inositol (PI) and phosphatidic acid (PA); more correctly the term anti-phospholipid antibodies indicate those antibodies directed against the complex.

Between β 2GPI and anionic phospholipids that can bind to the fifth domain of β 2GPI. Among these, Cardiolipin is the most commonly used phospholipid as an antigen for determining the aPL with ELISA method. Diagnostic laboratories measure the antibodies directed against the complex between β 2GPI and negatively charged phospholipids, as Phosphatidyl serine (PS) Phosphatidyl inositol (PI) and phosphatidic acid (PA). Some researchers suggest the use of PS instead of Cardiolipin in ELISA assays, for a more precise diagnosis. However, these antibodies against phospholipids are less commonly used, even if their use may increase the clinical sensitivity of patient samples with suspected Anti-phospholipid Syndrome (APS) but cannot replace the determination of autoantibodies anti-Cardiolipin.

3. PRINCIPLE OF THE METHOD

The Anti Phospholipid Screen assay is an indirect enzyme immunometric assay (ELISA) where patient samples, calibrators or controls are incubated on microtitre plates coated with antigenic complex between anionic phospholipids and β 2-Glycoprotein (solid phase). During the incubation, antibodies present in the test sample bind to the antigenic complex. After the incubation, the bound/free separation is performed by a simple solid phase washing.

A subsequent incubation occurs with either anti-human IgG or anti-human IgM conjugated with horseradish peroxidase (HRP), which binds to the immobilised antibodies. A further wash step is performed to remove excess conjugate. Then, a chromogenic substrate solution containing TMB is dispensed into the wells which reacts with the conjugated HRP and a blue colour develops that changes into yellow when the Stop Solution (H_2SO_4) is added. The level of colour is directly proportional to the concentration of IgG or IgM antibodies present in the original sample.

IgG or IgM concentration in the sample is calculated through a calibration curve.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (5 vial, 1.2 mL each)
Phosphate buffer 0.1M, NaN_3 <0.1%, BSA 3%
CAL0 REF DCE002/11406-0
CAL1 REF DCE002/11407-0
CAL2 REF DCE002/11408-0
CAL3 REF DCE002/11409-0
CAL4 REF DCE002/11410-0
2. Controls (2 vials, 1.2 mL each)
Phosphate buffer 0.1M, NaN_3 <0.1%, BSA 3%
Positive Control REF DCE045/11402-0
Negative Control REF DCE045/11401-0
3. Sample Diluent (1 vial, 100 mL)
Phosphate buffer 0.1M, NaN_3 <0.1%
REF DCE053-0
4. Conjugate IgG (1 vial, 15 mL)
Anti- human IgG conjugated with horseradish peroxidase (HRP). ProClin >0.0015% REF DCE002/11402-G
5. Conjugate IgM (1 vial, 15 mL)
Anti- human IgM conjugated with horseradish peroxidase (HRP). ProClin >0.0015% REF DCE002/11402-M
6. Coated Microplate (1 microplate breakable)
Antigenic phospholipid and β 2-Glycoprotein complexes antibodies adsorbed on microplate REF DCE002/11403-0
7. TMB Substrate (1 vial, 15 mL)
 H_2O_2 -TMB (0.26 g/L) (avoid any skin contact)
ProClin <0.0015% REF DCE004-0
8. Stop Solution (1 vial, 15 mL)
Sulphuric acid 0.15M (avoid any skin contact)
REF DCE005-0
9. 10X Conc. Wash Solution (1 vial, 50 mL)
Phosphate buffer 0.2M pH 7.4. ProClin >0.0015%
REF DCE054-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser
Precision Pipetting Devices
Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.

- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- Some reagents (calibrators, controls and sample diluent) contain small amounts of Sodium Azide (NaN_3) <0.1% which may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through with large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/ H_2O_2 to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 – 8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22 – 28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of

samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate

- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinisation and decontamination are effective in avoiding contamination of the conjugate; **this procedure is highly recommended when the kit is processed using analysers which are not equipped with disposable tips.**
For this purpose, Diametra supplies a separate decontamination reagent for cleaning needles.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the calibrators are stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2 – 8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum or plasma samples.

Sample Storage	Duration
2 – 8 °C	5 days
-20°C	6 months

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer.

The calibrators contain both IgG and IgM antibodies, are ready for use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄
AU/mL	0	5	10	20	80

The controls are ready to use; the concentration is printed on the label.

9.2. Preparation of the Conjugate

The conjugates are ready to use. Mix gently, for 5 minutes, on a roller mixer.

9.3. Preparation of the Wash Solution

Dilute the content of the vial “10X Conc. Wash Solution” with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

The determination of IgG or IgM class autoantibodies against anti- phospholipids can be performed in human serum or plasma samples.

All serum and plasma samples must be pre-diluted with sample diluent 1:100; for example, 10 µL of sample may be diluted with 990 µL of sample diluent.

Fasting samples are not necessary and no special sample preparations are required.

Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation. Testing of heat-inactivated sera is not recommended.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.5. Procedure

- **Allow all reagents to reach room temperature (22 – 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve

(C₀-C₄), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Calibrator C ₀ -C ₅	100 µL		
Control		100 µL	
Diluted sample		100 µL	

Incubate 60 minutes at room temperature (+22 – 28°C). Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

Conjugate (IgG or IgM)	100 µL	100 µL	
---------------------------	--------	--------	--

Incubate 30 minutes at room temperature (+22 – 28°C). Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Washing: follow the same indications of the previous point.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at 22 – 28°C for 15 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake gently the microplate.
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁵.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required**. A smoothed spline fit including Calibrator 0 can be used. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

The Positive Control are intended to monitor for substantial reagent failure and they will not ensure precision at the assay cut-off.

12. EXPECTED VALUES

In a normal range study with serum samples from healthy blood donors the following ranges have been established with the Anti Phospholipid Screen test:

	IgG (AU/mL)	IgM (AU/mL)
Normal	< 10	< 10
Elevated	≥ 10	≥ 10

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

13. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

13.1. Precision and reproducibility

Precision and reproducibility are evaluated by eight replicates of two positive samples by two different runs with two different lots. Dispensing and washing operations were performed manually by an operator.

The results in terms of standard deviation and coefficient of variation were below:

Sample	IgG			
	1		2	
	SD	CV%	SD	CV%
Intra-assay	1.03	5.9	1.31	7.4
Inter-assay	0.26	9.2	5.25	11.7

	IgM			
	1		2	
	SD	CV%	SD	CV%
Intra-assay	0.61	7.6	1.97	5.9
Inter-assay	0.15	7.1	2.98	6.6

13.2. Sensitivity

The clinical sensitivity of Anti Phospholipid Screen IgG assay is 92.3%.

The clinical sensitivity of Anti Phospholipid Screen IgM assay is 68.8%.

13.3. Specificity

The clinical specificity of Anti Phospholipid Screen IgG assay is 84.6%.

The clinical specificity of Anti Phospholipid Screen IgM assay is > 99.9%.

13.4. Detection Limit

The lowest concentration that can be distinguished from Calibrator zero is 0.03 AU/mL for IgG and 0.16 AU/mL for IgM.

14. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁶. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.
- The presence of immune complexes or other immunoglobulin aggregates in the patient sample may cause an increased level of non-specific binding and produce false positives in this assay.

15. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

16. BIBLIOGRAPHY

1. Hughes G.R.V., Harris E.N., Gharavi A.E.: The Anticardiolipin Syndrome. J. Rheumat. 13, 3: 486-489, 1986.
2. Harris E.N. et al.: Evaluation of the anti-Cardiolipin antibody test: report of an international Workshop held 4 April 1986. Clin. Exp. Immunol. 68: 215-222, 1987.
3. Domke N., Siebert G.: Phospholipidantikörper und ihre klinische Bedeutung. Zeitschrift für Klinische Medizin 16: 1399-1401, 1988.
4. Pengo V. et al.: Immunological Specificity and Mechanism of Action of IgG Lupus Anticoagulants. Blood, Vol. 70. N. 1: 69-76, 1987.
5. Basic QC Practices On-line Course; <http://www.Westgard.com>.
6. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. Clin Chem, 34, 1988, pp 27-33

17. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

18. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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