DCM099-8
Ed. 01/2024

ANA Screen

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ANA Screen

for routine analysis

Qualitative determination of IgG class antibodies against nuclear antigens in human serum or plasma

IVD

LOT

See external label

2°C 8°C

$\Sigma = 96$ tests

REF DKO099

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

The ANA Screen assay is a manual *in vitro* diagnostic device intended for the qualitative determination of IgG antibodies directed against Sm (Smith), RNP/Sm, Scl-70, SS-A (Ro) (52kDa e 60kDa), SS-B (La), Jo1, U1-SmRNP, CENP-B, dsDNA and Histones in human serum or plasma from an adult population. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of systemic rheumatic autoimmune diseases.

2. CLINICAL SIGNIFICANCE

Anti-nuclear autoantibodies (ANAs) are a heterogeneous group of autoantibodies that recognise nuclear macromolecules and their complexes; antibodies to certain cytoplasmic proteins might also fall into this category¹⁻⁶. These antibodies commonly occur in the sera of patients with systemic autoimmune rheumatic diseases (SARDs) and can bind to DNA, RNA and proteins, as well as complexes of nucleic acids with proteins that are essential for different intracellular functions such as replication and transcription.

ANAs' targets can be divided into two groups in the context of systemic rheumatic autoimmune diseases: antibodies that recognise DNA, histones and nucleosomes^{1,2,5,6}; and those that bind to complexes of RNA with RNA-binding proteins (RBPs) and small nuclear ribonucleoproteins (snRNPs)¹.

Autoantibodies directed to nuclear antigens are serological hallmarks of SARDs and their detection is used by clinicians to support the diagnosis of disorders such as SLE, Sjögren syndrome, systemic sclerosis, polymyositis and mixed connective tissue disease. It is well documented that these are complex disorders and characterised by the presence of a wide spectrum of autoantibodies. Furthermore, it has widely been shown that ANAs are also frequently found in the sera of patients with non-rheumatic diseases and healthy subjects that do not develop a rheumatic disease.

Due to the complexity of SARDs, ANA tests should always be considered together with the evaluation of the patient's history, clinical manifestations and other laboratory tests to diagnose these conditions.

3. PRINCIPLE OF THE METHOD

The ANA Screen test is an indirect enzyme immunometric assay (ELISA) based on the binding of present antibodies with Sm, RNP/Sm, SS-A (Ro), SS-B (La), Scl-70, Jo1, U1-SmRNP, CENP-B, dsDNA and Histones antigens coated on the microplates. Any antibodies present in calibrators, controls or prediluted patient samples bind to the inner surface of the wells. After a 30 minutes incubation the microplate is washed with wash buffer for removing non-reactive serum components.

An anti-human-IgG horseradish peroxidase conjugate solution recognises IgG class antibodies bound to the immobilised antigens. After a 30 minutes incubation any excess enzyme conjugate, which is not specifically bound is washed away with wash buffer.

A chromogenic substrate solution containing TMB is dispensed into the wells. After 15 minutes of incubation the colour development is stopped by adding the stop solution. The solution turns yellow at this point. The level of colour is directly proportional to the concentration of IgG antibodies present in the original sample.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrator (1 vial, 1.2 mL)

Phosphate buffer 0.1M, NaN₃ < 0.1%, BSA 3%

CAL1

REF DCE002/09906-0

2. Controls (2 vials, 1.2 mL each)

Phosphate buffer 0.1M, NaN₃ < 0.1%, BSA 3%

Negative Control

REF DCE045/09901-0

Positive Control

REF DCE045/09902-0

3. Sample Diluent (1 vial, 100 mL)

Phosphate buffer 0.1M, NaN₃ < 0.1%

REF DCE053-0

4. Conjugate (1 vial, 15 mL)

Anti h-IgG conjugated with horseradish peroxidase (HRP), BSA 0.1%, ProClin >0.0015%

REF DCE002/09902-0

5. Coated Microplate (1 breakable microplate)

Microplate coated with ANA antigens

REF DCE002/09903-0

6. TMB Substrate (1 vial, 15 mL)

H₂O₂ -TMB (0,26 g/L) (avoid any skin contact), ProClin <0.0015%

REF DCE004-0

7. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

8. 10X Conc. Wash Solution (1 vial, 50 mL)
Phosphate buffer 0.2M pH 7.4, ProClin <0.0015%
REF DCE054-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser.

Precision Pipetting Devices Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- **Classification according to Regulation (EC) No. 1272/2008 [CLP]**

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- Some reagents (calibrator, controls and sample diluent) contain small amounts of Sodium Azide (<0.1%) as preservative which may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately **validated for its intended use/purpose**
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- **WARNING: the conjugate reagent is designed to ensure maximum dose sensitivity and may be contaminated by external agents if not used properly;** therefore, it is recommended to use disposable consumables (tips, bottles, trays, etc.). For divided doses, take the exact amount of conjugate needed and do not re-introduce any waste product into the original bottle. In addition, **for doses dispensed with the aid of automatic and semi-automatic devices,** before using the conjugate, it is advisable to clean the fluid handling system, ensuring that the procedures of washing, deproteinization and decontamination are effective in avoiding contamination of the conjugate. This **procedure is highly recommended when the kit is processed using analysers which are not equipped with disposable tips.**
For this purpose, DiaMetra supplies a separate decontamination reagent for cleaning needles.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated in the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2 – 8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2-8°C	96 hours
Freeze/thaw cycles	3 cycles

9. PROCEDURE

9.1. Preparation of the Calibrator (C₁)

Calibrator is ready to use; the concentration of the Calibrator is printed on the label.

The Controls are ready to use.

9.2. Preparation of the Conjugate

The conjugate is ready to use. Mix gently, for 5 minutes, on a roller mixer.

9.3. Preparation of the Wash Solution

Dilute the contents of each vial of the buffered wash solution concentrate (10X) with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

For determination of ANA antibodies, human serum or plasma are the preferred sample matrixes.

All serum and plasma samples must be prediluted with sample diluent 1:100; for example, 10 µL of sample may be diluted with 990 µL of sample diluent.

Fasting samples are not necessary and no special preparations are required.

Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₁), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Controls	Blank
Calibrator C ₁	100 µL		
Controls		100 µL	
Diluted Sample		100 µL	

Incubate for 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

Conjugate	100 µL	100 µL	
Incubate for 30 minutes at room temperature (22-28°C). Remove the content from each well, wash the wells 3 times with 300 µL of diluted wash solution.			
Washing: follow the same indications of the previous point.			
TMB Substrate	100 µL	100 µL	100 µL

Incubate for 15 minutes in the dark at room temperature (22-28°C).

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake the microplate gently.

Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁷.

11. CALCULATION OF RESULTS

Determine the mean absorbance for each duplicate sample.

To obtain sample value: divide the mean absorbance of the sample by the absorbance of the Calibrator, then multiply by the concentration of the Calibrator, as shown below:

$$\text{Sample Conc.} = \frac{OD \text{ sample}}{OD \text{ C1}} \times \text{Conc. C1}$$

Reactivity is not connected in a linear way to the amount of antibodies present. Although an increase or a decrease in the concentration of antibodies results in increased or decreased responsiveness, the change is not in proportion (e.g., doubling the concentration of antibodies does not lead to a doubling of responsiveness).

To achieve greater accuracy in the determination of antibodies it is recommended to test serial dilutions of the sample. The final dilution that is positive in the test is the antibody concentration of the patient.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

12. METROLOGY AND TRACEABILITY

The ANA Screen assay is traceable to the Antinuclear Antibodies (ANA) Reference Standards by Autoantibody Standardisation Committee (ASC) - Center for Disease Control (CDC) and WHO anti-dsDNA 15/174.

13. INTERPRETATION OF RESULTS

ANA Screen (AU/ mL)	Interpretation
< 25	The sample should be considered negative
≥ 25	The sample should be considered positive

The above index and interpretation should be considered as guidelines only. Positive results should be verified concerning the entire clinical status of the patient. Also, every decision for therapy should be taken on an individual patient basis.

It is recommended that each laboratory establishes its own normal and pathological ranges of serum anti-ANA.

14. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

14.1. Diagnostic Sensitivity and Specificity

The sensitivity and specificity were determined with guidance from CLSI EP-24 "Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves" using 50 negative and 53 positive samples run on a single reagent lot.

		DKO099		Total
		Positive	Negative	
True state	Positive	53	0	53
	Negative	0	50	50
Total		53	50	103

Diagnostic sensitivity: 100%

Diagnostic specificity: 100%

14.2. Precision

Precision of the ANA Screen kit was determined by performing a complex precision study.

Repeatability: A total of A total of 3 serum samples were assayed using 2 lots of reagents in 10 replicates per sample, by 3 operators, once a day for 5 days. Data from one representative lot is shown below:

Repeatability results from one representative lot:

Sample	n	+ve results	-ve results	% +ve results	% -ve results
Cut off	150	82	68	54.7	45.3
Negative	150	1	149	0.7	99.3
Positive	150	150	0	100.0	0.0

Reproducibility results for the combined data from two lots is shown below:

Sample	n	+ve results	-ve results	% +ve results	% -ve results
Cut off	300	163	137	54.3	45.7
Negative	300	3	297	1.0	99.0
Positive	300	300	0	100.0	0.0

14.3. Serum-plasma study

The ANA Screen matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI GP34-A "Validation and Verification of Tubes for Venous and Capillary Blood. Specimen Collection". A total of 20 samples (16 native, 4 contrived) were evaluated.

Sample Type	No. Samples	Results vs Serum	
		False pos	False neg
SST	20	0	0
K2 EDTA	20	0	0
Lithium heparin	20	0	0
Sodium heparin	20	0	0

14.4. Interferents

The following substances do not interfere in the ANA Screen assay when the concentrations presented in the following table are below the stated threshold.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Haemoglobin	200 mg/dL
Total Protein	10 g/dL
Triglyceride	500 mg/dL

15. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁸. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

16. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

17. BIBLIOGRAPHY

1. Pisetsky DS. Antinuclear antibody testing - misunderstood or misbegotten? *Nat Rev Rheumatol*. 2017 Aug;13(8):495-502.
2. Pisetsky DS. Antinuclear antibodies in rheumatic disease: a proposal for a function-based classification. *Scand J Immunol*. 2012 Sep;76(3):223-8.
3. Fritzler MJ. Clinical relevance of autoantibodies in systemic rheumatic diseases. *Mol Biol Rep*. 1996;23(3-4):133-45.
4. von Mühlen CA, Tan EM. Autoantibodies in the diagnosis of systemic rheumatic diseases. *Semin Arthritis Rheum*. 1995 Apr;24(5):323-58.
5. Defendenti C, Atzeni F, Spina MF, Grosso S, Cereda A, Guercilena G, Bollani S, Saibeni S, Puttini PS. Clinical and laboratory aspects of Ro/SSA-52 autoantibodies. *Autoimmun Rev*. 2011 Jan;10(3):150-4.
6. Aggarwal A. Role of autoantibody testing. *Best Pract Res Clin Rheumatol*. 2014 Dec;28(6):907-20.
7. Basic QC Practices On-line Course; <http://www.Westgard.com>.
8. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

18. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

19. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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