DCM077-8
Ed. 03/2024

C-PEPTIDE ELISA

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BIBLIOGRAFIA

- Eastham R.D: Biochemical Values in Clinical Medicine, 7th Ed. Bristol. England. Jonh Wright & Sons, Ltd;. (1985).
- Gerbitz, V.K.D, J.Clin.Chem.Biochem. 18, 313-326 (1980)
- Boehm TM, et al Diabetes Care 479-490. (1079)
- National Committee for Clinical Laboratory Standards. Procedure for the collection of diagnostic blood specimens by venipuncture: approved standards. 4th Ed. NCCLS Document H3-A4, Wayne, PA(1988).
- Turkinton RW, et al Archive of Internal Med. 142

- Sacks BD Carbohydrates in Burtis, C.A. and Ashwood, AR (Eds) Tietz Textbook OF Clinical Chemistry.2nd Ed. Philadelphia w.B. Saunders Co. 1994
- Kahn CR, et al Diabetes Care 2, 283 – 295 (1979)

Ed. 03/2024

DCM077-8

DCM077-8
Ed. 03/2024

C-PEPTIDE ELISA

for routine analysis

Direct immunoenzymatic determination of C-Peptide level in human serum or plasma

IVD

LOT

See external label

2°C  8°C

$\Sigma = 96$ tests

REF DKO077

INTENDED USE

C-Peptide ELISA kit is a direct solid phase enzyme immunoassay for quantitative determination of C-Peptide in human serum or plasma.

C-Peptide ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

C-Peptide is the abbreviation for connecting peptide, it is a 31-aminoacid peptide. C-Peptide of insulin is the C-terminal cleavage product produced during processing of the insulin prohormone to the mature insulin molecule. Proinsulin is split when it is released from the pancreas into the blood - one C-Peptide for each insulin molecule. C-Peptide is devoid of any biological activity but appears to be necessary to maintain the structural integrity of Insulin.

In-vitro determination of Insulin and C-Peptide level help in differential diagnosis of liver disease, acromegaly, Cushing syndrome, familial glucose intolerance, Insulinemia, renal failure, ingestion of accidental oral hypoglycemic drugs or C-Peptide induced factitious hypoglycemia.

Newly diagnosed diabetes patient often get their C-Peptide levels measured, to find if they are type 1 diabetes or type 2 diabetes. The pancreas of patients with type 1 diabetes is unable to produce insulin and they will therefore usually have a decreased level of C-Peptide, while C-Peptide levels in type 2 patients is normal or higher than normal. Measuring C-Peptide in patients injecting insulin can help to determine how much of their own natural insulin these patients are still producing.

C-Peptide assays may be analytically more sensitive than insulin assays. Measurement of the C-Peptide may be useful in evaluating endogenous insulin secretion in a variety of clinical conditions. Insulin and C-Peptide are secreted into portal circulation in equimolar concentrations, fasting levels of C-Peptide are 5 – 10 fold higher than those of Insulin owing to the longer half-life of C-Peptide. The liver does not extract C-Peptide however; it is removed from the circulation by degradation in the kidneys with a fraction passing out unchanged in urine. Hence the


urine C-Peptide levels correlate well with fasting C-Peptide levels in serum.

2. PRINCIPLE

The essential reagents required for an immunoenzymometric assay include high affinity and specificity antibodies (enzyme-linked and immobilized), with different and distinct epitope recognition, in excess, and native antigen. In this procedure, the immobilization takes place during the assay at the surface of a microplate well through the interaction of the streptavidin coated on the well and exogenously added biotinylated monoclonal anti C-Peptide antibody (Ab).

Upon mixing the monoclonal biotinylated antibody, the enzyme-labelled antibody and a serum containing the native antigen (Ag), a reaction results between the native antigen and the antibodies, without competition or steric hindrance, to form a soluble sandwich complex.

The interaction is illustrated in the following equation:

BtAb(m) = biotinylated monoclonal antibody (excess quantity)

AgN = native antigen (variable quantity)

EAb = enzyme labelled antibody (excess quantity)


EAb – AgN – BtAb(m) = antigen-antibodies sandwich complex

k_a = rate constant of association

k_{-a} = rate constant of dissociation

Simultaneously, the complex is deposited into the well through the high affinity reaction of streptavidin and biotinylated antibody.

This interaction is illustrated below:

StreptavidinC.W. = Streptavidin immobilized on well

Immobilized complex = sandwich complex bound to the well

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by a washing step. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is quantitated by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 2 mL each, lyophilized), ProClin <0.0015%

CAL0	REF DCE002/7706-0
CAL1	REF DCE002/7707-0
CAL2	REF DCE002/7708-0
CAL3	REF DCE002/7709-0
CAL4	REF DCE002/7710-0
CAL5	REF DCE002/7711-0

2. Conjugate (1 vial, 13 mL)

Antibodies anti C-Peptide conjugated with horseradish peroxidase (HRP) and anti biotinylated C-Peptide, ProClin <0.0015% REF DCE002/7702-0

3. Coated Microplate (1 breakable microplate)

Microplate coated with streptavidin REF DCE002/7703-0

4. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB 0.26 g/L (avoid any skin contact), ProClin <0.0015% REF DCE004-0

5. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact) REF DCE005-0

6. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L; Tween-20 55 g/L, ProClin >0.0015% REF DCE006B-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Note

Store all reagents between 2-8°C in the dark.

Open the bag of reagent 3 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until the expiry date of the kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.

- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators should be handled in the same manner as potentially infectious material.
- Some reagents (wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- This method allows the quantitative determination of C-Peptide from 0.2 to 10.0 ng/mL.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.

- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of the Calibrator (C₀...C₅)

The Calibrators were calibrated using a reference preparation, which was assayed against the WHO 1st IRR 84/510.

The Calibrators have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
ng/mL	0	0,2	1,0	2,0	5,0	10,0

Reconstitute each Calibrator with 2 mL of distilled or deionised water.

Once reconstituted, the Calibrators are stable 7 days at 2-8°C.

In order to store for a longer period aliquot the reconstituted calibrators in vials and store them at -20°C (stable 6 months). Do not freeze and thaw more than once.

6.2. Preparation of Wash Solution

Dilute the content of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

6.3. Preparation of the Sample

Follow the good laboratory procedures for handling blood products.

For accurate comparison to established normal values, a fasting morning serum sample should be obtained.

To obtain the serum, the blood should be collected in a venipuncture tube without additives or anti-coagulants; allow the blood to clot; centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, they must be stored at -20°C for up to 30 days. Avoid repetitive freezing and thawing.

Patient specimens with C-Peptide concentrations above 10 ng/mL may be diluted (for example 1:10 or higher) with the Calibrator 0 and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.

6.4. procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample	Blank
Calibrator C ₀ -C ₅	50 µL		
Sample		50 µL	
Conjugate	100 µL	100 µL	
Incubate at room temperature (22-28°C) for 2 hours. Remove the content from each well. Wash the wells 3 times with 300 µL of diluted Wash Solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22-28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

7. QUALITY CONTROL

Each laboratory should assay controls at levels in the low, medium and high ranges of the dose response curve for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

8.1. Note

The optical densities (O.D.) of some calibrators and samples may be higher than 2.0, in such a case, they could be out of the measurement range of the microplate reader. It is therefore necessary, for O.D. higher than 2.0, to perform a reading at 405 nm (=wavelength of peak shoulder) in addition to 450 nm (peak wavelength) and 620 nm (reference filter for the subtraction of interferences due to the plastic). For microplate readers unable to read the plate at 3 wavelengths at the same time, it is advisable to proceed as follows:

- Read the microplate at 450 nm and at 620 nm.
- Read again the plate at 405 nm and 620 nm.
- Find out the wells whose OD at 450 nm are higher than 2.0

- Select the corresponding ODs read at 405 nm and multiply these values at 405 nm by the conversion factor 3.0 (where $OD_{450}/OD_{405} = 3.0$), that is: $OD_{450\text{ nm}} = OD_{405\text{ nm}} \times 3.0$.

Warning: The conversion factor 3.0 is suggested only. For better accuracy, the user is advised to calculate the conversion factor specific for his own reader.

8.2. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve and of each sample.

8.3. Calibration curve – Automatic method

Use the 4 parameters logistic – preferred – or the smoothed cubic spline function as calculation algorithm.

8.4. Calibration curve – Manual method

A dose response curve is used to ascertain the concentration of C-Peptide in unknown specimens. Record the OD obtained from the printout of the microplate reader.

Plot the OD for each duplicate calibrator versus the corresponding C-Peptide concentration in ng/mL on linear graph paper (do not average the duplicates of the calibrators before plotting).

Draw the best-fit curve through the plotted points.

To determine the concentration of C-Peptide for an unknown, locate the average absorbance of the duplicates for each unknown on the vertical axis of the graph, find the intersecting point on the curve, and read the concentration (in ng/mL) from the horizontal axis of the graph (the duplicates of the unknown may be averaged as indicated).

9. REFERENCE VALUES

C-Peptide values are consistently higher in plasma than in serum; thus, serum is preferred.

In concordance with the published literature the following ranges have been assigned.

	C-Peptide
Adult (Normal)	0,7 – 1,9 ng/mL

Please pay attention to the fact that the determination of a range of expected values for a “normal” population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay

Within run variation was determined by replicate measurements (16x) of three different control sera in one assay. The within assay variability is $\leq 6.2\%$.

10.1.2. Inter Assay

Between run variation was determined by replicate measurements (20x) of three different control sera in different lots. The between assay variability is $\leq 10.0\%$.

10.2. Sensitivity

The lowest detectable concentration of C-Peptide that can be distinguished from the Calibrator 0 is 0.01 ng/mL at the 95 % confidence limit.

10.3. Specificity

The cross reactivity of the C-Peptide ELISA method to selected substances was evaluated by adding the interfering substance. The cross reactivity was calculated by deriving a ratio between dose of interfering substance to dose of C -Peptide needed to produce the same absorbance.

Cross Reagent	Conc. Tested	Obtained	Cross Reactivity
C-Peptide	---	---	100 %
Insulin	10000 μ IU/mL	N.D.	Not Detected
Proinsulin	1000 ng/mL	N.D.	Not Detected

10.4. Correlation

Dia.Metra C-Peptide ELISA was compared to another commercially available C-Peptide assay. 194 Serum samples were analysed according in both test systems.

The linear regression curve was calculated

$$y = 1.012 x + 0.025$$

$$r^2 = 0.991$$

y = C-Peptide commercial kit

x = C-Peptide Dia.Metra kit

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

diagnostic blood specimens by venipuncture: approved standards. 4th Ed. NCCLS Document H3-A4, Wayne, PA(1988).

- Turkinton RW., Estkowski A., Link M. Secretion of Insulin dependence of connecting peptide; a predictor of insulin or dependence of obese diabetics. *Archive of Internal Med.* 142
- Sacks BD: Carbohydrates in Burtis, C.A. and Ashwood, AR (Eds) *Tietz Textbook OF Clinical Chemistry.* 2nd Ed. Philadelphia W. .B. Saunders Co. 1994
- Kahn CR, Rosenthal AS, Immunologic reactions to insulin, insulin allergy, insulin resistance and autoimmune insulin syndrome. *Diabetes Care* 2, 283 – 295 (1979)

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- Eastham R.D: Biochemical Values in Clinical Medicine, 7th Ed. Bristol. England. Jonh Wright & Sons, Ltd;. (1985).
- Gerbitz, V.K.D, J.Clin.Chem.Biochem. 18, 313-326 (1980)
- Boehm TM, et al Diabetes Care 479-490. (1079)
- National Committee for Clinical Laboratory Standards. Procedure for the collection of diagnostic blood specimens by venipuncture: approved standards. 4th Ed. NCCLS Document H3-A4, Wayne, PA(1988).
- Turkinton RW, et al Archive of Internal Med. 142
- Sacks BD Carbohydrates in Burtis, C.A. and Ashwood, AR (Eds) Tietz Textbook OF Clinical Chemistry. 2nd Ed. Philadelphia w.B. Saunders Co. 1994
- Kahn CR, et al Diabetes Care 2, 283 – 295 (1979)

Ed. 03/2024

DCM077-8

	DE ES FR EN IT PT	<i>In vitro</i> Diagnostikum Producto sanitario para diagnóstico <i>In vitro</i> Dispositif medical de diagnostic <i>in vitro</i> <i>In vitro</i> Diagnostic Medical Device Dispositivo medico-diagnostico <i>in vitro</i> Dispositivos medicos de diagnostico <i>in vitro</i>		DE ES FR EN IT PT	Hergestellt von Elaborado por Fabriqué par Manufacturer Produttore Produzido por
	DE ES FR EN IT PT	Achtung, Begleitdokumente Precaución, consulte los documentos adjuntos Attention, veuillez consulter les documents d'accompagnement Caution, consult accompanying documents Attenzione, consultare la documentazione allegata Atenção, consultar os documentos de acompanhamento	 yyyy-mm	DE ES FR EN IT PT	Herstellungs datum Fecha de fabricacion Date de fabrication Date of manufacture Data di produzione Data de produção
 yyyy-mm-dd	DE ES FR EN IT PT	Verwendbar bis Estable hasta (usar antes de último día del mes) Utiliser avant (dernier jour du mois indiqué) Use by (last day of the month) Utilizzare prima del (ultimo giorno del mese) Utilizar (antes ultimo día do mês)		DE ES FR EN IT PT	Biogefährdung Riesco biológico Risque biologique Biological risk Rischio biologico Risco biológico
	DE ES FR EN IT PT	Gebrauchsanweisung beachten Consultar las instrucciones Consulter le mode d'emploi Consult instructions for use Consultare le istruzioni per l'uso Consultar instruções para uso		DE ES FR EN IT PT	Chargenbezeichnung Codigo de lote Numero de lot Batch code Codice del lotto Codigo do lote
 $\Sigma = xx$	DE ES FR EN IT PT	Ausreichend für "n" Tests Contenido suficiente para "n" tests Contenu suffisant pour "n" tests Contains sufficient for "n" tests Contenuto sufficiente per "n" saggi Contém o suficiente para "n" testes		DE ES FR EN IT PT	Inhalt Contenido del estuche Contenu du coffret Contents of kit Contenuto del kit Conteúdo do kit
 Max Min	DE ES FR EN IT PT	Temperaturbereich Limitación de temperatura Limites de température de conservation Temperature limitation Limiti di temperatura Temperaturas limites de conservação		DE ES FR EN IT PT	Bestellnummer Número de catálogo Références du catalogue Catalogue number Numero di Catalogo Número do catálogo
	DE ES FR EN IT PT	Vor direkter Sonneneinstrahlung schützen Mantener alejado de la luz solar Tenir à l'écart de la lumière du soleil Keep away from sunlight Tenere lontano dalla luce solare Mantenha longe da luz solar			