DCM054-10  
Ed. 03/2024

## CA125 ELISA

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DCM054-10  
Ed. 03/2024

# CA125 ELISA

Direct immunoenzymatic determination of CA125 in human serum or plasma.

for routine analysis

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO054

## INTENDED USE

Dia.Metra CA125 ELISA kit is an immunoenzymatic colorimetric method for the quantitative determination of the tumoral antigen CA125 concentration in human serum or plasma.

CA125 ELISA kit is intended for laboratory use only.

## 1. CLINICAL SIGNIFICANCE

CA125 is a tumour marker, it is a glycoprotein as high molecular weight entity (MW > 200,000).

It is best known as a marker for ovarian cancer, but it may also be elevated in other malignant cancers, including those originating in the endometrium, fallopian tubes, lungs, breast and gastrointestinal tract. CA125 may also be elevated in a number of relatively benign conditions, such as endometriosis, several diseases of the ovary, and pregnancy. It also tends to be elevated in the presence of any inflammatory condition in the abdominal area, both cancerous and benign. Thus, CA125 is not perfectly specific for cancer nor is it perfectly sensitive since not every patient with cancer will have elevated levels of CA125 in the blood.

CA125 is clinically approved for following the response to treatment and predicting prognosis after treatment. It is especially useful for detecting the recurrence of ovarian cancer. Its potential role for the early detection of ovarian cancer is controversial and has not yet been adopted for widespread screening efforts in asymptomatic women.

Elevated levels in post-menopausal women are usually an indication that further screening is necessary. In pre-menopausal women, the test is less reliable as values are often elevated due to a number of non-cancerous causes.

## 2. PRINCIPLE


The essential reagents required for an immunoenzymatic assay include high affinity and specificity antibodies (enzyme and immobilised) with different and distinct epitope recognition, in excess, and native antigen.

In this method, CA125 calibrators, patient specimens and/or controls containing the native antigen are first added to streptavidin coated wells. Biotinylated monoclonal and enzyme labeled antibodies are

added and the reactants mixed: these antibodies have high affinity and specificity and are detected against distinct and different epitopes of CA125.

Reaction between the various CA125 antibodies and native CA125 occurs in the microwells without competition or steric hindrance forming a soluble sandwich complex.

The interaction is illustrated by the following equation:

**BAb(m)** = biotinylated monoclonal antibody (excess quantity)

**AgCA125** = native antigen (variable quantity)


**EAb** = enzyme labelled antibody (excess quantity)

**EAb-AgCA125-BtAb(m)** = antigen-antibodies sandwich complex

**ka** = rate constant of association

**k-a** = rate constant of dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

**Streptavidin C.W.** = Streptavidin immobilized on well  
**Immobilized complex** = sandwich complex bound to the well.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. By utilizing several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

### 3. REAGENTS, MATERIALS AND INSTRUMENTATION

#### 3.1. Reagents and materials supplied in the kit

##### 1. Calibrators (6 vials, 1 mL each)

ProClin <0.0015%

|      |                   |
|------|-------------------|
| CAL0 | REF DCE002/5406-0 |
| CAL1 | REF DCE002/5407-0 |
| CAL2 | REF DCE002/5408-0 |
| CAL3 | REF DCE002/5409-0 |
| CAL4 | REF DCE002/5410-0 |
| CAL5 | REF DCE002/5411-0 |

##### 2. Conjugate (1 vial, 13 mL)

Antibody anti CA125 conjugated with horseradish peroxidase (HRP) and antibody anti CA125 biotinylated, ProClin <0.0015%

REF DCE002/5402-0

##### 3. Coated Microplate (1 microplate breakable)

Microplate coated with Streptavidin

REF DCE002/5403-0

##### 4. TMB Substrate (1 vial, 15 mL)

H<sub>2</sub>O<sub>2</sub> -TMB (0.26 g/L) (avoid any skin contact), ProClin <0.0015%

REF DCE004-0

##### 5. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

REF DCE005-0

##### 6. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L, Tween-20 55 g/L, ProClin >0.0015%

REF DCE006B-0

#### 3.2. Reagents necessary not supplied

Distilled water.

#### 3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

#### Note

Store all reagents at 2÷8°C in the dark.

Open the bag of reagent 3 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, it is stable until the expiry date of the kit.

### 4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators

should be handled in the same manner as potentially infectious material.

- Some reagents (wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- **Classification according to Regulation (EC) No. 1272/2008 [CLP]**

Skin sensitivity, Category 1

Contains: ProClin 300

Warning

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H<sub>2</sub>O<sub>2</sub> to directed sunlight, metals or oxidants. Do not freeze the solution.
- CA125 has a low clinical sensitivity and specificity as a tumour marker. Clinically an elevated CA125 value alone is not of diagnostic value as a test for cancer and should only be used in conjunction with other clinical manifestations (observations) and diagnostic parameters.

### 5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately tested. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

## 6. PROCEDURE

### 6.1. Preparation of the Calibrators (C<sub>0</sub>...C<sub>5</sub>)

The Calibrators are made using a >99% pure affinity purified preparation of CA125. The preparation is calibrated against Centocor CA125 IRMA test.

The references CA125 Antigen has the following concentrations:

|      | C <sub>0</sub> | C <sub>1</sub> | C <sub>2</sub> | C <sub>3</sub> | C <sub>4</sub> | C <sub>5</sub> |
|------|----------------|----------------|----------------|----------------|----------------|----------------|
| U/mL | 0              | 15             | 50             | 100            | 200            | 400            |

Once opened, the Calibrators are stable 6 months at 2-8°C. The Calibrators are ready to use.

### 6.2. Preparation of Wash Solution

Dilute the contents of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2-8°C.

### 6.3. Preparation of the Sample

The specimens can be human serum or plasma; the usual precautions in the collection of venipuncture samples should be observed. For accurate comparison to established normal values, a fasting morning serum sample should be obtained.

To obtain the serum, the blood should be collected in a plain redtop venipuncture tube without additives or anti-coagulants; allow the blood to clot; centrifuge the specimen to separate the serum from the cells.

Samples may be refrigerated at 2-8°C for a maximum period of 5 days. If the specimens cannot be assayed within this time, they may be stored at temperatures of -20°C for up to 30 days. Avoid repetitive freezing and thawing.

Patient specimens with CA125 concentrations above 400 U/mL may be diluted (for example 1:10 or higher) with normal male serum (CA125 < 5 U/mL) and re-assayed. The sample's concentration is obtained by multiplying the result by the dilution factor.

When assayed in duplicate, 0.050 mL of the specimen is required.

## 6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>), two for each Control, two for each sample, one for Blank.

| Reagent                                                                                                                                                                                                                                                                                                                                                                                                                        | Calibrator | Sample | Blank  |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------|--------|
| Calibrator<br>C <sub>0</sub> -C <sub>5</sub>                                                                                                                                                                                                                                                                                                                                                                                   | 25 µL      |        |        |
| Sample                                                                                                                                                                                                                                                                                                                                                                                                                         |            | 25 µL  |        |
| Conjugate                                                                                                                                                                                                                                                                                                                                                                                                                      | 100 µL     | 100 µL |        |
| Incubate at room temperature (22±28°C) for 1 hour. Remove the content from each well; wash the wells 3 times with 300 µL of diluted Wash Solution.<br><b>Important note:</b> during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.<br><b>Automatic washer:</b> if you use automated equipment, wash the wells at least 5 times. |            |        |        |
| TMB<br>Substrate                                                                                                                                                                                                                                                                                                                                                                                                               | 100 µL     | 100 µL | 100 µL |
| Incubate at room temperature (22±28°C) for 15 minutes in the dark.                                                                                                                                                                                                                                                                                                                                                             |            |        |        |
| Stop<br>Solution                                                                                                                                                                                                                                                                                                                                                                                                               | 100 µL     | 100 µL | 100 µL |
| Shake the microplate gently.<br>Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.                                                                                                                                                                                                                                                                              |            |        |        |

## 7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of CA125 for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits.. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

## 8. RESULTS

### 8.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C<sub>0</sub>-C<sub>5</sub>) and of each sample.

### 8.2. Calibration curve

Plot the mean value of absorbance (Em) of the calibrators (C<sub>0</sub>-C<sub>5</sub>) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

### 8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in U/mL.

## 9. REFERENCE VALUES

|                                   | CA125     |
|-----------------------------------|-----------|
| Healthy and non-pregnant subjects | < 35 U/mL |

The serum CA125 is elevated in 1% of normal healthy women, 3% of normal healthy women with benign ovarian diseases, 6% of patients with non-neoplastic conditions (including but not limited to first trimester pregnancy, menstruation, endometriosis uterine fibrosis, acute salpingitis, hepatic diseases and inflammation of peritoneum or pericardium).

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

## 10. PERFORMANCE AND CHARACTERISTICS

### 10.1. Precision

#### 10.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (18x) of three different control sera in one assay. The within assay variability is ≤ 10.0%.

#### 10.1.2. Inter Assay Variation

Between run variation was determined by replicate measurements (20x) of three different control sera in different lots. The between assay variability is ≤ 9.9%.

### 10.2. Sensitivity

The lowest detectable concentration of CA125 that can be distinguished from the Calibrator zero is 1.7 U/mL at the 95% confidence limit.

### 10.3. Specificity

In order to test the specificity of the antibody pair used massive concentrations of possible cross-reactants were added to known serum pools and assayed in parallel with the base sera.

| Substance | U.M.  | Concentration | Cross Reactivity |
|-----------|-------|---------------|------------------|
| CA125     | U/mL  | ---           | 100.0%           |
| PSA       | mg/mL | 1.0           | N.D.             |
| β-HCG     | ng/mL | 250           | N.D.             |
| HCG       | U/mL  | 5000          | N.D.             |
| LH        | U/mL  | 1000          | N.D.             |
| CA-19-9   | U/mL  | 1000          | N.D.             |
| CA-15-3   | U/mL  | 1000          | N.D.             |

According to data, the antibody pair was found to be highly specific for CA125 only.

### 10.4. Correlation with RIA

Dia.Metra CA125 ELISA kit was compared to another commercially available CA125 assay. 142 Serum samples were analysed according in both test systems.

The linear regression curve was calculated:

$$y = 0.987x + 1.019$$

$$r^2 = 0.984$$

$$y = \text{CA125 commercial kit}$$

$$x = \text{CA125 Dia.Metra kit}$$

## 11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

## BIBLIOGRAPHY

- 1) Zamcheck N, Adv Intern Med, 19, 413 (1974)
- 2) Rayncao G, Chu TM, JAMA, 220, 381 (1972)
- 3) Harrison, Principles of Internal Medicine., McGraw Hill Book Company, New York. 12th. Ed.
- 4) Wild D, The Immunoassay Handbook., Stockton Press (1994)

- 5) Hasholzner U, Tumor Diagnosis & Ther 15:114-117 (1994)
- 6) Hasholzner U, Int J Cancer 69:329-34 (1996)
- 7) Ovarian Cancer–NIH Conference. JAMA 273:491-497 (1995)
- 8) Daoud E. Clin Chem 37:1968-74 (1991)
- 9) De Bruijn HWA, Curr Opin Gynecol 9:8-13 (1997)
- 10) Sikorska H, Cancer Detection Preview 12:321-355 (1988)
- 11) NIH CA125 Ann.Inter.Med. 94:407-409 (1981)

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