DCM076-10

Ed. 03/2024

INSULIN ELISA

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18. BIBLIOGRAFIA

1. Khan RMM, Chua ZJY, Tan JC, Yang Y, Liao Z, Zhao Y. From Pre-Diabetes to Diabetes: Diagnosis, Treatments and Translational Research. *Medicina* (Kaunas). 2019 Aug 29;55(9):546
2. Pourmotabbed G, Kitabchi AE. Hypoglycemia. *Obstet Gynecol Clin North Am.* 2001 Jun;28(2):383-400.
3. Göke B. Implications of blood glucose, insulin resistance and beta-cell function in impaired glucose tolerance. *Diabetes Res Clin Pract.* 1998 Jul;40 Suppl:S15-20.
4. Rett K, Wicklmayr M, Mehnert H. What is the clinical significance of insulin resistance? *J Cardiovasc Pharmacol.* 1992;20 Suppl 11:S22-6.
5. Behboudi-Gandevani, S., Amiri, M., Bidhendi Yarandi, R. et al. The impact of diagnostic criteria for gestational diabetes on its prevalence: a systematic review and meta-analysis. *Diabetol Metab Syndr* 11, 11 (2019).
6. John Service F. Hyperinsulinemic hypoglycemia in adults. *Ann Endocrinol (Paris).* 2004 Feb;65(1):88-95.
7. Service FJ. Hypoglycemias. *West J Med.* 1991 Apr;154(4):442-54.
8. Service FJ. Hypoglycemia. *Endocrinol Metab Clin North Am.* 1997 Dec;26(4):937-55.
9. Field JB. Hypoglycemia. Definition, clinical presentations, classification, and laboratory tests. *Endocrinol Metab Clin North Am.* 1989 Mar;18(1):27-43.
10. Basic QC Practices On-line Course; <http://www.Westgard.com>.
11. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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INSULIN ELISA

for routine analysis

Direct immunoenzymatic detection of insulin level in human serum or plasma

IVD

LOT

See external label

2°C  8°C $\Sigma = 96$ tests

REF DKO076

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use
For Laboratory Professional Use

Insulin ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of insulin in human serum or plasma from an adult population. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of diabetes mellitus and hypoglycaemia.

2. CLINICAL SIGNIFICANCE

All body organs and tissues derive their energy primarily from glucose. The sources of body glucose are dietary intake, endogenous glycogen and glucose derived from gluconeogenesis.

During a fasting state, the liver becomes the first source of glucose through glycogenolysis during the first 8 to 12 hours of fasting. During fasting, glucose and thus insulin levels decrease while counterregulatory hormone levels, such as epinephrine which inhibits insulin secretion, begin to rise to help maintain normoglycemia.

During the fed state, insulin is the main hormone in controlling glucose metabolism. After a meal containing carbohydrate, glucose levels increase and stimulate insulin secretion.

Insulin is a hormone produced by beta cells of the pancreas that stabilise blood glucose levels by directing the cells to take up glucose, promoting hepatic storage as glycogen, glucose uptake by muscles and adipose tissue and by inhibiting hepatic glucose production.

In the context of diabetes, beta cell functions are usually impaired as well as insulin secretion and thus blood glucose levels^{1,2}. Insulin resistance leads at first to an increased insulin production to maintain a near-normoglycaemia. Diabetes then results from a further deterioration of insulin resistance of beta cell function, leading to a gradual deterioration of insulin-secreting beta cells and to a complete loss of blood glucose level control^{3,4}. Affected individuals usually pass through a phase of impaired glucose tolerance and increased fasting plasma glucose levels before developing overt type 2 diabetes³.

Insulin secretion is also impaired in the context of gestational diabetes mellitus (GDM). GDM is a heterogenous disorder that occurs in around 7 – 10% of all pregnancies⁵.

insulin levels are also altered in the context of hypoglycaemia, a condition characterised by a blood glucose level lower than normal and abnormally high insulin levels.

Hypoglycaemia is usually caused by diabetes-related treatments, for example an excessively high insulin dose, an inadequate food intake after taking diabetes medication or, in rare cases, by insulinoma, a pancreatic cancer where beta cells secrete abnormal amounts of insulin. Measurements of insulin levels are considered, together with other clinical and laboratory data, to help in the diagnosis of hypoglycaemia in patients that show typical clinical manifestations^{2,6-9}.

3. PRINCIPLE OF THE METHOD

The Insulin ELISA is a sandwich assay based on simultaneous binding of human insulin to two monoclonal antibodies, one immobilised on the microwell plate and the other conjugated with horseradish peroxidase (HRP).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate (H_2O_2) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H_2SO_4) is added. The colour intensity is proportional to the insulin concentration in the sample.

Insulin concentration in the sample is calculated through a calibration curve.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each)

ProClin >0.0015%, BSA 1%

CAL0	REF DCE002/7606-0
CAL1	REF DCE002/7607-0
CAL2	REF DCE002/7608-0
CAL3	REF DCE002/7609-0
CAL4	REF DCE002/7610-0
CAL5	REF DCE002/7611-0

2. Control (1 vial, 1 mL)

Concentration of Control is indicated on the Certificate of Analysis. ProClin >0.0015%, BSA 1%

REF DCE045/7603-0

3. Conjugate (1 vial, 13 mL)

Monoclonal anti- Insulin conjugated with horseradish peroxidase (HRP). ProClin >0.0015%, BSA 0.1%

REF DCE002/7602-0

4. Coated Microplate (1 microplate breakable)

Monoclonal anti- Insulin adsorbed on microplate

REF DCE002/7603-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB (0.26 g/L) (avoid any skin contact)

ProClin <0.0015%

REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

7. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L; Tween-20 55 g/L. ProClin >0.0015%

REF DCE006-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, controls, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.

- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 – 8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22 – 28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2 – 8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	96 hours
Freeze/thaw cycles	3 cycles

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer. The calibrators are ready for use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
µIU/mL	0	3	10	30	80	200

The control is ready to use; the concentration is printed on the label.

9.2. Preparation of the Conjugate

The conjugate is ready to use.

9.3. Preparation of the Wash Solution

Dilute the content of the vial "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated

wash solution to 1000 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

The determination of insulin can be performed in human serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

For accurate comparison to establish normal values, a fasting morning serum sample should be obtained.

Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.5. Procedure

- **Allow all reagents to reach room temperature (22 – 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for the Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C ₀ -C ₅	100 µL		
Sample/Control		100 µL	
Conjugate	100 µL	100 µL	

Incubate 2 h at room temperature (22 – 28 °C). Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at 22 – 28°C for 15 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake gently the microplate.
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories¹⁰.

11. CALCULATION OF RESULTS

The optical densities (ODs) of some calibrators and samples may be higher than 3.0 which may be outside of the measurement range of the microplate reader. Where ODs are higher than 3.0, perform a reading at 405 nm in addition to 450 nm and 620 nm (reference filter for the subtraction of interferences due to the plastic).

For microplate readers unable to read the plate at 3 wavelengths at the same time it is advisable to proceed as follows:

- Read the microplate at 450 nm and at 620 nm.
- Read again the plate at 405 nm and 620 nm.
- Find out the wells whose ODs at 450 nm are higher than 2.0
- Select the corresponding ODs read at 405 nm and multiply these values at 405 nm by the conversion factor

3.0 (where $OD\ 450/OD\ 405 = 3.0$), that is: $OD\ 450\ nm = OD\ 405\ nm \times 3.0$.

Warning: The conversion factor 3.0 is suggested only. For better accuracy, the user is advised to calculate the conversion factor specific for his own reader.

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A smoothed cubic spline curve fit, **including Calibrator 0 is required**. A 4-parameter logistic (4-PL) fit including Calibrator 0 can be used. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

12. MEASURING RANGE

The assay measuring range (AMR) is 2.1 – 200 µIU/mL. Any value that reads below 2.1 µIU/mL should be reported as “< 2.1 µIU/mL”.

13. METROLOGY AND TRACEABILITY

The calibrators of this kit are traceable to the WHO International Standard for insulin (IS-83/500).

14. EXPECTED VALUES

The following ranges were determined using the Insulin ELISA and are provided for information only. The 95 % reference interval for apparently healthy adults collected in a fasting state were calculated by a non-parametric method following guidance from CLSI C28-A3 “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

	n	Median (µIU/mL)	Reference interval (µIU/mL)
Adults	119	8.8	2.1 – 27.7

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	1.2 µIU/mL
Limit of Detection (LoD)	1.8 µIU/mL
Limit of Quantitation (LoQ)	2.1 µIU/mL

15.2. Trueness

Trueness has been demonstrated through a recovery test for the Insulin ELISA with the WHO International Standard for Insulin (IS-83/500).

15.3. Precision

Precision of the Insulin ELISA was determined by performing a complex precision study.

Repeatability: A total of 6 serum samples were assayed in five replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (µIU/mL)	Within run (Repeatability)	
			SD	CV%
1	75	9.2	0.5	5.3%
2	75	22.7	0.9	4.0%
3	75	36.6	1.3	3.6%
4	75	84.7	4.1	4.8%
5	75	117.3	6.8	5.8%
6	75	167.8	7.6	4.5%

Reproducibility: A total of 6 serum samples were assayed in five replicates, once a day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (µIU/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	9.3	1.2	12.6%
2	150	22.5	1.5	6.5%
3	150	36.0	2.9	8.2%
4	150	84.0	5.8	6.9%
5	150	117.8	10.9	9.2%
6	150	169.3	17.2	10.1%

15.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For insulin concentration by Insulin ELISA, the measurement procedure shows linearity for the interval from 1.7 to 251.9 µIU/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$.

15.5. Method comparison

The Insulin ELISA was compared against a commercially available quantitative manual assay following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 41 samples, selected to represent a wide range of insulin concentrations, was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
41	1.0 [0.952 – 1.093]	-0.6 [-1.56 – 0.19]	0.98

15.6. Analytical Specificity

The specificity was assessed with the following cross-reactants.

Cross-reactant	Concentration tested (unit)	Mean %Cross reactivity
C-Peptide	100 ng/mL	-0.1%
Proinsulin	1000 ng/mL	0.1%
Insulin Glargine	1000 µIU/mL	23.7%

The following substances do not interfere with a bias of $> \pm 15\%$ in the Insulin ELISA assay when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	15 mg/dL
Hemoglobin	200 mg/dL
Total Protein	8.3 g/dL
Triglycerides	500 mg/dL

15.7. Serum-plasma study

The Insulin ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI EP9-A3 guidelines. A total of 20 samples (16 native, 4 spiked) to cover the assay range were evaluated. Linear regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	1.0 [0.92 to 1.02]	0.3 [-1.41 to 2.01]	0.99
Lithium Heparin	1.0 [0.93 to 1.09]	-1.8 [-4.56 to 0.87]	0.99
Sodium Heparin	1.0 [0.90 to 1.04]	-1.4 [-3.86 to 1.13]	0.99
EDTA	1.1 [1.04 to 1.18]	1.3 [-1.28 to 3.80]	0.99

15.8. Hook effect

The hook effect was tested using concentrations of insulin up to 25000 µIU/mL. No hook effect was observed.

16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays¹¹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

18. BIBLIOGRAPHY

1. Khan RMM, Chua ZJY, Tan JC, Yang Y, Liao Z, Zhao Y. From Pre-Diabetes to Diabetes: Diagnosis, Treatments and Translational Research. *Medicina* (Kaunas). 2019 Aug 29;55(9):546
2. Pourmotabbed G, Kitabchi AE. Hypoglycemia. *Obstet Gynecol Clin North Am.* 2001 Jun;28(2):383-400.
3. Göke B. Implications of blood glucose, insulin resistance and beta-cell function in impaired glucose tolerance. *Diabetes Res Clin Pract.* 1998 Jul;40 Suppl:S15-20.

4. Rett K, Wicklmayr M, Mehnert H. What is the clinical significance of insulin resistance? *J Cardiovasc Pharmacol.* 1992;20 Suppl 11:S22-6.
5. Behboudi-Gandevani, S., Amiri, M., Bidhendi Yarandi, R. et al. The impact of diagnostic criteria for gestational diabetes on its prevalence: a systematic review and meta-analysis. *Diabetol Metab Syndr* 11, 11 (2019).
6. John Service F. Hyperinsulinemic hypoglycemia in adults. *Ann Endocrinol (Paris).* 2004 Feb;65(1):88-95.
7. Service FJ. Hypoglycemias. *West J Med.* 1991 Apr;154(4):442-54.
8. Service FJ. Hypoglycemia. *Endocrinol Metab Clin North Am.* 1997 Dec;26(4):937-55.
9. Field JB. Hypoglycemia. Definition, clinical presentations, classification, and laboratory tests. *Endocrinol Metab Clin North Am.* 1989 Mar;18(1):27-43.
10. Basic QC Practices On-line Course; <http://www.Westgard.com>.
11. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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