DCM157-4
Ed. 03/2025

Intact PTH ELISA

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BIBLIOGRAFIA

1. Segre, G.V., Niall H.D., Habener J.F. et. al. : Metabolism of parathyroid hormone: physiological and clinical significance. *Am. J. Med.* 56: 774, 1974.
2. Mallette, L.E., Gagel, R.F.: Parathyroid Hormone and Calcitonin. In: Murray
3. J.F. (ed) *Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism.* American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 65-69, 1990.
4. Bilezikian, J.P.: Primary Hyperparathyroidism. In: Murray J.F. (ed) *Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism.* American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 109-111, 1990.
5. Stewart, A.F.: Humoral Hypercalcemia of Malignancy. In: Murray J.F. (ed) *Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism.* American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 115-118, 1990.
6. Mallette, L.E.: The parathyroid polyhormones: New concepts in the spectrum of peptide hormone action. *Endocrin. Rev.* 12:110-117, 1991.
7. Kruger, L., Rosenblum, S., Zaazra, J. and Wong, J. Intact PTH is stable in unfrozen EDTA plasma for 48 hours prior to laboratory Analysis. *Clin. Chem.* 41:6: page S47, 1995.
8. Raisz, L.G., Yajnik, C.H., Bockman, R.S., and Bower, B.B.: Comparison of commercially available parathyroid hormone immunoassay in the differential diagnosis of hypercalcemia due to primary hyperparathyroidism or malignancy. *Ann. Intern. Med.* 91:739-740, 1979.
9. Endres, D., Brickman, A., Goodman, W., Maloney, D., and Sherrard, D.: N-Terminal PTH radioimmunoassays in assessment of renal osteodystrophy. *Kidney International.* 21:132, 1982.
10. Dambacher, M.A., Fischer, J.A., Hunziker, W.H., et.al.: Distribution of circulating immunoreactive components of parathyroid hormone in normal subjects and in patients with primary and secondary hyperparathyroidism: the role of kidney and of the serum calcium concentration. *Clin. Sci.* 57:435, 1979.
11. Kao, P.C., Jiang, N.S., Klee, G.G., and Purnell, D.C.: Development and validation of a new radioimmunoassay for parathyrin (PTH). *Clin. Chem.* 28:69, 1982.

Ed. 03/2025

DCM157-4

DCM157-4
Ed. 03/2025

Intact PTH ELISA

for routine analysis

Quantitative immunoenzymatic determination of intact Parathyroid Hormone (PTH) in human serum and plasma

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO157

INTENDED USE

Colorimetric immunoenzymatic method for quantitative determination of intact Parathyroid Hormone (PTH) concentration in human serum and plasma.

Intact PTH ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

Parathyroid Hormone (PTH) is a hormone synthesized in the parathyroid gland in a 115 aminoacids form, that, after post-translational modifications, reaches the biologically active form of 84 amino acids (Intact PTH). PTH is released in the blood, where is rapidly degraded (half-life about 4 minutes); intact PTH concentration in the blood is regulated by a negative feedback mechanism that involves calcium in serum. PTH is degraded inside the parathyroid glands, in liver, kidneys and bones. A high concentration of intact PTH in the blood may indicate kidney dysfunction.

The measurement of the concentration of intact PTH in the blood is an effective way to diagnose the primary hyperparathyroidism and distinguish it from other diseases not related to parathyroid: in fact, intact PTH is elevated in 90% of cases of hyperparathyroidism.

2. PRINCIPLE OF THE METHOD

Dia.Metra Intact PTH ELISA measures the concentration of PTH in its biologically active form of 84 amino acids (intact PTH). The assay uses two polyclonal antibodies directed against specific and different regions of human PTH; one antibody is biotinylated (and allows the binding of intact PTH to the microplate, through the binding of biotin with streptavidin coated in the wells of the microplate), the other antibody is linked to horseradish peroxidase (HRP) (and allows the detection of the signal derived from the molecules of captured intact PTH). This method has been set to avoid interference from inactive fragments of PTH: therefore, only the form of 84 amino acids is measured by the assay.

From a procedural point of view, Calibrators, Controls and the serum samples are incubated with both antibodies simultaneously in the wells of the microplate. After a period of incubation (where intact PTH binds to the wells of the microplate through the bridge formed by the two antibodies), the microplate wells are washed with a wash solution to remove unbound material. The detection of the signal that quantifies the associated intact PTH is accomplished by incubation with TMB (tetramethylbenzidine), that reacts with the peroxidase and produces a colorimetric signal. This reaction is stopped by the addition of Stop Solution, a solution of sulphuric acid. The intensity of the color is directly proportional to the amount of intact PTH present in the wells; through a calibration curve, it is possible to determine the exact concentration of intact PTH bound.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 0.5 mL each)

CAL0	REF DCE002/15706-0
CAL1	REF DCE002/15707-0
CAL2	REF DCE002/15708-0
CAL3	REF DCE002/15709-0
CAL4	REF DCE002/15710-0
CAL5	REF DCE002/15711-0

2. Controls (2 vial, 0.5 mL each)

See concentration on the Certificate of Analysis

Control A	REF DCE045/15703A-0
Control B	REF DCE045/15703B-0

3. HRP Conjugate (1 vial, 7 mL)

Sheep anti PTH polyclonal antibody, conjugated to horseradish peroxidase (HRP) REF DCE002/15702-0

4. Biotin Conjugate (1 vial, 7 mL)

Sheep anti PTH polyclonal antibody, biotinylated
REF DCE019/15719-0

5. Coated Microplate (1 breakable microplate)
Microplate coated with streptavidin
REF DCE002/15703-0
6. TMB Substrate (1 vial, 12 mL)
H₂O₂-TMB (avoid any skin contact)
REF DCE004/15704-0
7. Stop Solution (1 vial, 12 mL)
Sulphuric acid 1N (avoid any skin contact)
REF DCE005/15705-0
8. 20X Conc. Wash Solution (1 vial, 25 mL)
Buffered solution
REF DCE007/15707-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 630 nm)

Notes

Store all reagents at 2÷8°C in the dark.

Open the bag of Reagent 5 (Coated Microplate) only when it is at room temperature and close it immediately after use.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the reagents should be handled in the same manner as potentially infectious material.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, controls and wash solution) contain small amounts of ProClin™ 300 as preservative. Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P280 - Wear protective gloves/ protective clothing / eye protection / face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

Stop solution

Eye irritation (Category 2A); Skin irritation

(Category 2)

Warning

Hazard statements:

H315 - Causes skin irritation

H319 - Causes a serious eye irritation

Precautionary statements:

P264 Wash exposed skin thoroughly after handling.

P280 Wear eye protection/face protection. Wear protective gloves.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P332 + P313 If skin irritation occurs: Get medical advice/attention.

P337 + P313 If eye irritation persists: Get medical advice/attention.

P362 + P364 Take off contaminated clothing and wash before reuse

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels

must be observed. Do not use any kit component beyond their expiry date.

- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.

6. PROCEDURE

6.1. Preparation of Calibrators and Controls

Calibrators and Controls are ready to use

The Calibrators have approximate concentrations of:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
pg/mL	0	15	30	100	300	1000

Exact levels are claimed on the Certificate of Analysis and labels.

Once opened, the Calibrators are stable 6 months at 2-8°C.

6.2. Preparation of the Sample

The determination of intact PTH with this kit can be performed in human serum or plasma.

The samples are stable at 2-8°C for 5 days; store at -20°C for longer time (samples are stable at -20°C for 1 month).

Avoid freezing and thawing cycles.

Do not use Sodium Azide as preservative, because it could inhibit the HRP activity.

6.3. Preparation of Wash Solution

Dilute the content of each vial of the "20X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:20 dilution ratio. The diluted wash solution is stable for 90 days at room temperature (22-28°C).

6.4. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature..
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample.

Reagent	Calibrator	Sample/Control
Calibrator C ₀ -C ₅	25 µL	
Sample/ Controls		25 µL
Biotin Conjugate	50 µL	50 µL
HRP Conjugate	50 µL	50 µL
Cover the plate and incubate at room temperature (22-28°C) for 90 minutes gently shaking (500-600 rpm). Remove the reaction mixture. Wash the wells 4 times with 300 µL of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: in case you use an automatic washer, it is advised to do 6 washing steps.		
TMB Substrate	100 µL	100 µL
Incubate 15 minutes in the dark at room temperature (22±28°C).		
Stop Solution	50 µL	50 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 630 nm within 5 minutes.		

7. QUALITY CONTROL

Each laboratory should assay controls at normal, high and low levels range of intact PTH for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied

reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

8. RESULTS

1. Calculate the mean of the absorbance (Em) for each point of the calibration curve (C₀-C₅) and of each sample.
2. Plot the mean value of absorbance (Em) of the calibrators (C₀-C₅) against concentration.
3. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).
4. Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in pg/mL.
5. For Intact PTH concentrations greater than the Cal 5 concentration (1000 pg/mL) dilute the sample 1:4 with Cal 0 (for example: to obtain 100 µL of diluted sample mix 25 µL of original sample with 75 µL of Cal 0) and test again.

9. REFERENCE VALUES

The following range can be used as initial guideline:

Expected range: 9 - 94 pg/mL

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay

Within run variation was determined by replicate (16x) the measurement of two different sera in one assay. The within assay variability is ≤ 6.3%.

10.1.2. Inter Assay

Between run variation was determined by replicate (16x) the measurement of three different sera in different lots. The between assay variability is ≤ 8.2%.

10.2. Specificity

PTH antibodies used in this method have been selected so that Dia.Metra Intact PTH ELISA kit recognises N-terminal PTH fragments, with the specificity assessed as below:

Fragment	Tested concentration	Cross-reactivity
PTH 1-84	1000 pg/mL	100%
PTH 7-84	1000 pg/mL	107%
PTH 1-13	100 ng/mL	0%
PTH 1-34	100 ng/mL	0%
PTH 39-84	1000 ng/mL	0.01%

10.3. Sensitivity

The lowest detectable concentration of intact PTH that can be distinguished from the Calibrator zero is 0.49 pg/mL at the 95% confidence limit.

10.4. Correlation

Dia.Metra Intact PTH ELISA kit was compared to a commercially available kit. 110 serum samples were analysed.

The linear regression curve was calculated:

$$Y = 1.0124 \cdot X + 2.7005$$

$$r = 0.996$$

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

BIBLIOGRAPHY

1. Segre, G.V., Niall H.D., Habener J.F. et. al. : Metabolism of parathyroid hormone: physiological and clinical significance. Am. J. Med. 56: 774,1974.
2. Mallette, L.E., Gagel, R.F.: Parathyroid Hormone and Calcitonin. In: Murray
3. J.F. (ed) Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 65-69, 1990.
4. Bilezikian, J.P.: Primary Hyperparathyroidism. In: Murray J.F. (ed) Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 109-111, 1990.
5. Stewart, A.F.: Humoral Hypercalcemia of Malignancy. In: Murray J.F. (ed) Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 115-118, 1990.
6. Mallette, L.E.: The parathyroid polyhormones: New concepts in the spectrum of peptide hormone action. Endocrin. Rev. 12:110-117, 1991.
7. Kruger, L., Rosenblum, S., Zaazra, J. and Wong,

- J. Intact PTH is stable in unfrozen EDTA plasma for 48 hours prior to laboratory Analysis. Clin. Chem. 41:6: page S47, 1995.
8. Raisz, L.G., Yajnik, C.H., Bockman, R.S., and Bower, B.B.: Comparison of commercially available parathyroid hormone immunoassay in the differential diagnosis of hypercalcemia due to primary hyperparathyroidism or malignancy. Ann. Intern. Med. 91:739-740, 1979.
 9. Endres, D., Brickman, A., Goodman, W., Maloney, D., and Sherrard, D.: N-Terminal PTH radioimmunoassays in assessment of renal osteodystrophy. Kidney International. 21:132, 1982.
 10. Dambacher, M.A., Fischer, J.A., Hunziker, W.H., et.al.: Distribution of circulating immunoreactive components of parathyroid hormone in normal subjects and in patients with primary and secondary hyperparathyroidism: the role of kidney and of the serum calcium concentration. Clin. Sci. 57:435, 1979.
 11. Kao, P.C., Jiang, N.S., Klee, G.G., and Purnell, D.C.: Development and validation of a new radioimmunoassay for parathyrin (PTH). Clin. Chem. 28:69, 1982.
 12. Endres, D.B., Villanueva, R., Sharp, C.F. Jr, Singer, F.R.: Measurement of parathyroid hormone. Endocrinol Metab. Clin. North Am. 18:611-629, 1989.
 13. Kao, P.C., van Heerden, J.A., Grant, C.S., Klee, G.G., Khosla S: Clinical performance of parathyroid hormone immunometric assays. Mayo Clin. Proc. 67:637-645, 1992
 14. Marcus, R.: Laboratory diagnosis of primary hyperparathyroidism. Endocrinol Metab. Clin. North Am. 18:647-658, 1989.

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primary hyperparathyroidism or malignancy. Ann. Intern. Med. 91:739-740, 1979.

9. Endres, D., Brickman, A., Goodman, W., Maloney, D., and Sherrard, D.: N-Terminal PTH radioimmunoassays in assessment of renal osteodystrophy. Kidney International. 21:132, 1982.
10. Dambacher, M.A., Fischer, J.A., Hunziker, W.H., et.al.: Distribution of circulating immunoreactive components of parathyroid hormone in normal subjects and in patients with primary and secondary hyperparathyroidism: the role of kidney and of the serum calcium concentration. Clin. Sci. 57:435, 1979.
11. Kao, P.C., Jiang, N.S., Klee, G.G., and Purnell, D.C.: Development and validation of a new radioimmunoassay for parathyrin (PTH). Clin. Chem. 28:69, 1982.
12. Endres, D.B., Villanueva, R., Sharp, C.F. Jr, Singer, F.R.: Measurement of parathyroid hormone. Endocrinol Metab. Clin. North Am. 18:611-629, 1989.
13. Kao, P.C., van Heerden, J.A., Grant, C.S., Klee, G.G., Khosla S: Clinical performance of parathyroid hormone immunometric assays. Mayo Clin. Proc. 67:637-645, 1992
14. Marcus, R.: Laboratory diagnosis of primary hyperparathyroidism. Endocrinol Metab. Clin. North Am. 18:647-658, 1989.

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1. Segre, G.V., Niall H.D., Habener J.F. et. al. : Metabolism of parathyroid hormone: physiological and clinical significance. Am. J. Med. 56: 774, 1974.
2. Mallette, L.E., Gagel, R.F.: Parathyroid Hormone and Calcitonin. In: Murray
3. J.F. (ed) Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 65-69, 1990.
4. Bilezikian, J.P.: Primary Hyperparathyroidism. In: Murray J.F. (ed) Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 109-111, 1990.
5. Stewart, A.F.: Humoral Hypercalcemia of Malignancy. In: Murray J.F. (ed) Primer on the Metabolic Bone Diseases and Disorders of Mineral Metabolism. American Society for Bone and Mineral Research, Kelseyville; William Byrd Press, Richmond, pp. 115-118, 1990.
6. Mallette, L.E.: The parathyroid polyhormones: New concepts in the spectrum of peptide hormone action. Endocrin. Rev. 12:110-117, 1991.
7. Kruger, L., Rosenblum, S., Zaazra, J. and Wong, J. Intact PTH is stable in unfrozen EDTA plasma for 48 hours prior to laboratory Analysis. Clin. Chem. 41:6: page S47, 1995.
8. Raisz, L.G., Yajnik, C.H., Bockman, R.S., and Bower, B.B.: Comparison of commercially available parathyroid hormone immunoassay in the differential diagnosis of hypercalcemia due to

Ed. 03/2025

DCM157-4

	DE <i>In vitro</i> Diagnostikum ES Producto sanitario para diagnóstico <i>In vitro</i> FR Dispositif médical de diagnostic <i>in vitro</i> EN <i>In vitro</i> Diagnostic Medical Device IT Dispositivo medico-diagnostico <i>in vitro</i> PT Dispositivos medicos de diagnostico <i>in vitro</i>		DE Hergestellt von ES Elaborado por FR Fabriqué par EN Manufacturer IT Produttore PT Produzido por
	DE Achtung, Begleitdokumente ES Precaución, consulte los documentos adjuntos FR Attention, veuillez consulter les documents d'accompagnement EN Caution, consult accompanying documents IT Attenzione, consultare la documentazione allegata PT Atenção, consultar os documentos de acompanhamento	 yyyy-mm	DE Herstellungs datum ES Fecha de fabricacion FR Date de fabrication EN Date of manufacture IT Data di produzione PT Data de produção
 yyyy-mm-dd	DE Verwendbar bis ES Estable hasta (usar antes de último día del mes) FR Utiliser avant (dernier jour du mois indiqué) EN Use by (last day of the month) IT Utilizzare prima del (ultimo giorno del mese) PT Utilizar (antes ultimo dia do mês)		DE Biogefährdung ES Riesco biológico FR Risque biologique EN Biological risk IT Rischio biologico PT Risco biológico
	DE Gebrauchsanweisung beachten ES Consultar las instrucciones FR Consulter le mode d'emploi EN Consult instructions for use IT Consultare le istruzioni per l'uso PT Consultar instruções para uso		DE Chargenbezeichnung ES Código de lote FR Numero de lot EN Batch code IT Codice del lotto PT Código do lote
 $\Sigma = xx$	DE Ausreichend für "n" Tests ES Contenido suficiente para "n" tests FR Contenu suffisant pour "n" tests EN Contains sufficient for "n" tests IT Contenuto sufficiente per "n" saggi PT Contém o suficiente para "n" testes		DE Inhalt ES Contenido del estuche FR Contenu du coffret EN Contents of kit IT Contenuto del kit PT Conteúdo do kit
 Max Min	DE Temperaturbereich ES Limitación de temperatura FR Limites de température de conservation EN Temperature limitation IT Limiti di temperatura PT Temperaturas limites de conservação		DE Bestellnummer ES Número de catálogo FR Références du catalogue EN Catalogue number IT Numero di Catalogo PT Número do catálogo
	DE Vor direkter Sonneneinstrahlung schützen ES Mantener alejado de la luz solar FR Tenir à l'écart de la lumière du soleil EN Keep away from sunlight IT Tenere lontano dalla luce solare PT Mantenha longe da luz solar		