DCM004-13

Ed. 01/2024

17OH PROGESTERONE ELISA

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18. BIBLIOGRAFIA

1. Honour JW. 17-Hydroxyprogesterone in children, adolescents and adults. *Ann Clin Biochem.* 2014 Jul;51(Pt 4):424-40.
2. Speiser PW, Arlt W, Auchus RJ, Baskin LS, Conway GS, Merke DP, Meyer-Bahlburg HFL, Miller WL, Murad MH, Oberfield SE, White PC. Congenital Adrenal Hyperplasia Due to Steroid 21-Hydroxylase Deficiency: An Endocrine Society Clinical Practice Guideline. *J Clin Endocrinol Metab.* 2018 Nov 1;103(11):4043-4088.
3. Bornstein SR, Allolio B, Arlt W, Barthel A, Don-Wauchope A, Hammer GD, Husebye ES, Merke DP, Murad MH, Stratakis CA, Torpy DJ. Diagnosis and Treatment of Primary Adrenal Insufficiency: An Endocrine Society Clinical Practice Guideline. *J Clin Endocrinol Metab.* 2016 Feb;101(2):364-89.
4. Sharma R, Seth A. Congenital adrenal hyperplasia: issues in diagnosis and treatment in children. *Indian J Pediatr.* 2014 Feb;81(2):178-85.
5. Legro RS, Arslanian SA, Ehrmann DA, Hoeger KM, Murad MH, Pasquali R, Welt CK; Endocrine Society. Diagnosis and treatment of polycystic ovary syndrome: an Endocrine Society clinical practice guideline. *J Clin Endocrinol Metab.* 2013 Dec;98(12):4565-92.
6. Grandone A, Marzuillo P, Luongo C, Toraldo R, Mariani M, Miraglia Del Giudice E, Perrone L. Basal levels of 17-hydroxyprogesterone can distinguish children with isolated precocious pubarche. *Pediatr Res.* 2018 Oct;84(4):533-536.
7. Rachoń D. Differential diagnosis of hyperandrogenism in women with polycystic ovary syndrome. *Exp Clin Endocrinol Diabetes.* 2012 Apr;120(4):205-9.
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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17OH PROGESTERONE ELISA

for routine analysis

Direct immunoenzymatic determination of 17OH Progesterone in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

 $\Sigma = 96$ tests

REF DKO004

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use For Laboratory Professional Use

17OH Progesterone ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of 17OH Progesterone in human serum or plasma from adult and paediatric populations. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis and monitoring of various disorders of the adrenal glands.

2. CLINICAL SIGNIFICANCE

17-OH Progesterone (17-OH P) is an important intermediate of steroid biosynthesis. Elevated levels of 17-OH P are linked to dysfunction of steroid biosynthetic pathway that leads to an androgen excess. Therefore, measurements of 17-OH P are useful for the differential diagnosis of clinical conditions related to an hyperandrogenic phenotype and linked to adrenal disfunctions in both sexes.

Congenital Adrenal Hyperplasia (CAH) is a group of autosomal recessive disorders characterised by impaired cortisol synthesis^{1,4}. The predominant cause of the condition (95% of cases) is mutations of the *CYP21A2* gene that encodes the adrenal steroid 21-hydroxylase. 21-hydroxylase is responsible for conversion of 17-OHP to 11-deoxycortisol and progesterone to deoxycorticosterone. Assessment of 17-OH P levels provides supporting information for the differential diagnosis of CAH caused by mutations of the genes involved in cortisol synthesis from *CYP21A2* gene mutations.

Measurement of 17-OH P is also considered useful in supporting diagnosis of other conditions characterised by hyperandrogenism such as polycystic ovarian syndrome (PCOS)⁵ and hirsutism in women, precocious puberty⁶ and precocious adrenarche. Determination of 17-OH P levels are used in such cases to exclude a diagnosis of CAH and non-classical CAH.

Patients affected by CAH shall be treated with glucocorticoids and mineral-corticoids. Treatment with glucocorticoids are of particular importance to avoid adrenal crisis and virilization due to excessive androgen

synthesis^{2,4,7}. Since the goal of glucocorticoid treatments is not the complete suppression of 17-OH P synthesis, it is important to ensure appropriate monitoring of treatment to evaluate the correct dosage of the drugs. This is managed by the regular assessment of 17-OH P (and androstenedione) levels as traditional indicators of the adequacy of glucocorticoids treatment in CAH⁵.

3. PRINCIPLE OF THE METHOD

The 17OH Progesterone ELISA is a competitive enzyme immunometric assay (ELISA) where 17OH Progesterone (antigen) in the sample competes with the antigenic 17OH Progesterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti 17OH Progesterone coated on the microplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in the bound fraction reacts with the Substrate (H_2O_2) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H_2SO_4) is added. The colour intensity is inversely proportional to the 17OH Progesterone concentration of in the sample.

17OH Progesterone concentration in the sample is calculated through a calibration curve.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each)

Human serum, ProClin >0.0015%, NaN_3 <0.1%

CAL0	REF DCE002/0406-0
CAL1	REF DCE002/0407-0
CAL2	REF DCE002/0408-0
CAL3	REF DCE002/0409-0
CAL4	REF DCE002/0410-0
CAL5	REF DCE002/0411-0

2. Control (2 vials, 1 mL each)

Human serum, ProClin >0.0015%, NaN_3 <0.1%

Control A	REF DCE045/0403A-0
Control B	REF DCE045/0403B-0

Control concentrations are indicated on the Certificate of Analysis

3. Conjugate (1 vial, 22 mL)
ProClin >0.0015%, BSA 0.1%
17OH Progesterone conjugated with Horseradish peroxidase (HRP) **REF DCE002/0402-0**
4. Coated Microplate (1 breakable microplate)
Microplate coated with anti 17OH Progesterone antibody **REF DCE002/0403-0**
5. TMB Substrate (1 vial, 15 mL)
H₂O₂-TMB 0.26 g/L (*avoid any skin contact*)
ProClin <0.0015% **REF DCE004-0**
6. Stop Solution (1 vial, 15 mL)
Sulphuric acid 0.15 mol/L (*avoid any skin contact*) **REF DCE005-0**
7. 10X Conc. Wash Solution (1 vial, 50 mL)
0.2M phosphate buffer pH 7.4 ProClin >0.0015% **REF DCE054-0**

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser
Precision Pipetting Devices
Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
-  All human source material used in the preparation of the reagents has been tested and found negative for antibody to HIV 1&2, HbsAg, and HCV. No test method however can offer complete assurance that HIV, HBV, HCV or other infectious agents are absent. Therefore, the Calibrators and the Controls should be handled in the same manner as potentially infectious material.
- Some reagents (calibrators, controls, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- **Classification according to Regulation (EC) No. 1272/2008 [CLP]**
Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
P321 - Specific treatment (see supplemental first aid instruction on this label).
P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.
P362+P364 - Take off contaminated clothing and wash it before reuse.

- Some reagents (calibrators and controls) contain small amounts of Sodium Azide (NaN₃) <0.1% as preservative. Sodium Azide may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through with large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate.
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2-8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Sample Storage	Duration
2 – 8 °C	96 hours
Freeze/thaw cycles	4 cycles

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

The Calibrators are ready to use and have the following concentration of 17OH Progesterone:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
ng/mL	0	0.2	0.6	2	6	16

The Controls are ready to use.

9.2. Preparation of the Conjugate

The conjugate is ready to use.

9.3. Preparation of the Wash Solution

Dilute the content of the vial "Conc. Wash Solution 10X" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

The determination of 17OH Progesterone can be performed in human serum (standard sampling tubes or tubes containing serum separating gel) or plasma (lithium heparin, sodium heparin or potassium EDTA) samples.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2-8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C ₀ -C ₅	25 µL		
Sample/Control		25 µL	
Conjugate	200 µL	200 µL	

Incubate for 1 hour at 37°C (± 0.5 °C).

Remove the contents from each well; wash the wells 6 times with 300 µL of diluted Wash Solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

In the case you use an automatic washer, after the last washing cycle remove the excess solution at the same manner.

TMB Substrate	100 µL	100 µL	100 µL
Incubate at room temperature (22±28°C) for 15 minutes in the dark.			
Stop Solution	100 µL	100 µL	100 µL

Shake the microplate gently.

Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁸.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

12. MEASURING RANGE

The assay measuring range (AMR) is 0.26 – 16 ng/mL. Any value that reads below 0.26 ng/mL should be reported as “< 0.26 ng/mL”. Any value that reads above 16 ng/mL should be reported as “> 16 ng/mL”.

13. METROLOGY AND TRACEABILITY

The 17OH Progesterone ELISA has been standardised against internal reference standards (serum matrix) which have been value assigned to a commercially available LC-MS/MS method.

14. EXPECTED VALUES

The following ranges were determined using the 17OH Progesterone ELISA and are provided for information only. The 95 % reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A “Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory”.

	n	Median (ng/mL)	Reference interval (ng/mL)
Women			
Luteal phase	123	0.45	<LoQ – 2.19
Follicular phase	122	0.32	<LoQ – 2.10
Post-menopausal	122	<LoQ	<LoQ – 0.95
Men	122	0.75	<LoQ – 1.97
Children (3 – 18 years)	124	<LoQ	<LoQ – 1.48

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

15. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

15.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A, “Protocols for Determination of Limits of Detection and Limits of Quantitation” using 6 blanks and 6 low level samples.

Sensitivity	Concentration
Limit of Blank (LoB)	0.02 ng/mL
Limit of Detection (LoD)	0.11 ng/mL
Limit of Quantitation (LoQ)	0.26 ng/mL

15.2. Trueness

Trueness has been demonstrated through method comparison of the 17OH Progesterone ELISA to a commercially available LC-MS/MS using native donor samples – refer to section 15.5.

15.3. Precision

Precision of the 17OH Progesterone ELISA was determined by performing a complex precision study.

Repeatability: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within run (Repeatability)	
			SD	CV%
1	75	0.90	0.10	11.4%
2	75	1.45	0.15	10.0%
3	75	2.10	0.20	9.3%
4	75	5.14	0.49	9.5%
5	75	8.94	0.74	8.3%
6	75	14.50	1.11	7.6%

Reproducibility: A total of 6 serum samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Data from one representative lot is shown below:

Sample	n	Mean Conc. (ng/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	75	0.90	0.12	13.0%
2	75	1.45	0.17	11.6%
3	75	2.10	0.23	11.1%
4	75	5.14	0.55	10.7%
5	75	8.94	0.82	9.2%
6	75	14.50	1.29	8.9%

15.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". For 17OH Progesterone concentration by 17OH Progesterone ELISA, the measurement procedure shows linearity for the interval from 0.2 to 17.22 ng/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$.

15.5. Method comparison

The 17OH Progesterone ELISA was compared against a commercially available LC-MS/MS method following CLSI EP-9A, "Method Comparison and Bias Estimation Using Patient Samples". A total of 120 samples, selected to represent a wide range of 17OH Progesterone concentrations was assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
120	0.91 [0.87 to 0.97]	0.21 [0.16 to 0.28]	0.99

15.6. Analytical Specificity

The specificity was assessed with the following cross-reactants.

Cross-reactant	Concentration tested (unit)	Mean %Cross reactivity
11-deoxycortisol	80 ng/mL	2.15%
Progesterone	2000 ng/L	0.47%
Pregnenolone	2000 ng/mL	0.15%
Testosterone	2000 ng/mL	0.07%
17 β -Estradiol	2000 ng/mL	0.01%
Aldosterone	2000 ng/mL	0.00%
Estriol	2000 ng/mL	0.00%
Estrone-3-Sulphate	2000 ng/mL	0.00%
Spironolactone	2000 ng/mL	0.00%
Androstenedione	2000 ng/mL	0.05%
Androsterone	2000 ng/mL	0.04%
Corticosterone	2000 ng/mL	0.09%
Cortisol	2000 ng/mL	0.10%
Cortisone	2000 ng/mL	0.07%
DHEA	2000 ng/mL	0.07%
DHEA-S	20000 ng/mL	0.01%
DHT	2000 ng/mL	0.04%
Prednisolone	3000 ng/mL	0.02%
Prednisone	2000 ng/mL	0.01%

The following substances do not interfere with a bias of $> \pm 15\%$ in the 17OH Progesterone ELISA when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	20 mg/dL
Bilirubin, unconjugated	20 mg/dL
Haemoglobin	200 mg/dL
Triglyceride	600 mg/dL

15.7. Serum-plasma study

The 17OH Progesterone ELISA matrix comparison study was performed to evaluate the difference across tube types (serum separator tubes (SST), lithium heparin plasma, sodium heparin plasma and K2 EDTA plasma) versus the control samples (red top serum, without additive) following CLSI (EP9-A) guidelines. A total of 27 samples (23 native, 4 spiked) to cover the assay range were evaluated. Passing-

Bablok regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (ng/mL) [95% CI]	Correlation coefficient (r)
SST	1.02 [0,91 to 1,13]	-0.03 [-0,21 to 0,11]	0.99
Lithium Heparin	1.00 [0,87 to 1,04]	0.00 [-0,07 to 0,08]	0.99
Sodium Heparin	1.00 [0,94 to 1,10]	-0.02 [-0,14 to 0,03]	0.99
EDTA	1.04 [0,90 to 1,12]	-0.03 [-0,19 to 0,09]	0.99

16. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

17. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

18. BIBLIOGRAPHY

1. Honour JW. 17-Hydroxyprogesterone in children, adolescents and adults. *Ann Clin Biochem*. 2014 Jul;51(Pt 4):424-40.
2. Speiser PW, Arlt W, Auchus RJ, Baskin LS, Conway GS, Merke DP, Meyer-Bahlburg HFL, Miller WL, Murad MH, Oberfield SE, White PC. Congenital Adrenal Hyperplasia Due to Steroid 21-Hydroxylase Deficiency: An Endocrine Society Clinical Practice Guideline. *J Clin Endocrinol Metab*. 2018 Nov 1;103(11):4043-4088.
3. Bornstein SR, Allolio B, Arlt W, Barthel A, Don-Wauchope A, Hammer GD, Husebye ES, Merke DP, Murad MH, Stratakis CA, Torpy DJ. Diagnosis and Treatment of Primary Adrenal Insufficiency: An

Endocrine Society Clinical Practice Guideline. *J Clin Endocrinol Metab*. 2016 Feb;101(2):364-89.

4. Sharma R, Seth A. Congenital adrenal hyperplasia: issues in diagnosis and treatment in children. *Indian J Pediatr*. 2014 Feb;81(2):178-85.
5. Legro RS, Arslanian SA, Ehrmann DA, Hoeger KM, Murad MH, Pasquali R, Welt CK; Endocrine Society. Diagnosis and treatment of polycystic ovary syndrome: an Endocrine Society clinical practice guideline. *J Clin Endocrinol Metab*. 2013 Dec;98(12):4565-92.
6. Grandone A, Marzuillo P, Luongo C, Toraldo R, Mariani M, Miraglia Del Giudice E, Perrone L. Basal levels of 17-hydroxyprogesterone can distinguish children with isolated precocious pubarche. *Pediatr Res*. 2018 Oct;84(4):533-536.
7. Rachoń D. Differential diagnosis of hyperandrogenism in women with polycystic ovary syndrome. *Exp Clin Endocrinol Diabetes*. 2012 Apr;120(4):205-9.
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

19. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

20. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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1. Honour JW. 17-Hydroxyprogesterone in children, adolescents and adults. Ann Clin Biochem. 2014 Jul;51(Pt 4):424-40.
2. Speiser PW, Arlt W, Auchus RJ, Baskin LS, Conway GS, Merke DP, Meyer-Bahlburg HFL, Miller WL, Murad MH, Oberfield SE, White PC. Congenital Adrenal Hyperplasia Due to Steroid 21-Hydroxylase Deficiency: An Endocrine Society Clinical Practice Guideline. J Clin Endocrinol Metab. 2018 Nov 1;103(11):4043-4088.
3. Bornstein SR, Allolio B, Arlt W, Barthel A, Don-Wauchope A, Hammer GD, Husebye ES, Merke DP, Murad MH, Stratakis CA, Torpy DJ. Diagnosis and Treatment of Primary Adrenal Insufficiency: An Endocrine Society Clinical Practice Guideline. J Clin Endocrinol Metab. 2016 Feb;101(2):364-89.
4. Sharma R, Seth A. Congenital adrenal hyperplasia: issues in diagnosis and treatment in children. Indian J Pediatr. 2014 Feb;81(2):178-85.
5. Legro RS, Arslanian SA, Ehrmann DA, Hoeger KM, Murad MH, Pasquali R, Welt CK; Endocrine

Society. Diagnosis and treatment of polycystic ovary syndrome: an Endocrine Society clinical practice guideline. *J Clin Endocrinol Metab.* 2013 Dec;98(12):4565-92.

6. Grandone A, Marzuillo P, Luongo C, Toraldo R, Mariani M, Miraglia Del Giudice E, Perrone L. Basal levels of 17-hydroxyprogesterone can distinguish children with isolated precocious pubarche. *Pediatr Res.* 2018 Oct;84(4):533-536.
7. Rachoń D. Differential diagnosis of hyperandrogenism in women with polycystic ovary syndrome. *Exp Clin Endocrinol Diabetes.* 2012 Apr;120(4):205-9.
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

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