DCM020-13

Ed. 03/2025

CORTISOL SALIVA ELISA

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CORTISOL SALIVA ELISA

for routine analysis

Direct immunoenzymatic determination of Cortisol in saliva.

IVD

LOT

See external label

2°C  8°C $\Sigma = 96$ tests

REF DKO020

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

Cortisol Saliva ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of cortisol in human saliva from an adult population.

2. CLINICAL SIGNIFICANCE

Cortisol is a steroid hormone released from the adrenal cortex in response to a hormone called ACTH (produced by the pituitary gland), it is involved in the response to stress; it increases blood pressure, blood sugar levels, may cause infertility in women, and suppresses the immune system.

Cortisol acts through specific intracellular receptors and has effects in numerous physiologic systems, including immune function, glucose-counter regulation, vascular tone, substrate utilization and bone metabolism. Cortisol is excreted primarily in urine in an unbound (free) form.

The majority of cortisol in saliva is unbound and enters the saliva via intracellular mechanisms. Salivary cortisol levels are unaffected by salivary flow rate or salivary enzymes. There is a high correlation between serum and saliva cortisol levels.

These normal endogenous functions are the basis for the physiological consequences of chronic stress - prolonged cortisol secretion causes muscle wastage, hyperglycaemia, and suppresses immune / inflammatory responses. The same consequences arise from long-term use of glucocorticoid drugs.

3. PRINCIPLE OF THE METHOD

The Cortisol Saliva ELISA is a competitive enzyme immunometric assay (ELISA) where cortisol (antigen) in the sample competes with the antigenic cortisol conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti- cortisol coated on the microplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in

the bound fraction reacts with the Substrate (H_2O_2) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H_2SO_4) is added. The colour intensity is inversely proportional to the cortisol concentration in the sample.

Cortisol concentration in the sample is calculated through a calibration curve.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (7 vials, 1 mL each)

ProClin >0.0015, BSA 0.1%

CAL0	REF DCE002/2006-0
CAL1	REF DCE002/2007-0
CAL2	REF DCE002/2008-0
CAL3	REF DCE002/2009-0
CAL4	REF DCE002/2010-0
CAL5	REF DCE002/2011-0
CAL6	REF DCE002/2012-0

2. Control (2 vials, 1 mL each)

Concentration of Control is indicated on the Certificate of Analysis. ProClin >0.0015, BSA 0.1%

Control A	REF DCE045/2003A-0
Control B	REF DCE045/2003B-0

3. Incubation Buffer (1 vial, 30 mL)

Phosphate buffer 50 mM pH 7.4, BSA 1 g/L.

Proclin >0.0015% REF DCE001-0

4. Conjugate (1 vial, 1 mL)

Cortisol conjugated with horseradish peroxidase (HRP). ProClin >0.0015, BSA 0.1% REF DCE002/2002-0

5. Coated Microplate (1 microplate breakable)

Anti- Cortisol antibodies adsorbed on microplate

REF DCE002/2003-0

6. TMB Substrate (1 vial, 15 mL)

 H_2O_2 -TMB (0.26 g/L) (avoid any skin contact)

ProClin <0.0015% REF DCE004-0

7. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

8. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M pH 7.4. ProClin >0.0015%

REF DCE054-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

Sarstedt Salivette

REF 51.1534.500

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, controls, incubation buffer, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 – 8°C in their original container. Any exceptions are clearly indicated.

- Allow all kit components and specimens to reach room temperature (22 – 28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the calibrators are stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2 – 8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION

The assay should be performed using saliva samples.

It is recommended to collect saliva samples with a centrifuge glass tube and a plastic straw, with the "Salivette" (Sarstedt, Ref. 511534500). Other commercially available sample collector devices have not been tested.

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer.

The calibrators are ready for use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆
ng/mL	0	0.5	1	5	10	20	100

The controls are ready to use; the concentration is printed on the label.

9.2. Preparation of the Conjugate

Prepare immediately before use.

Add 10 µL of Conjugate (reagent 3) to 1 mL of Incubation Buffer (reagent 2).

Mix gently. Stable for 3 hours at room temperature (22-28°C).

9.3. Preparation of the Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

9.4.1. Saliva Processing Instructions with Salivette Sardstedt

1. Remove the swab from the suspended insert of the Salivette
2. Gently chewing the swab for 1 minute produces a sufficient quantity of saliva.
3. Replace the swab into the Salivette and firmly close the tube using the stopper.
4. Centrifuge the Salivette for 2 minutes at 1000g (rcf) for saliva generation.
5. Remove the insert complete with the swab from the centrifuge vessel and discard. The clear saliva is now ready for analysis (at least 1 mL of saliva should be recovered with this method).

9.4.2. Method and Limitations

Collect saliva samples at the times indicated.

If no specific instructions have been given, saliva samples may be collected at any time; however the following should be noted:

- a) If saliva is collected in the morning ensure that this is carried out prior to brushing teeth
- b) During the day allow 1 hour after a meal, oral intake of pharmaceutical drugs or tooth cleaning.
- c) It is very important that a good clear sample is received – i.e. no contamination with food, lipstick, blood (bleeding gums) or other such extraneous materials.

9.5. Procedure

- **Allow all reagents to reach room temperature (22 – 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₆), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C ₀ -C ₆	25 µL		
Sample/Control		25 µL	
Conjugate	200 µL	200 µL	

Incubate 1h at +37°C ±0.5°C.

Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 6 times.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at 22 – 28°C for 15 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake gently the microplate.

Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁸.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required.** Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

Conversion of units

To convert results to SI units:
 $\text{nmol/L} = \text{ng/mL} \times 2.76$

To convert results to mass units:
 $\text{ng/mL} = \text{nmol/L} \times 0.36$

12. EXPECTED VALUES

The following values can be used as preliminary guideline until each laboratory established its own normal range.

	Reference value (ng/mL)
AM	3 – 10 ng/mL
PM	0.6 – 2.5 ng/mL

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

13. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

13.1. Precision

13.1.1. Intra Assay Variation

Within run variation was determined by replicate (20x) the measurement of three different saliva samples in one assay. The within assay variability is $\leq 10\%$.

13.1.2. Inter Assay Variation

Between run variation was determined by replicate (10x) the measurement of three different saliva samples in different lots of kits. The between assay variability is $\leq 8.3\%$.

13.2. Accuracy

The recovery of 6.25 - 12.5 - 25 - 50 ng/mL of Cortisol added to a saliva sample gave an average value ($\pm$ SD) of $95.42\% \pm 9.11\%$ with reference to the original concentrations.

13.3. Sensitivity

The lowest detectable concentration of Cortisol that can be distinguished from the Calibrator 0 is 0.12 ng/mL at the 95% confidence limit.

13.4. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

Cortisol	100 %
Prednisolone	46.2 %
11-Deoxycortisol	4 %
Cortisone	3.69 %
Prednisone	3.10 %
11 α OH Progesterone	1 %
Progesterone	< 0.1 %
Aldosterone	< 0.1 %
Pregnenolone	< 0.1 %
17b Estradiol	< 0.1 %
Estrone 3-solfato	< 0.1 %
Estriol	< 0.1 %
Testosterone	< 0.1 %
Spironolactone	< 0.1 %
DHEA	< 0.1 %
DHEA-S	< 0.1 %
Androstenedione	< 0.1 %
Androsterone	< 0.1 %
DHT	< 0.1 %
Danazol	< 0.1 %
Cholesterol	< 0.1 %
Dexamethasone	< 0.1 %

13.5. Method comparison

The Cortisol Saliva ELISA kit was compared to a previous version of the DiaMetra Cortisol Saliva ELISA kit. 35 saliva samples were analysed.

The linear regression curve was calculated:

n	Slope	Intercept (ng/mL)	R ²
35	0.94	-0.11	0.799

14. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

15. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

16. BIBLIOGRAPHY

1. Foster, L. B. and Dunn, R.T. Clin. Chem: 20/3, 365 (1974)
2. De Lacerda L, et al J. Clin. Endocr. and Metab: 36, 227 (1973)
3. Roller, E., et al Clin chim Acta 66 319 (1976)
4. Kobayashi, Y., et al Steroids, 32 no. 1 (1978)
5. Arakawa, H., et al Anal. Biochem. 97 248 (1979)
6. Kirschbaum C. et al, Psychoneuroendocrinology, 19: 313-333
7. Cristina Mihaela Ghiciuc C.M et al., Neuroendocrinol Lett 2011; 32(4):475-480
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. Clin Chem, 34, 1988, pp 27-33

17. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

18. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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