DCM053-13

Ed. 03/2024

ALDOSTERONE ELISA

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19. BIBLIOGRAFIA

1. Williams Textbook of Endocrinology, 10th Edition, Editors Larsen, Kronenberg, Melmed, Polonski, Willson
2. Duncan J. Campbell, Juerg Nussberger, Michael Stowasser, A.H. Jan Danser, Alberto Morganti, Erik Frandsen, and Joël Menard, Activity Assays and Immunoassays for Plasma Renin and Prorenin: Information Provided and Precautions Necessary for Accurate Measurement. Clinical Chemistry 55:5, 867–877 (2009).
3. John W. Funder et al. Case Detection, Diagnosis and Treatment of Patients with Primary

Aldosteronism: An Endocrine Society Clinical Practice Guideline. J. Clin Endocrinol. Metab. September 2008 93 (9) 2266-3281.

4. Caroline Schirpenbach, Lysann, Seiler, Christiane Maser-Gluth, Felix Beuschlein, Martin Reincke and Martin Bidlingmaier Automated Chemiluminescence-Immunoassay for Aldosterone during Dynamic Testing: Comparison to Radioimmunoassays with and without Extraction Steps. Clinical Chemistry. 52-9 1749-1755 (2006).
5. Basic QC Practices On-line Course; <http://www.Westgard.com>.
6. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. Clin Chem, 34, 1988, pp 27-33

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ALDOSTERONE ELISA

for routine analysis

Immunoenzymatic determination of Aldosterone in human serum, human plasma or urine.

IVD

LOT

See external label

2°C 8°C

Σ = 96 tests

REF DKO053

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

Aldosterone ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of Aldosterone in human serum, plasma or urine from an adult population. Results are to be used in conjunction with other clinical and laboratory data as an aid in the diagnosis of hypertension related syndromes.

2. CLINICAL SIGNIFICANCE

Aldosterone is a steroid hormone produced by the zona glomerulosa of the adrenal cortex in the adrenal gland. It plays a primary role in the regulation of sodium and potassium balance in the blood¹. Aldosterone secretion is stimulated primarily via the renin angiotensin-aldosterone system (RAAS). The production of renin is stimulated when there is a reduction in renal perfusion pressure, reduction in plasma volume and a negative sodium balance. Renin cleaves Angiotensinogen into Angiotensin I which ultimately leads to the production of Angiotensin II via angiotensin converting enzyme². Angiotensin II acts on the vascular system causing vasoconstriction as well as stimulating the adrenal cortex to secrete aldosterone. Aldosterone acts on mineralocorticoid receptors within the principal cells of the distal tubule to increase reabsorption of sodium and chloride and decrease reabsorption of potassium and hydrogen. Aldosterone measurements are used in the diagnosis and treatment of primary aldosteronism^{3,4} (a disorder caused by excessive secretion of aldosterone by the adrenal gland), hypertension caused by primary aldosteronism, selective hyperaldosteronism, edematous states and other conditions of electrolyte balance.

3. PRINCIPLE OF THE METHOD

The Aldosterone ELISA is a competitive enzyme immunometric assay (ELISA) where aldosterone (antigen) in the sample competes with the antigenic aldosterone conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti- aldosterone coated on the microplate (solid phase).

After the incubation, the bound/free separation is performed by a simple solid phase washing. Then, the enzyme HRP in

the bound fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blue colour that changes into yellow when the Stop Solution (H₂SO₄) is added. The colour intensity is inversely proportional to the aldosterone concentration in the sample.

Aldosterone concentration in the sample is calculated through a calibration curve.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each)

Contains ProClin > 0.0015% and BSA

CAL0	REF DCE002/5306-0
CAL1	REF DCE002/5307-0
CAL2	REF DCE002/5308-0
CAL3	REF DCE002/5309-0
CAL4	REF DCE002/5310-0
CAL5	REF DCE002/5311-0

2. Controls (2 vial, 1 mL each)

Contains ProClin > 0.0015% and BSA

Control A	REF DCE045/5303A-0
Control B	REF DCE045/5303B-0

Control concentration is indicated on the Certificate of Analysis

3. Conjugate (1 vial, 15 mL)

Aldosterone conjugated with horseradish peroxidase (HRP).

Contains ProClin < 0.0015% and BSA

REF DCE002/5302-0

4. Coated Microplate (1 breakable microplate)

Microplate coated with anti- aldosterone antibody

REF DCE002/5303-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB (0.26 g/L) (avoid any skin contact)

REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)

REF DCE005-0

7. 10X Conc. Wash Solution (1 vial, 50 mL)

Contains ProClin > 0.0015%

0.2M Phosphate buffer pH 7,4

REF DCE054-0

4.2. Materials required but not provided

Distilled water

PBS buffer pH 7.4

TRIS buffer 1M pH 8.5

4.3. Auxiliary materials and instrumentation

Polypropylene tubes (for urine samples)

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents contain small amounts of ProClin™ 300 as preservative. Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the kit is stable at 2 – 8°C for 6 months*.
- The diluted wash solution is stable for 30 days at 2-8°C.

* in use stability data supports reagent stability when used two times within this period.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum (standard sampling tubes or tubes containing serum separating gel), plasma (lithium heparin, sodium heparin or potassium EDTA) or urine samples.

Sample Storage	Duration
2 – 8 °C (serum, plasma & urine)	96 hours
Freeze/thaw cycles (serum, plasma & urine)	3 cycles
-20°C (plasma & urine)	12 months

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

The calibrators are ready to use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
pg/mL	0	40	80	300	800	1260

Mix gently on a roller mixer for at least 5 minutes before using.

The Controls are ready to use.

9.2. Preparation of the Conjugate

Ready to use. Mix gently for 5 minutes on a roller mixer.

9.3. Preparation of the Wash Solution

Dilute the content of the vial "10X Conc. Wash Solution" with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

The determination of Aldosterone can be performed in human serum (standard sampling tubes or tubes containing serum separating gel), plasma (lithium heparin, sodium heparin or potassium EDTA) or urine samples. Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.4.1. Serum or plasma samples

For serum or plasma samples proceed to section 9.5

9.4.2. Urine samples

To obtain urine samples use the following procedure:

- Pipette 25 µL of each urine sample into a polypropylene tube

- Add 25 µL of 1M HCl to every tube and mix thoroughly
- Cap and incubate at 2 – 8°C for 24h in the dark
- Add 25 µL of TRIS pH 8.5 to every tube and mix thoroughly
- Add 400 µL of PBS buffer
- Use 50 µL of this solution to perform the assay as stated in section 9.5

Note: this pretreatment leads to a 1:19 dilution. The dilution factor 19 has to be taken into account for the calculation of the final concentration of the urine sample. Please refer to section 11

Store the sample at -20°C if the determination is not performed on the same day of the sample collection.

9.5. Procedure

- Allow all reagents to reach room temperature (22 – 28°C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Samples/ Control	Blank
Calibrator C ₀ -C ₅	50 µL		
Samples/ Controls		50 µL	

Incubate at 22 – 28°C for 30 minutes.

Conjugate	100 µL	100 µL	
-----------	--------	--------	--

Incubate at 22 – 28°C for 1 hour.

Remove the contents from each well; wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at 22 – 28°C for 20 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
Shake the microplate gently. Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.			

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The controls provided in the kit should be tested as unknowns and are intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁵.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in pg/mL.

Conversion of units:

pg/mL / 10 = ng/dL

pg/mL x 2.774 = pmol/L

For urine samples the concentration in pg/mL read from the calibration curve has to be multiplied to the **dilution factor 19**.

To convert the Aldosterone concentration in "mcg Aldosterone/24h" (µg Aldosterone/24h), calculate as above and correct for total volume of urine collected in 24 hours:

$$\frac{\text{pg/mL} \times \text{Vol (mL) urine in 24hr}}{1,000,000} = \text{mcg Aldosterone/24hr}$$

12. MEASURING RANGE

The assay measuring range (AMR) for plasma samples is 28.0 – 1260 pg/mL.

Any value that reads below 28.0 pg/mL should be reported as "< 28.0 pg/mL". Any value that reads above 1260 pg/mL should be reported as "> 1260 pg/mL".

The assay measuring range (AMR) for urine samples is 38.3 – 1260 pg/mL.

Any value that reads below 38.3 pg/mL should be reported as "< 38.3 pg/mL". Any value that reads above 1260 pg/mL should be reported as "> 1260 pg/mL".

13. METROLOGY AND TRACEABILITY

The Aldosterone ELISA has been standardised against in-house reference standards which have been value assigned by the Isotope dilution mass spectrometry (ID/GC/MS) from DGKC.

14. Uncertainty of Measurement

The uncertainty associated with the measurement of the assay calibrators (C₀ – C₅) was determined following CLSI EP29-A "Expression of Measurement Uncertainty in Laboratory Medicine; Approved Guideline". The uncertainty for each calibrator level was assessed over a minimum of 5 days (1 run per day) by 3 different operators. Result from one representative lot:

Calibrator	Mean conc. (pg/mL)	95% CI (pg/mL)		Calibrator uncertainty Uc (pg/mL)
CAL0	<LoB	N/A		N/A
CAL1	39.6	28.2	51.0	4.6
CAL2	84.3	64.1	104.5	6.0
CAL3	299.0	212.6	385.3	23.0
CAL4	792.2	737.9	846.4	43.9
CAL5	1271.4	1103.1	1439.7	37.0

15. EXPECTED VALUES

The following ranges were determined using the Aldosterone ELISA and are provided for information only. The 95% reference interval for apparently healthy adults were calculated by a non-parametric method following guidance from CLSI C28-A3 "Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory".

	n	Median	Normal range
Plasma samples (pg/mL)			
Upright	60	37.4	<28.0 – 221.8
Supine	60	45.5	<28.0 – 141.5
Urine samples (mcg/24 hr)			
	120	10.5	3.4 – 23.2

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

16. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

16.1. Detection Capability

The limit of blank (LoB), limit of detection (LoD) and limit of quantitation (LoQ) were determined with guidance from CLSI EP17-A2, "Protocols for Determination of Limits of Detection and Limits of Quantitation" using 6 blanks and 6 low level samples.

Sensitivity	Concentration	
	Plasma	Urine
Limit of Blank (LoB)	0.9	4.2
Limit of Detection (LoD)	10.1	16.1
Limit of Quantitation (LoQ)	28.0	38.3

16.2. Trueness

Trueness has been demonstrated through comparison of the Aldosterone ELISA to the isotope dilution mass spectrometry (ID/GC/MS) method from DGKC with a panel of native plasma samples.

	pg/mL
Linear regression	
Slope	1.01
Intercept	25.80
Correlation Coefficient, r	0.99

16.3. Precision

Precision of the Aldosterone ELISA was determined by performing a complex precision study.

Plasma samples

Repeatability: A total of 6 EDTA plasma samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (pg/mL)	Within run (Repeatability)	
			SD	CV%
1	75	87.2	8.8	10.0%
2	75	233.3	19.6	8.4%
3	75	351.8	23.1	6.6%
4	75	564.1	42.7	7.6%
5	75	1028.4	59.0	5.7%
6	75	1136.6	51.3	4.5%

Reproducibility: A total of 6 EDTA plasma samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (pg/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	86.9	8.5	9.8%
2	150	235.5	24.7	10.5%
3	150	353.7	37.1	10.5%
4	150	570.4	50.4	8.8%
5	150	1035.6	71.6	6.9%
6	150	1136.8	64.7	5.7%

Urine samples

Repeatability: A total of 6 urine samples were assayed in 5 replicates, once a day for 5 days by 3 operators.

Data from one representative lot is shown below:

Sample	n	Mean Conc. (pg/mL)	Within run (Repeatability)	
			SD	CV%
1	75	122.9	12.7	10.3%
2	75	218.9	16.8	7.7%
3	75	323.4	20.1	6.2%
4	75	563.1	27.6	4.9%
5	75	1027.3	54.0	5.3%
6	75	1187.9	42.1	3.5%

Reproducibility: A total of 6 urine samples were assayed in 5 replicates, once a day for 5 days by 3 operators. Results for the combined data from two lots is shown below:

Sample	n	Mean Conc. (pg/mL)	Within Laboratory (Reproducibility)	
			SD	CV%
1	150	123.5	18.9	15.3%
2	150	219.3	24.4	11.1%
3	150	322.0	31.0	9.6%
4	150	555.8	48.5	8.7%
5	150	1003.4	93.8	9.4%
6	150	1170.3	76.9	6.6%

16.4. Linearity

Linearity was evaluated based on CLSI EP-06, "Evaluation of the Linearity of Quantitative Measurement Procedures". A serial dilution of high and low plasma / urine samples was prepared to create a panel from above the upper limit of linearity (ULL) to below the lower limit of quantitation (LLOQ).

For aldosterone concentration measured in plasma samples by Aldosterone ELISA, the measurement procedure shows linearity for the interval from 15.4 to 1394 pg/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$ or $\leq \pm 28$ pg/mL below LoQ.

For aldosterone concentration measured in urine samples by Aldosterone ELISA, the measurement procedure shows linearity for the interval from 20.6 to 1635.8 pg/mL within the allowable deviation of linearity (ADL) of $\pm 15\%$ or $\leq \pm 38.3$ pg/mL below LoQ.

16.5. Method comparison

The Aldosterone ELISA was compared against a commercially available quantitative manual ELISA method, following CLSI EP-9C, "Method Comparison and Bias Estimation Using Patient Samples". A total of 106 plasma samples selected to represent a wide range of aldosterone concentrations, were assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
106	1.1 [1.0 to 1.2]	-3.9 [-13.2 to 5.2]	0.98

A total of 114 urine samples, selected to represent a wide range of aldosterone concentrations, were assayed by each method. Passing-Bablok regression analysis was performed on the comparative data:

n	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
114	0.95 [0.9 to 1.0]	18.9 [10.7 to 27.8]	1.0

16.6. Analytical Specificity

The specificity was assessed with the following cross-reactants – assessed separately in both plasma and urine.

Cross-reactant	Concentration tested (unit)	Mean Cross reactivity	
		Plasma	Urine
Corticosterone	0.1 mg/mL	0.1%	0.1%
Androsterone	5 mg/mL	0.0%	0.0%
Cortisone	0.1 mg/mL	0.0%	0.1%
Dihydrotestosterone (DHT)	0.5 mg/mL	0.0%	0.0%
Estradiol	10 mg/mL	0.0%	0.0%
Estrone	10 mg/mL	0.0%	0.0%
Estriol	10 mg/mL	0.0%	0.0%
Testosterone	0.2 mg/mL	0.0%	0.0%
DHEA-S	10 mg/mL	0.0%	0.0%

Plasma samples

The following substances do not interfere with a bias of $> \pm 15\%$ in the Aldosterone ELISA assay when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Bilirubin, conjugated	15 mg/dL
Bilirubin, unconjugated	10 mg/dL
Haemoglobin	200 mg/dL
Total Protein	8.3 g/dL
Triglyceride	500 mg/dL

Urine samples

The following substances do not interfere with a bias of $> \pm 15\%$ in the Aldosterone ELISA assay when the concentrations are below the stated threshold presented in the following table.

Potentially Interfering Reagent	Threshold Concentration
Protein	80 mg/dL
Sodium Chloride	0.5 M
Urea	700 mM
Creatinine	10 mM
Glucose	4 mM
Boric Acid	20 g/L

16.7. Serum-plasma study

The Aldosterone ELISA matrix comparison study was performed to evaluate the difference across tube types (serum, serum separator tubes (SST), lithium heparin plasma and sodium heparin plasma) versus the control samples (K2 EDTA) following CLSI EP35-Ed1 guidelines. A total of 22 samples (18 native, 4 spiked or diluted) were evaluated. Passing-Bablok regression analysis was performed on the comparative data:

Sample type	Slope [95% CI]	Intercept (pg/mL) [95% CI]	Correlation coefficient (r)
Serum	0.98 [0.9 to 1.0]	-12.25 [-30.1 to 5.6]	1.0
SST	1.03 [1.0 to 1.1]	-11.6 [-37 to 13.8]	1.0
Lithium Heparin	0.98 [0.9 to 1.0]	-3.4 [-28.5 to 21.7]	1.0
Sodium Heparin	0.94 [0.9 to 1.0]	-6.4 [-22.4 to 9.5]	1.0

17. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁶. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

18. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

19. BIBLIOGRAPHY

1. Williams Textbook of Endocrinology, 10th Edition, Editors Larsen, Kronenberg, Melmed, Polonski, Willson
2. Duncan J. Campbell, Juerg Nussberger, Michael Stowasser, A.H. Jan Danser, Alberto Morganti, Erik Frandsen, and Joël Menard, Activity Assays and Immunoassays for Plasma Renin and Prorenin: Information Provided and Precautions Necessary for Accurate Measurement. *Clinical Chemistry* 55:5, 867–877 (2009).
3. John W. Funder et al. Case Detection, Diagnosis and Treatment of Patients with Primary Aldosteronism: An Endocrine Society Clinical Practice Guideline. *J. Clin Endocrinol. Metab.* September 2008 93 (9) 2266-3281.
4. Caroline Schirpenbach, Lysann, Seiler, Christiane Maser-Gluth, Felix Beuschlein, Martin Reincke and Martin Bidlingmaier Automated Chemiluminescence-Immunoassay for Aldosterone during Dynamic Testing: Comparison to Radioimmunoassays with and without Extraction Steps. *Clinical Chemistry*. 52:9 1749-1755 (2006).
5. Basic QC Practices On-line Course; <http://www.Westgard.com>.
6. Boscatto, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

20. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

21. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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1. Williams Textbook of Endocrinology , 10th Edition, Editors Larsen, Kronenberg, Melmed, Polonski, Willson
2. Duncan J. Campbell, Juerg Nussberger, Michael Stowasser, A.H. Jan Danser, Alberto Morganti, Erik Frandsen, and Joël Menard, Activity Assays and Immunoassays for Plasma Renin and Prorenin: Information Provided and Precautions Necessary for Accurate Measurement. Clinical Chemistry 55:5, 867–877 (2009).

3. John W. Funder et al. Case Detection, Diagnosis and Treatment of Patients with Primary Aldosteronism: An Endocrine Society Clinical Practice Guideline. *J. Clin Endocrinol. Metab.* September 2008 93 (9) 2266-3281.
4. Caroline Schirpenbach, Lysann, Seiler, Christiane Maser-Gluth, Felix Beuschlein, Martin Reincke and Martin Bidlingmaier Automated Chemiluminescence-Immunoassay for Aldosterone during Dynamic Testing: Comparison to Radioimmunoassays with and without Extraction Steps. *Clinical Chemistry*. 52-9 1749-1755 (2006).
5. Basic QC Practices On-line Course; <http://www.Westgard.com>.
6. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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