DCM009-13

Ed. 03/2024

LH ELISA

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16. BIBLIOGRAFIA

1. Kosasa T.S., "Measurement of Human Luteinizing Hormone." *Journal of Reproductive Medicine*, 26 (1981) pg. 201-6.
2. Danzer H., Braunstein G.D., et al., "Maternal Serum Human Chorionic Gonadotropic Concentrations and Fetal Sex Predictions." *Fertility and Sterility*, 34 (1980) pg. 336-40.
3. Braunstein G.D., et al., "Serum Human Luteinizing Hormone Levels through Normal Pregnancy ", *American Journal of Obstetrics and Gynecology*, 126 (1976) pg. 678-81.
4. Goldstein D.P., and Kosasa T.S., "The Subunit Radioimmunoassay for LH Clinical Application." *Gynecology*, 6 (1975) pg. 145-84.
5. Batzer F., "Hormonal Evaluation of Early Pregnancy", *Fertility and Sterility*, 34 (1980) pg. 1-12.
6. Braunstein, G.D., et al., "First-Trimester Luteinizing Hormone Measurements as an Aid to the Diagnosis of Early Pregnancy Disorders", *American Journal of Obstetrics and Gynecology*, 131 (1978) pg. 25-32.
7. Lenton E., Neal L. and Sulaiman R., "Plasma Concentrations of Human Gonadotropin from the time of Implantation until the Second Week of Pregnancy ", *Fertility and Sterility*, 37 (1982) pg. 773-78.
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.

9. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

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LH ELISA

for routine analysis

Direct immunoenzymatic determination of luteinising hormone (LH) in human serum or plasma

IVD

LOT

See external label

2°C  8°C $\Sigma = 96$ tests

REF DKO009

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use
For Laboratory Professional Use

LH ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of luteinising hormone (LH) in human serum or plasma from an adult population.

2. CLINICAL SIGNIFICANCE


Luteinising hormone (LH) is a glycoprotein consisting of two subunits with a molecular mass of 30,000 Daltons. The α -subunit is similar to other pituitary hormones [follicle stimulating hormone (FSH), thyroid stimulating hormone (TSH) and chorionic gonadotropin (HCG)] while the β -subunit is unique. The β -subunit confers the biological activity to the molecule.

The α -subunit consists of 89 amino acid residues while the β -subunit contains 129 amino acids. The carbohydrate content is between 15% and 30%. The clinical usefulness of the measurement of luteinising hormone (LH) in ascertaining the homeostasis of fertility regulation via the hypothalamic - pituitary - gonadal axis has been well established^{1,2}. In addition, the advent of *in vitro* fertilisation (IVF) technology to overcome infertility associated problems has provided the impetus for rapid improvement in LH assay methodology from the technically demanding bioassay³ to the procedurally simple and rapid immunoenzymometric assays.

3. PRINCIPLE OF THE METHOD

In this method, the LH calibrators, the patient specimens and/or the controls (containing the native antigen) are first added to streptavidin coated wells. Then monoclonal biotinylated and enzyme labelled antibodies are added and the reactants mixed: these antibodies have high affinity and specificity and are directed against distinct and different epitopes of LH. A reaction between the various LH antibodies and native LH occurs in the microwells without competition or steric hindrance forming a soluble sandwich complex.

The interaction is illustrated in the following equation:

BtnAb(m) = biotinylated monoclonal antibody (excess quantity)

AgLH = native antigen (variable quantity)


EnzAb = enzyme labeled antibody (excess quantity)

EnzAb-AgLH-BtnAb(m) = antigen-antibodies sandwich complex

k_a = rate constant of association

k_{-a} = rate constant of dissociation

Simultaneously, the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

Streptavidin C.W. = streptavidin immobilised on well.

Immobilised complex = antibodies-antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by a washing step. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well quantitated by reaction with a suitable substrate to produce colour. By utilising several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each)

ProClin >0.0015%, BSA 1%

CAL0	REF DCE002/0906-0
CAL1	REF DCE002/0907-0
CAL2	REF DCE002/0908-0
CAL3	REF DCE002/0909-0
CAL4	REF DCE002/0910-0
CAL5	REF DCE002/0911-0

2. Control (1 vial, 1 mL)

Concentration of Control is indicated on the Certificate of Analysis. ProClin >0.0015%, BSA 1%

REF DCE045/0903-0

3. Conjugate (1 vial, 12 mL)

Anti- LH antibodies conjugated with horseradish peroxidase (HRP) and biotinylated anti- LH. ProClin >0.0015%, BSA 0.1%

REF DCE002/0902-0

4. Coated Microplate (1 microplate breakable)

Streptavidin adsorbed on microplate

REF DCE002/0903-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB (0.26 g/L) (avoid any skin contact)

ProClin <0.0015% REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

7. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L; Tween-20 55 g/L, contains ProClin >0.015%

REF DCE006-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, control, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.

- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 – 8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22 – 28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.

- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the calibrators are stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2 – 8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum or plasma samples.

Sample Storage	Duration
2 – 8 °C	5 days
Freeze/thaw cycles	1 cycle

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer. The Calibrators are ready to use, are calibrated against the international reference WHO 1st IRP 68/40 and have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
mIU/mL	0	5	25	50	100	200

The control is ready to use; the concentration is printed on the label.

9.2. Preparation of the Conjugate

The conjugate is ready to use.

9.3. Preparation of the Wash Solution

Dilute the content of the vial “50X Conc. Wash Solution” with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated

wash solution to 1000 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

The determination of LH can be performed in human serum or plasma samples.

For accurate comparison to established normal values, a fasting morning serum or plasma sample should be obtained.

Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.5. Procedure

- **Allow all reagents to reach room temperature (22 – 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for the Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/Control	Blank
Calibrator C ₀ -C ₅	20 µL		
Sample/Control		20 µL	
Conjugate	100 µL	100 µL	

Incubate 1 h at room temperature +22 – 28°C. Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at 22 – 28°C for 15 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake gently the microplate.
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of each control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁸.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, **including Calibrator 0 is required**. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

12. EXPECTED VALUES

Reference ranges are hereby reported:

	Expected values (mIU/mL)
Males	0.7 – 7.4
Females	
Follicular phase	0.5 – 10.5
Ovulation phase	18.4 – 61.2
Luteal phase	0.5 – 10.5
Postmenopausal	8.2 – 40.8

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

13. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

13.1. Precision

13.1.1. Intra Assay Variation

Within run variation was determined by replicate (20x) the measurement of three different control sera in one assay. The within assay variability is $\leq 9.21\%$.

13.1.2. Inter Assay Variation

Between run variation was determined by replicate (10x) the measurements of three different control sera in different lots. The between assay variability is $\leq 7.91\%$.

13.2. Method comparison

LH ELISA was compared to a commercially available LH kit. 36 serum samples were tested.

The regression curve is:

n	Slope	Intercept (mIU/mL)	R ²
36	0.91	0.05	0.971

The new LH kit was compared to the old DiaMetra LH kit. 36 serum samples were tested.

The regression curve is:

n	Slope	Intercept (mIU/mL)	R ²
36	1.08	-1.22	0.981

13.3. Sensitivity

The lowest detectable concentration of LH that can be distinguished from the Calibrator 0 is 0.22 mIU/mL.

13.4. Specificity

The cross-reactivity of the LH ELISA to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between dose of interfering substance to dose of luteinising hormone needed to produce the same absorbance.

Substance tested	Cross Reactivity
LH	100 %
HCG	0.007 %
Prolactin	N.D.
FSH	N.D.
TSH	N.D.

13.5. Accuracy

The recovery test performed on three different samples, enriched with 5.63 - 11.25 - 22.5 - 45 - 90 mIU/mL of LH, gave a average value ($\pm$ SD) of 97.17% $\pm$ 4.00%.

In the dilution test three different samples were diluted 2, 4, 8 and 16 times with Calibrator 0; the average value ($\pm$ SD) obtained is 99.13% $\pm$ 7.37%.

13.6. Hook Effect

LH ELISA shows no Hook Effect up to 400 mIU/mL.

14. LIMITATIONS OF USE

- LH is suppressed by estrogen but in woman taking oral contraceptives the level may be low or normal. Excessive dieting and weight loss may lead to low gonadotropin concentrations. Luteinising hormone is dependent upon diverse factors other than pituitary homeostasis; thus, the determination alone is not sufficient to assess clinical status.
- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁹. Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

15. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

16. BIBLIOGRAPHY

1. Kosasa T.S., "Measurement of Human Luteinizing Hormone." *Journal of Reproductive Medicine*, 26 (1981) pg. 201-6.
2. Danzer H., Braunstein G.D., et al., "Maternal Serum Human Chorionic Gonadotropic Concentrations and Fetal Sex Predictions." *Fertility and Sterility*, 34 (1980) pg. 336-40.
3. Braunstein G.D., et al., "Serum Human Luteinizing Hormone Levels through Normal Pregnancy", *American Journal of Obstetrics and Gynecology*, 126 (1976) pg. 678-81.
4. Goldstein D.P., and Kosasa T.S., "The Subunit Radioimmunoassay for LH Clinical Application." *Gynecology*, 6 (1975) pg. 145-84.
5. Batzer F., "Hormonal Evaluation of Early Pregnancy", *Fertility and Sterility*, 34 (1980) pg. 1-12.
6. Braunstein, G.D., et al., "First-Trimester Luteinizing Hormone Measurements as an Aid to the Diagnosis of Early Pregnancy Disorders", *American Journal of Obstetrics and Gynecology*, 131 (1978) pg. 25-32.
7. Lenton E., Neal L. and Sulaiman R., "Plasma Concentrations of Human Gonadotropin from the time of Implantation until the Second Week of Pregnancy", *Fertility and Sterility*, 37 (1982) pg. 773-78.
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27-33

17. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

18. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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16. BIBLIOGRAFÍA

1. Kosasa T.S., "Measurement of Human Luteinizing Hormone." Journal of Reproductive Medicine, 26 (1981) pg. 201-6.
2. Danzer H., Braunstein G.D., et al., "Maternal Serum Human Chorionic Gonadotropic Concentrations and Fetal Sex Predictions." Fertility and Sterility, 34 (1980) pg. 336-40.
3. Braunstein G.D., et al., "Serum Human Luteinizing Hormone Levels through Normal Pregnancy", American Journal of Obstetrics and Gynecology, 126 (1976) pg. 678-81.
4. Goldstein D.P., and Kosasa T.S., "The Subunit Radioimmunoassay for LH Clinical Application." Gynecology, 6 (1975) pg. 145-84.
5. Batzer F., "Hormonal Evaluation of Early Pregnancy", Fertility and Sterility, 34 (1980) pg. 1-12.
6. Braunstein, G.D., et al., "First-Trimester Luteinizing Hormone Measurements as an Aid to the Diagnosis of Early Pregnancy Disorders", American Journal of Obstetrics and Gynecology, 131 (1978) pg. 25-32.
7. Lenton E., Neal L. and Sulaiman R., "Plasma Concentrations of Human Gonadotropin from the time of Implantation until the Second Week of Pregnancy", Fertility and Sterility, 37 (1982) pg. 773-78.
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. Clin Chem, 34, 1988, pp 27-33

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