DCM011-14

Ed. 03/2025

PROLACTIN ELISA

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PROLACTIN ELISA

for routine analysis

Direct immunoenzymatic determination of Prolactin in human serum.

IVD

LOT

See external label

2°C  8°C $\Sigma = 96$ tests

REF DKO011

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

Prolactin ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of Prolactin in human serum from an adult population.

2. CLINICAL SIGNIFICANCE

Prolactin is a polypeptide hormone synthesized and secreted by the Adenohypophysis (anterior Pituitary gland) and the placenta. It is also produced in other tissues including the breast and the decidua. Pituitary prolactin secretion is regulated by neuroendocrine neurons in the hypothalamus, most importantly by neurosecretory dopamine neurons of the arcuate nucleus, which inhibit prolactin secretion.

Prolactin is present in several body fluids, including blood plasma, amniotic fluid, milk, mucosal secretions and cerebrospinal fluid.

Prolactin has many effects, the most important of which is to stimulate the mammary glands to produce milk (lactation). Other possible functions of prolactin include the surfactant synthesis of the fetal lungs at the end of the pregnancy and immune tolerance of the foetus by the maternal organism during pregnancy.

Prolactin may also have inhibitory effects on gonadal function when present in high concentrations. There is a diurnal cycle in prolactin secretion.

During pregnancy, high circulating concentrations of estrogen promote prolactin production. The resulting high levels of prolactin secretion cause maturation of the mammary glands, preparing them for lactation. After childbirth, prolactin levels fall as the internal stimulus for them is removed.

High prolactin levels also tend to suppress the ovulatory cycle by inhibiting the secretion of both FSH and GnRH. Prolactin levels may be checked as part of a sex hormone workup, as elevated prolactin secretion can suppress the

secretion of FSH and GnRH, leading to hypogonadism, and sometimes causing erectile dysfunction in men.

Elevations in plasma prolactin concentrations occur during ovulation, pregnancy, nursing and stress. Abnormal elevations in plasma prolactin levels (hyperprolactinemia) can occur as a result of pituitary adenomas, other anatomic and traumatic abnormalities, in response to certain pharmacologic agents and in hypothyroidism.

Hyperprolactinemia (low prolactin levels) are observed in cases of hypopituitarism.


3. PRINCIPLE OF THE METHOD

The essential reagents required for an immunoenzymatic assay include high affinity and specificity antibodies (enzyme and immobilised) with different and distinct epitope recognition, in excess, and native antigen.

In this procedure the immobilization takes place during the assay at the surface of a microplate well through the interaction of streptavidin coated on the well and exogenously added biotinylated monoclonal anti-PRL antibody.

Upon mixing monoclonal biotinylated antibody, the enzyme labelled antibody and a serum containing the native antigen, a reaction results between the native antigen and the antibodies without competition or steric hindrance to form a soluble sandwich complex.

The interaction is illustrated by the following equation:

$\text{B}^{\text{tn}}\text{Ab}_{(m)}$ = Biotinylated Monoclonal Antibody (Excess quantity)

Ag_{PRL} = Native PRL antigen (variable quantity)

$\text{EnzAb}_{(p)}$ = Enzyme labeled polyclonal antibody (Excess quantity)


$\text{EnzAb}_{(p)} - \text{Ag}_{\text{PRL}} - \text{B}^{\text{tn}}\text{Ab}_{(m)}$ = Antigen-Antibodies Sandwich complex

K_a = Rate constant of association

K_{-a} = Rate constant of disassociation

Simultaneously the complex is deposited to the well through the high affinity reaction of streptavidin and biotinylated antibody.

This interaction is illustrated below:

Streptavidin_{cw} = Streptavidin immobilized on well

Immobilized complex = Antibodies-Antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody-bound fraction is directly proportional to the native antigen concentration. By using several different serum references of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each)

ProClin >0.0015%, Sodium Azide <0.1% and BSA 8%

CAL0	REF DCE002/1106-0
CAL1	REF DCE002/1107-0
CAL2	REF DCE002/1108-0
CAL3	REF DCE002/1109-0
CAL4	REF DCE002/1110-0
CAL5	REF DCE002/1111-0

2. Control (1 vial, 1 mL)

Concentration of Control is indicated on the Certificate of Analysis. ProClin >0.0015%, Sodium Azide <0.1% and BSA 8%

REF DCE045/1103-0

3. Conjugate (1 vial, 12 mL)

Anti- Prolactin antibody conjugated with horseradish peroxidase (HRP) and Anti Prolactin biotinylated antibody. ProClin >0.0015% and BSA 0.1%

REF DCE002/1102-0

4. Coated Microplate (1 microplate breakable)

Streptavidin adsorbed on microplate

REF DCE002/1103-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB (0.26 g/L) (avoid any skin contact)

ProClin <0.0015%

REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact)

REF DCE005-0

7. 10X Conc. Wash Solution (1 vial, 50 mL)

Phosphate buffer 0.2M pH 7.4. ProClin >0.0015%

REF DCE054-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.

- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.

-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.

- Some reagents (calibrators, control, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.

- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- Some reagents (calibrators and control) contain small amounts of Sodium Azide (NaN₃) <0.1% which may be toxic if ingested or absorbed through the skin or eyes; moreover, it may react with lead or copper plumbing to form potentially explosive metal azides. If you use a sink to remove the reagents, wash through with large amounts of water to prevent azide build-up.
- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 – 8°C in their original container. Any exceptions are clearly indicated.
- Allow all kit components and specimens to reach room temperature (22 – 28°C) and mix well prior to use.

- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the calibrators are stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2 – 8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum samples.

Sample Storage	Duration
2 – 8 °C	2 days
Freeze/thaw cycles	1 cycle

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer.

The Calibrators are ready to use, are calibrated against WHO 4th IS 83/573 and have the following concentrations:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
ng/mL	0	5	10	25	50	100

The control is ready to use; the concentration is printed on the label.

9.2. Preparation of the Conjugate

The conjugate is ready to use. Mix gently, for 5 minutes, on a roller mixer.

9.3. Preparation of the Wash Solution

Dilute the content of the vial “10X Conc. Wash Solution” with distilled water to a final volume of 500 mL prior to use. For smaller volumes respect the 1:10 dilution ratio.

It is possible to observe the presence of crystals within the concentrated wash solution; in this case mix at room temperature until the complete dissolution of crystals. For greater accuracy, dilute the whole bottle of concentrated wash solution to 500 mL, taking care also to transfer crystals completely by rinsing of the bottle, then mix until crystals are completely dissolved.

9.4. Preparation of Samples

The determination of Prolactin can be performed in human serum samples. **Do not use plasma for this assay – plasma sample may lead to false results.**

Fasting samples are not necessary and no special sample preparations are required.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.5. Procedure

- **Allow all reagents to reach room temperature (22 – 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for the Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Calibrator C ₀ -C ₅	50 µL		
Sample/ Control		50 µL	
Conjugate	100 µL	100 µL	

Incubate 1h at room temperature (22 – 28 °C).
Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at 22 – 28°C for 15 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake gently the microplate.
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of the control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories¹¹.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) or Sigmoid curve fit, **including Calibrator 0 is required.**

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

NOTE: If the absorbance of a sample exceeds that of C₅, the sample should be diluted 1:2 in C₀ and re-analysed.

12. EXPECTED VALUES

The serum Prolactin values are comprised in the following intervals:

Adults	Expected values (ng/mL)
Males	1.8 – 17.0
Females	
Premenopausal	1.2 – 19.5
Postmenopausal	1.5 – 19.5

Some of the female population tested in this group were probably using oral contraceptives which may affect results.

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

13. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

13.1. Precision

13.1.1. Intra Assay Variation

Within run variation was determined by replicate (20x) the measurement of three different control sera in one assay.

Sample	N	X	σ	C.V.
Level 1	20	5.33	0.15	2.78%
Level 2	20	18.21	0.73	4.03%
Level 3	20	37.20	1.38	3.71%

13.1.2. Inter Assay Variation

Between run variation was determined by replicate (10x) the measurement of three different control sera with kits of different lots.

Sample	N	X	σ	C.V.
Level 1	10	5.46	0.30	5.49%
Level 2	10	17.72	0.91	5.16%
Level 3	10	36.29	1.67	4.60%

13.2. Sensitivity

The lowest detectable concentration of Prolactin that can be distinguished from the Calibrator 0 is 0.12 ng/mL at the 95 % confidence limit.

13.3. Accuracy

The recovery of 3,13 - 6,25 - 12,50 - 25,00 - 50,00 ng/mL of Prolactin added to sample gave an average value ($\pm$ SD) of 102.52% $\pm$ 9.75% with reference to the original concentrations.

The dilution test performed on three sera diluted 2 - 4 - 8 - 16 times gave an average value ($\pm$ SD) of 102.19% $\pm$ 9.80%.

13.4. Specificity

The cross reaction of the antibody calculated at 50% according to Abraham are shown in the table:

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13.5. Correlation

Dia.Metra Prolactin ELISA kit (y) was compared to another commercially available Prolactin assay (x). Serum samples of 37 subjects were analysed.

The linear regression curve was calculated:

n	Slope	Intercept (ng/mL)	R ²
37	1.01	1.94	0.957

The new Dia.Metra Prolactin ELISA kit (y) was compared to the old Dia.Metra Prolactin ELISA kit (x). Serum samples of 37 subjects were analysed.

The linear regression curve was calculated:

n	Slope	Intercept (ng/mL)	R ²
37	0.85	2.58	0.937

13.6. Hook Effect

Dia.Metra Prolactin assay shows no Hook effect up to 200 ng/mL.

14. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays¹². Patients routinely exposed to animals or to animal serum products can be prone to this interference and anomalous values may be observed.

15. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

16. BIBLIOGRAPHY

1. Wisdom, G.B. Clin. Chem. 22/8 1243 - 1255 (1976)
2. Shome, B. and Parlow, A.F.,
3. J. Clin. Endocr. Metab.,39, 199 – 202 (1974).
4. Shome, B. and Parlow, A.F.,
5. J. Clin. Endocr. Metab.,39, 203 – 205 (1974).
6. Uotila, M. Ruoslahti, E. Engvall, E.,
7. J. Immunol. Methods, 42, 11-15 (1981)
8. Acosta, A.A.M.D. and Wright, G.L.,
9. Journal of Clinical Immunoassays, 6, 41 (1983).
10. Cohen, K. L., Metabolism, 26 1165-1177 (1977)

11. Basic QC Practices On-line Course;
<http://www.Westgard.com>.
12. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

17. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

18. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its

authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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