DCM045-11
Ed. 02/2024

T4 ELISA

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BIBLIOGRAFIA

1. Barker, S.B., "Determination of Protein Bound Iodine" Journal Biological Chemistry, 173, 175 (1948).
2. Chopra, I.J., Solomon, D.H., and Ho, R.S., "A Radioimmunoassay of Thyroxine," J.Clinical Endocrinol, 33, 865 (1971).
3. Young, D.S., Pestaner, L.C., and Gilberman, U., "Effects of Drugs on clinical Laboratory Tests." Clinical Chemistry, 21, 3660 (1975).
4. Sterling, L., Diagnosis and Treatment of Thyroid Disease, Cleveland CRC press, P.19-51 (1975).

Ed. 02/2024

DCM045-11

DCM045-11
Ed. 02/2024

T4 ELISA

for routine analysis

Direct immunoenzymatic determination of thyroxine (T4) in human serum or plasma.

IVD

LOT

See external label

2°C 8°C

$\Sigma = 96$ test

REF DKO045

INTENDED USE

Immunoenzymatic colorimetric method for quantitative determination of thyroxine (T4) concentration in human serum and plasma from an adult population.

T4 ELISA kit is intended for laboratory use only.

1. CLINICAL SIGNIFICANCE

The thyroid hormone, thyroxine (T4) is produced by the thyroid gland. An important component in the synthesis is iodine. The major form of thyroid hormone in the blood is thyroxine (T4). Thyroxine is converted to the active T₃ (three to four times more potent than T4) within cells by deiodinases (5'-iodinase).

Thyroxine-binding globulin (TGB) is the major carrier protein for circulating thyroid hormone. Only a very small fraction of the circulating hormone is free (unbound) - T4 0.03%.

The thyronines act on the body to increase the basal metabolic rate, affect protein synthesis and increase the body's sensitivity to catecholamines (such as adrenaline) by permissiveness. The thyroid hormones are essential to proper development and differentiation of all cells of the human body. These hormones also regulate protein, fat, and carbohydrate metabolism, affecting how human cells use energetic compounds. Numerous physiological and pathological stimuli influence thyroid hormone synthesis.

Thyroid hormones act through an unknown mechanism to inhibit neuronal activity; one of the effects is the reduction of the body temperature.


Thyrotoxicosis or hyperthyroidism is the clinical syndrome caused by an excess of circulating free thyroxine, free triiodothyronine, or both.

Both T3 and T4 are used to treat thyroid hormone deficiency (hypothyroidism).

2. PRINCIPLE

The T4 (antigen) in the sample competes with the antigenic T4 conjugated with horseradish peroxidase (HRP) for binding to the limited number of antibodies anti T4 coated on the microplate (solid phase) (the enzyme conjugate should have no measurable binding to serum proteins especially TBG and albumin).

The interaction is illustrated through the following equation:

Ab_{cw}: monospecific immobilised antibody (constant quantity)

Ag: native antigen (variable quantity)

EnzAg: antigen conjugated to enzyme HRP (constant quantity)

Ag Ab_{cw}: antigen-antibody complex

EnzAg Ab_{cw}: antigen-HRP-antibody complex

K_a: rate constant of association

K_{-a}: rate constant of disassociation

K = k_a / k_{-a}: equilibrium constant

After incubation, the bound/free separation is performed by a simple solid-phase washing.

Then, the enzyme HRP in the bound-fraction reacts with the Substrate (H₂O₂) and the TMB Substrate and develops a blue color that changes into yellow when the Stop Solution (H₂SO₄) is added.

The colour intensity is inversely proportional to the T4 concentration of in the sample. T4 concentration in the sample is calculated through a calibration curve.

3. REAGENTS, MATERIALS AND INSTRUMENTATION

3.1. Reagents and materials supplied in the kit

1. Calibrators (6 vials, 1 mL each)

CAL0	REF DCE002/4506-0
CAL1	REF DCE002/4507-0
CAL2	REF DCE002/4508-0
CAL3	REF DCE002/4509-0
CAL4	REF DCE002/4510-0
CAL5	REF DCE002/4511-0

2. Conjugate Buffer (1 vial, 12.5 mL)

Buffered solution with preservative and binding protein inhibitors
REF DCE002/4501-0

3. Conjugate (1 vial, 1.4 mL)

T4 conjugated with horseradish peroxidase (HRP)
(keep away from sunlight) REF DCE002/4502-0

4. Coated Microplate (1 microplates breakable)

Anti T4 antibody adsorbed on microplates
REF DCE002/4503-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB (0.26 g/L) (avoid any skin contact)
ProClin <0.0015% REF DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15 mol/L (avoid any skin contact)
REF DCE005-0

7. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L, Tween-20 55 g/L
ProClin >0.0015% REF DCE006-0

3.2. Reagents necessary not supplied

Distilled water.

3.3. Auxiliary materials and instrumentation

Automatic dispenser.

Microplates reader (450 nm, 620-630 nm)

Note

Store the reagents at 2-8°C in the dark. Open the bag of reagent 4 (Coated Microplate) only when it is at room temperature and close it immediately after use; once opened, the microplate is stable until the expiry date of kit.

4. WARNINGS

- This kit is intended for in vitro use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
- Some reagents (wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid the contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The TMB Substrate contains an irritant, which may be harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous and corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to directed sunlight, metals or oxidants. Do not freeze the solution.
- The concentration of total thyroxine in serum depends on several factors: the function of the thyroid gland and its regulation, thyroxine binding globulin (TBG) and thyroxine binding to TBG. However, the concentration of total thyroxine alone is not sufficient to monitor the clinical status.
- Total serum thyroxine values may increase during pregnancy or administration of oral contraceptives. The table of drugs and interfering conditions in which the values of total thyroxine are affected have been compiled by the Journal of the American Association of Clinical Chemists.
- Not intended for newborn screening

5. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2-8°C in their original container. Any exceptions are clearly indicated. The reagents are stable until the expiry date when stored and handled as indicated.
- Allow all kit components and specimens to reach room temperature (22-28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems, it is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

6. PROCEDURE

6.1. Preparation of the Calibrators (C₀...C₅)

The calibrators are ready to use, are calibrated against Human Serum Reference for thyroxine and have approximate concentrations of:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅
µg/dL	0	2.0	5.0	10.0	15.0	25.0

Calibrators concentration is lot-specific; the exact values are stated on labels and Certificate of Analysis for each lot. Once opened, the Calibrators are stable 6 months at 2-8°C.

6.2. Preparation of diluted Conjugate

Dilute the Conjugate 1:11 with Conjugate Buffer in a suitable container. For example, dilute 160 μ L of conjugate in 1.6 mL of conjugate buffer. This reagent should be used within 24 hours for maximum performance of the assay. Store at 2÷8°C.

IMPORTANT NOTE: prolonged exposure of the Conjugate to sunlight may affect the functional characteristics of the assay; therefore do not expose Conjugate (and diluted Conjugate) to direct sunlight.

6.3. Preparation of Wash Solution

Dilute the contents of each vial of the "50X Conc. Wash Solution" with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio. The diluted wash solution is stable for 30 days at 2÷8°C.

6.4. Preparation of the Sample

The determination of T4 should be performed in human serum or plasma.

Specimens may be refrigerated at 2÷8°C (for a maximum period of 48 hours). If the sample cannot be assayed within 48 hours, it may be stored at -20°C up to 30 days. Avoid repetitive freezing and thawing of samples. When assayed in duplicate, 0.050 mL of the specimen is required.

6.5. Procedure

- **Allow all reagents to reach room temperature (22-28°C) for at least 30 minutes.** At the end of the assay, store immediately the reagents at 2-8°C: avoid long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2-8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₅), two for each Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Sample/ Control		25 μ L	
Calibrator C ₀ -C ₅	25 μ L		
Diluted Conjugate	100 μ L	100 μ L	
Incubate 1 h at room temperature (22÷28°C). Remove the contents from each well, wash the wells three times with 300 μ L of diluted wash solution. Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel. Automatic washer: if you use automated equipment, wash the wells at least 5 times.			
TMB Substrate	100 μ L	100 μ L	100 μ L

Incubate 15 minutes in the dark at room temperature (22÷28°C).

Stop Solution

100 μ L

100 μ L

100 μ L

Shake the microplate gently.

Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

7. QUALITY CONTROL

Each laboratory should assay controls at levels in the hypothyroid, euthyroid and hyperthyroid range for monitoring assay performance. These controls should be treated as unknowns and values determined in every test procedure performed. Quality control charts should be maintained to follow the performance of the supplied reagents. Pertinent statistical methods should be employed to ascertain trends. The individual laboratory should set acceptable assay performance limits. Other parameters that should be monitored include the 80, 50 and 20% intercepts of the calibration curve for run-to-run reproducibility. In addition, maximum absorbance should be consistent with past experience. Significant deviation from established performance can indicate unnoticed change in experimental conditions or degradation of kit reagents. Fresh reagents should be used to determine the reason for the variations.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

8. RESULTS

8.1. Mean Absorbance

Calculate the mean of the absorbance (Em) for each point of the calibration curve (C₀-C₅) and of each sample.

8.2. Calibration curve

Plot the mean value of absorbance (Em) of the Calibrators (C₀-C₅) against concentration. Draw the best-fit curve through the plotted points. (es: Four Parameter Logistic).

8.3. Calculation of Results

Interpolate the values of the samples on the calibration curve to obtain the corresponding values of the concentrations expressed in μ g/dL.

9. REFERENCE VALUES

A study of euthyroid adult population was undertaken to determine expected values for T4 ELISA.

	Mean (μ g/dL)	SD	Range (μ g/dL)
Normal Values	7,6	1,6	4,4 - 10,8

Please pay attention to the fact that the determination of a range of expected values for a "normal" population in a given method is dependent on many factors, such as specificity and sensitivity of the method used and type of population under investigation. Therefore each laboratory should consider the range given by the Manufacturer as a general indication and produce their own range of expected values

based on the indigenous population where the laboratory works.

10. PERFORMANCE AND CHARACTERISTICS

10.1. Precision

10.1.1. Intra Assay Variation

Within run variation was determined by replicate (20x) the measurement of three different control sera in one assay. The within assay variability is $\leq 8.16\%$.

10.1.2. Inter Assay Variation

Between run variations was determined by replicate (10x) the measurement of three different control sera in different lots of kit. The between assay variability is $\leq 8.42\%$.

10.2. Accuracy

The recovery of 1.25 - 2.50 - 5.00 - 10.00 $\mu\text{g/dL}$ of T4 added to sample gave an average value ($\pm\text{SD}$) of $97.8\% \pm 5\%$ with reference to the original concentrations.

10.3. Sensitivity

The lowest detectable concentration of T4 that can be distinguished from the Calibrator 0 is 0.4 $\mu\text{g/dL}$ at the 95% confidence limit.

10.4. Specificity

The cross-reactivity of the triiodothyronine antibody to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by deriving a ratio between doses of interfering substance to dose of Triiodothyronine needed to displace the same amount of tracer.

Substance	Cross Reactivity	Concentration
I-Thyroxine	1.0000	-
d-Thyroxine	0.9800	10 $\mu\text{g/dL}$
I-Triiodo-thyronine	0.0300	100 $\mu\text{g/dL}$
d-Triiodo-thyronine	0.0150	100 $\mu\text{g/dL}$
Monoiodo-Tyrosine	N/D	100 $\mu\text{g/mL}$
Diiodo-Tyrosine	N/D	100 $\mu\text{g/mL}$
Triiodothyroacetic Acid	N/D	100 $\mu\text{g/mL}$
Tetraiodothyroacetic Acid	0.0001	100 $\mu\text{g/mL}$

10.5. Correlation with RIA

Dia.Metra T4 ELISA was compared to another commercially available T4 assay. 114 serum samples were analysed according in both test systems.

The linear regression curve was calculated:

$$(\text{T4 Dia.Metra}) = 0.959 \cdot (\text{T4 RIA}) + 0.339$$

$$r^2 = 0.936$$

11. WASTE MANAGEMENT

Reagents must be disposed off in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

BIBLIOGRAPHY

1. Barker, S.B., "Determination of Protein Bound Iodine" *Journal Biological Chemistry*, **173**, 175 (1948).
2. Chopra, I.J., Solomon, D.H., and Ho, R.S., "A Radioimmunoassay of Thyroxine," *J. Clinical Endocrinol.*, **33**, 865 (1971).
3. Young, D.S., Pestaner, L.C., and Gilberman, U., "Effects of Drugs on clinical Laboratory Tests." *Clinical Chemistry*, **21**, 3660 (1975).
4. Sterling, L., *Diagnosis and Treatment of Thyroid Disease*, Cleveland CRC press, P.19-51 (1975).

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