DCM013-13

Ed. 05/2025

TSH ELISA

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TSH ELISA

for routine analysis

Direct immunoenzymatic determination of TSH in human serum or plasma

IVD

LOT

See external label

2°C  8°C $\Sigma = 96$ tests

REF DKO013

1. INTENDED PURPOSE

For *In Vitro* Diagnostic Use

For Laboratory Professional Use

TSH ELISA is a manual *in vitro* diagnostic device intended for the quantitative determination of thyroid stimulating hormone in human serum or plasma from an adult population.

2. CLINICAL SIGNIFICANCE

Thyroid-stimulating hormone (TSH or thyrotropin) is a glycopolypeptide hormone synthesised and secreted by the anterior pituitary gland which regulates the endocrine function of the thyroid gland. TSH consists of two subunits - *alpha* and *beta*. The α subunit is identical to that of human chorionic gonadotropin (HCG), luteinising hormone (LH) and follicle-stimulating hormone (FSH). The β (*beta*) subunit is unique to TSH, and therefore determines its function.

TSH stimulates the thyroid gland to secrete the hormones thyroxine (T4) and triiodothyronine (T3).

TSH production is controlled by a Thyrotropin Releasing Hormone (TRH). TRH is manufactured in the hypothalamus and transported to the pituitary gland, where it increases TSH production and release. Somatostatin is also produced by the hypothalamus and has an opposite effect on the pituitary production of TSH, decreasing or inhibiting its release.

Release of TSH is regulated by the circulating free fraction of thyroid hormones in the blood. TSH levels are depressed when peripheral concentrations of the free fraction of thyroid hormones are high. Conversely, TSH levels are high when peripheral concentrations of thyroid hormones are low. This effect creates a regulatory negative feedback loop.

TSH levels are tested in the blood of patients suspected of suffering from excess or deficiency of thyroid hormone hyperthyroidism or hypothyroidism respectively). Generally, a normal range for TSH is between 0.3 and 3.0 mIU/L, but the interpretation depends also on what the blood levels of thyroid hormones (T3 and T4) are.

Higher than normal levels of TSH combined with high levels of thyroid hormone (T3 and T4) may indicate dysfunction of the hypothalamus and pituitary gland. In this case, a high TSH is often produced by a benign tumour of the pituitary (adenoma). Conversely, low levels of TSH, while blood levels of T3 and T4 are also low, indicates abnormally low function of the pituitary, known as hypopituitarism.

Due to the negative feedback detailed above, abnormally high levels of Thyroid hormone, caused by overproduction in the thyroid, results in low TSH levels. This occurs in diseases such as hyperthyroidism or Grave's disease. Conversely, an underproduction of T3 and T4 caused by diseases such as congenital hypothyroidism (cretinism), hypothyroidism or thyroid hormone resistance, gives rise to an increase in the measured TSH.


Clearly both TSH and T3 and T4 should be measured to ascertain where a specific thyroid disfunction is caused by primary pituitary or by a primary thyroid disease.

3. PRINCIPLE OF THE METHOD

The essential reagents required for an immunoenzymatic assay include high affinity and specificity antibodies (enzyme and immobilised) with different and distinct epitope recognition, in excess, and native antigen.

In this method, TSH calibrators, patient specimens and/or controls containing the native antigen are first added to streptavidin coated wells. Biotinylated monoclonal and enzyme labelled antibodies are added and the reactants mixed: these antibodies have high affinity and specificity and are detected against distinct and different epitopes of TSH. Reaction between the various TSH antibodies and native TSH occurs in the microwells without competition or steric hindrance forming a soluble sandwich complex.

The interaction is illustrated by the following equation:

BtAb(m) = Biotinylated Monoclonal Antibody (Excess Quantity)

AgTSH = Native Antigen (Variable Quantity)

EAb(p) = Enzyme labeled Antibody (Excess Quantity)

EAb(p)-AgTSH-BtAb(m) = Antigen-Antibodies Sandwich Complex

K_a = Rate Constant of Association

K_{-a} = Rate Constant of Dissociation

Simultaneously, the complex is deposited to the well though the high affinity reaction of streptavidin and biotinylated antibody. This interaction is illustrated below:

EAb(p)-AgTSH-BtAb(m) + Streptavidin CW. ⇨ Immobilised Complex

Streptavidin CW. = Streptavidin immobilized on well
Immobilised Complex = Antibodies-Antigen sandwich bound.

After equilibrium is attained, the antibody-bound fraction is separated from unbound antigen by decantation or aspiration. The enzyme activity in the antibody bound fraction is directly proportional to the native antigen concentration. The activity of the enzyme present on the surface of the well is quantified by reaction with a suitable substrate to produce colour. By utilizing several different calibrators of known antigen values, a dose response curve can be generated from which the antigen concentration of an unknown can be ascertained.

4. REAGENTS, MATERIALS AND INSTRUMENTATION

4.1. Reagents and materials supplied in the kit

1. Calibrators (7 vials, 1 mL each)

ProClin >0.0015%

CAL0	REF DCE002/1306-0
CAL1	REF DCE002/1307-0
CAL2	REF DCE002/1308-0
CAL3	REF DCE002/1309-0
CAL4	REF DCE002/1310-0
CAL5	REF DCE002/1311-0
CAL6	REF DCE002/1312-0

2. Control (1 vial, 1 mL each)

Concentration of Control is indicated on the Certificate of Analysis. ProClin >0.0015% **REF** DCE045/1303-0

3. Conjugate (1 vial, 13 mL)

Goat anti- TSH conjugated with horseradish peroxidase (HRP); mouse anti- TSH biotinylated. ProClin >0.0015% BSA 0.1% **REF** DCE002/1302-0

4. Coated Microplate (1 microplate breakable)

Streptavidin adsorbed on microplate **REF** DCE002/1303-0

5. TMB Substrate (1 vial, 15 mL)

H₂O₂-TMB (0.26 g/L) (avoid any skin contact)

ProClin <0.0015% **REF** DCE004-0

6. Stop Solution (1 vial, 15 mL)

Sulphuric acid 0.15M (avoid any skin contact) **REF** DCE005-0

7. 50X Conc. Wash Solution (1 vial, 20 mL)

NaCl 45 g/L, Tween-20 55 g/L. ProClin >0.0015% **REF** DCE006-0

4.2. Materials required but not provided

Distilled water

4.3. Auxiliary materials and instrumentation

Automatic dispenser

Precision Pipetting Devices

Microplate reader (450 nm, 620-630 nm)

5. WARNINGS

- This kit is intended for *in vitro* use by professional persons only. Not for internal or external use in Humans or Animals.
- Use appropriate personal protective equipment while working with the reagents provided.
- Follow Good Laboratory Practice (GLP) for handling blood products.
-  Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy and the bovine protein has been obtained from countries not infected by BSE, but these materials should be handled as potentially infectious.
- Some reagents (calibrators, control, conjugate and wash solution) contain small amounts of ProClin™ 300 (>0.0015%, <0.06%) as preservative. Avoid contact with skin or mucosa.
- Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin sensitivity, Category 1

Warning

Contains: ProClin 300

Hazard statements:

H317 - May cause an allergic skin reaction.

Precautionary statements:

P261 - Avoid breathing dust / fume / gas / mist / vapours / spray.

P280 - Wear protective gloves/ protective clothing / eye protection / face protection / hearing protection.

P321 - Specific treatment (see supplemental first aid instruction on this label).

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 - Take off contaminated clothing and wash it before reuse.

- The TMB Substrate contains an irritant, which harmful if inhaled, ingested or absorbed through the skin. To prevent injury, avoid inhalation, ingestion or contact with skin and eyes.
- The Stop Solution consists of a diluted sulphuric acid solution. Sulphuric acid is poisonous, corrosive and can be toxic if ingested. To prevent chemical burns, avoid contact with skin and eyes.
- Avoid the exposure of reagent TMB/H₂O₂ to direct sunlight, metals or oxidants. Do not freeze the solution.

6. PRECAUTIONS

- Please adhere strictly to the sequence of pipetting steps provided in this protocol. The performance data represented here were obtained using specific reagents listed in this Instruction For Use.
- All reagents should be stored refrigerated at 2 – 8°C in their original container. Any exceptions are clearly indicated.

- Allow all kit components and specimens to reach room temperature (22 – 28°C) and mix well prior to use.
- Do not interchange kit components from different lots. The expiry date printed on box and vials labels must be observed. Do not use any kit component beyond their expiry date.
- If you use automated equipment, the user has the responsibility to make sure that the kit has been appropriately validated for its intended use/purpose.
- The incomplete or inaccurate liquid removal from the wells could influence the assay precision and/or increase the background. To improve the performance of the kit on automatic systems is recommended to increase the number of washes.
- It is important that the time of reaction in each well is held constant for reproducible results. Pipetting of samples should not extend beyond ten minutes to avoid assay drift. If more than 10 minutes are needed, follow the same order of dispensation. If more than one plate is used, it is recommended to repeat the dose response curve in each plate
- Addition of the TMB Substrate solution initiates a kinetic reaction, which is terminated by the addition of the Stop Solution. Therefore, the TMB Substrate and the Stop Solution should be added in the same sequence to eliminate any time deviation during the reaction.
- Observe the guidelines for performing quality control in medical laboratories by assaying controls and/or pooled sera.
- Maximum precision is required for reconstitution and dispensation of the reagents.
- Samples microbiologically contaminated, highly lipemic, icteric or haemolysed should not be used in the assay.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Fresh disposable tips must be used when pipetting assay reagents including samples, calibrators and controls to mitigate the risk of carryover contamination. Failure to do so may lead to invalid results.

7. REAGENT STORAGE AND STABILITY

Store the kit at 2 – 8°C in the dark.

- The kit is stable at 2 – 8°C until the expiry date stated on the external kit label.
- Once opened, the calibrators are stable at 2 – 8°C for 6 months.
- The diluted wash solution is stable for 30 days at 2 – 8°C.

Important note: open the bag containing the Coated Microplate only when it is at room temperature and close it immediately after use.

8. SAMPLE COLLECTION AND STORAGE

The assay should be performed using serum or plasma samples.

Sample Storage	Duration
2 – 8 °C	24 hours
Freeze/thaw cycles	1 cycle

9. PROCEDURE

9.1. Preparation of Calibrators and Controls

Before use, mix for 5 minutes with a rotating mixer.

The calibrators are ready for use and have the following concentration:

	C ₀	C ₁	C ₂	C ₃	C ₄	C ₅	C ₆
mIU/L	0	0.2	0.5	2.5	5.0	10	20

The control is ready to use; the concentration is printed on the label.

9.2. Preparation of the Conjugate

The conjugate is ready to use.

9.3. Preparation of the Wash Solution

Dilute the content of the vial “50X Conc. Wash Solution” with distilled water to a final volume of 1000 mL prior to use. For smaller volumes respect the 1:50 dilution ratio.

9.4. Preparation of Samples

The determination of TSH can be performed in human serum or plasma samples.

Fasting samples are not necessary and no special sample preparations are required.

Collect blood by venepuncture into vacutainers and separate serum (after clot formation) or plasma from the cells by centrifugation.

Store the sample at -20°C if the determination is not performed on the same day of the sample collection. Before using, mix gently, for 5 minutes, with a roller mixer.

9.5. Procedure

- **Allow all reagents to reach room temperature (22 – 28 °C) for at least 30 minutes.** At the end of the assay, immediately store the reagents at 2 – 8°C: avoiding long exposure to room temperature.
- Unused coated microwell strips should be released securely in the foil pouch containing desiccant and stored at 2 – 8°C.
- To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.
- As it is necessary to perform the determination in duplicate in order to improve accuracy of the test results, prepare two wells for each point of the calibration curve (C₀-C₆), two for the Control, two for each sample, one for Blank.

Reagent	Calibrator	Sample/ Control	Blank
Calibrator C ₀ -C ₆	50 µL		
Sample/ Control		50 µL	
Conjugate	100 µL	100 µL	

Incubate 1 h at room temperature (22 – 28 °C).
To increase the sensibility read the absorbance after 120 min of incubation.

Remove the contents from each well, wash the wells 3 times with 300 µL of diluted wash solution.

Important note: during each washing step, gently shake the plate for 5 seconds and remove excess solution by tapping the inverted plate on an absorbent paper towel.

Automatic washer: if you use automated equipment, wash the wells at least 5 times.

TMB Substrate	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Incubate at 22 – 28°C for 20 minutes in the dark.

Stop Solution	100 µL	100 µL	100 µL
---------------	--------	--------	--------

Shake gently the microplate.
Read the absorbance (E) at 450 nm against a reference wavelength of 620-630 nm or against Blank within 5 minutes.

10. QUALITY CONTROL

Good Laboratory Practice (GLP) requires the use of quality control specimens in each series of assays in order to check the performance of the assay. Controls should be treated as unknown samples, and the results analysed with appropriate statistical methods.

The kit control provided in the kit should be tested as unknown and is intended to assist in assessing the validity of results obtained with each assay plate.

The mean concentration of the control level is documented in the QC report included with each kit. These mean concentration levels are determined over several assays which are run in duplicate in multiple locations across each plate.

DiaMetra recommends the users to maintain graphic records of the control values generated with each assay run, including the running means, SDs and %CVs. This information will facilitate the controls trending analysis relating to the performance of current and historical control lots relative to the supplied Quality Control data. The trending will assist in the identification of assays which give

control values significantly different from their average range.

When interpreting control data, users should note that this product was designed and developed as a manual product. The range stated on the QC certificate should be appropriate for assays that are performed manually and with strict adherence to the Assay Procedure described above. It is recognised by Quality Control professionals, that as a result of differences in conditions and practices, there will always be variability in the mean values and precision of control measurements between different laboratories⁸.

11. CALCULATION OF RESULTS

A variety of data reduction software packages are available, which may be employed to generate the mean calibration curve and to calculate the mean concentrations of unknown samples and controls. A 4-parameter logistic (4PL) curve fit, or cubic spline **including Calibrator 0 is required**. Other curve fitting algorithms are not recommended.

Alternatively, a calibration curve may be prepared on semi-log graph paper by plotting mean absorbance on the Y-axis against concentration of analyte on the X-axis. Calibrator 0 should be included in the calibration curve. Read the mean absorbance value of each unknown sample off the curve.

In order for the assay results to be considered valid the kit calibrators and control must fall within the specifications detailed in the lot specific certificate of analysis.

If a control is out of its specified range, the associated test results are invalid and samples must be retested.

NOTE: If the absorbance of a sample exceeds that of C₆, the sample should be diluted in C₀ and re-analysed. Correct the result using an appropriate dilution factor.

12. EXPECTED VALUES

Serum samples of apparently healthy women and men were assayed using the TSH ELISA test, with the following results:

	Mean (mIU/L)	Reference interval (mIU/L)
Adults	1.85	0.39 – 6.16

The above ranges should be considered as guidelines only; it is recommended that each laboratory establish its own expected range based upon its own patient population.

13. PERFORMANCE CHARACTERISTICS

Representative performance data are shown. Results obtained at individual laboratories may vary.

13.1. Precision

13.1.1. Intra Assay Variation

Within run variation was determined by replicate measurements (16x) of two different control sera in one assay. The within assay variability is ≤ 4.6%.

13.1.2. Inter Assay Variation

Between run variation was determined by replicate measurements (12x) of three different control sera in different lots of kit. The between assay variability is $\leq 10.8\%$.

13.2. Accuracy

The recovery of 0.25 – 2.5 – 10 mIU/L of TSH added to sample gave an average value ($\pm$ SD) of $85.85\% \pm 2.62\%$ with reference to the original concentrations.

13.3. Sensitivity

The lowest detectable concentration of TSH that can be distinguished from the Calibrator 0 is 0.01 mIU/L with a confidence limit of 95%.

13.4. Specificity

The cross-reactivity of the thyrotropin ELISA method to selected substances was evaluated by adding the interfering substance to a serum matrix at various concentrations. The cross-reactivity was calculated by the ratio between the dose of the interfering substance and the amount of TSH needed to obtain the same absorbance.

Thyrotropin (hTSH)	100 %
Thyrotropin - alpha (hTSH- α)	9.3 %
Thyrotropin - beta (hTSH- β)	152 %
Follicle-Stimulating Hormone (hFSH)	1.0 %
Luteinising hormone (hLH)	1.0 %
Chorionic Gonadotropin (hCG)	<0.1 %

13.5. Method comparison

TSH ELISA was compared with two other commercially available TSH assays. 54 serum samples were tested. The two linear regression curves were calculated:

Assay type	Slope	Intercept (mIU/L)	R ²
CLIA	0.96	0.06	0.972
ELISA	0.87	0.29	0.971

13.6. Hook effect

The hook effect was tested using concentrations of analyte up to 500 mIU/L. No hook effect was observed.

14. LIMITATIONS OF USE

- As in the case of any diagnostic procedure, results must be interpreted in conjunction with the patient's clinical presentation and other information available to the physician.
- The performance characteristics of this assay have not been established in a paediatric population.
- Heterophilic antibodies in human serum can react with reagent immunoglobulins, interfering with *in vitro* immunoassays⁹. Patients routinely exposed to animals or to animal serum products can be

prone to this interference and anomalous values may be observed.

15. WASTE MANAGEMENT

Reagents must be disposed of in accordance with local regulations.

All materials that have come into contact with samples and reagents must be disposed of in accordance with country, state and local regulations.

16. BIBLIOGRAPHY

1. Hopton M. R, et al Clinical Chemistry, 32, 691 (1986)
2. Caldwell, G. et al. Lancet, I, 1117, (1985)
3. Young, D. S, et al Clinical Chemistry 21, 3660 (1975)
4. Spencer, CA, et al Clinical Chemistry 41, 367 (1995)
5. Beck-Peccoz P, et al Eur. J. Endocrinol 131: 331-340 (1994)
6. Bravemann, L. E, Clinical Chemistry 42: 174-181 (1996)
7. Fisher, D. A., Clinical Chemistry 42: 135-139 (1996)
8. Basic QC Practices On-line Course; <http://www.Westgard.com>.
9. Boscato, LM. and Stuart, MC., 'Heterophilic antibodies: a problem for all immunoassays'. *Clin Chem*, 34, 1988, pp 27–33

17. REVISION IDENTIFIER

Additions or changes to the IFU are indicated by grey highlighting.

18. PRODUCT COMPLAINTS AND TECHNICAL SUPPORT

For a patient/user/third party in the European Union and in countries with similar regulatory regime (Regulation 2017/746/EU on IVD Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorised representative and to your national regulatory authority.

The manufacturer can be contacted through their customer service or technical support team. The contact details can be found below and on the company website: www.diametra.com.

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