
STANDARD F
C. difficile Toxin A/B FIA
STANDARD™ F *C. difficile* Toxin A/B FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

KIT CONTENTS

MATERIALS REQUIRED BUT NOT PROVIDED

- STANDARD F Analyzer

SPECIMEN COLLECTION AND PREPARATION



■ Feces

- Fecal specimen should be collected in a clean and dry container and collected at any time of the day may be used.
- Urine should be excluded from the fecal specimen.
- Fresh fecal specimen should be used in this test. Do not use any transport media to store and transport the fecal specimen.
- Fecal specimen may be stored in a refrigerator at 2-8°C for up to 24 hours. For prolonged storage, specimens may be frozen under -20°C.
- It should be brought to room temperature prior to use.

PREPARATION AND TEST PROCEDURE

■ Preparation

- Allow kit components and collected sample to room temperature at least 30 minutes before starting the test.
- Carefully read instructions for STANDARD F *C. difficile* Toxin A/B FIA.
- Check the expiry date at the back of the foil pouch. Use another lot if expiry date has passed.
- Open the foil pouch, and check the test device and the silica gel pack inside the foil pouch.

- Do not write on the barcode or damage the barcode of the test device.
- If the color of silica gel is changed from yellow to green, please do not use the test device.

■ Collection of sample

- Poke the 4 different areas of the fecal sample using the sterile swab to collect a portion of feces (about 40-70mg). In case of liquid sample, soak the sterile swab completely in the liquid fecal sample.

- The amount of fecal swab may affect the results. It is required to follow the swab amount of feces as shown in the picture above. Excessive fecal amount may induce a false positive result and slow migration.

- Take off the top of an extraction buffer tube and insert the swab into the tube.
- Swirl the swab at least ten times into the extraction buffer tube to dissolve the sample.
- Remove the swab and dispose of the used swab in accordance with your biohazard waste disposal protocol.
- Tightly press down the filter cap onto the tube.

- Do not use the extraction buffer tube without the filter cap. It affects the result value.

■ Analysis of specimen


STANDARD TEST™ mode STANDARD F100, F200 and F2400 analyzer

- Prepare a STANDARD F Analyzer and select the 'Standard Test' mode according to the analyzer's manual.

STANDARD F2400 analyzer	Workplace -> 'Run Test' -> Insert patient ID and / or operator ID on the analyzer
STANDARD F100 and F200 analyzer	'Standard Test' mode -> Insert patient ID and / or operator ID on the analyzer

- Take the test device out of the foil pouch and place it on a flat and dry surface. Write patient information on the label of test device.

- Insert the test device to the test slot of the analyzer. When inserting the test device to the analyzer, the analyzer will read the barcode data, and check the test device is valid. Select the sample type on the sample selection screen.

- Apply 3 drops of mixed sample to the sample well in the test device holding the prepared extraction buffer tube upside down.

- After applying the buffer, immediately press the 'TEST START' button.

- The analyzer will automatically display the test result within 15 minutes.

- Do not put bubbles in the sample well of the test device.
- Discard the used test kit in proper container, according to your healthcare professional's recommendation.
- If you feel the test result is inaccurate or do not agree your test result, do not change the treatment and contact your healthcare professional.

INTERPRETATION OF TEST RESULT

Result	COI (Cutoff Index) value	Interpretation
C. difficile Toxin A Positive	COI ≥ 1.0	Detection of C. difficile Toxin A
C. difficile Toxin A Negative	COI < 1.0	No detectable C. difficile Toxin A detected.
C. difficile Toxin B Positive	COI ≥ 1.0	Detection of C. difficile Toxin B
C. difficile Toxin B Negative	COI < 1.0	No detectable C. difficile Toxin B detected.
Invalid	N/A	Retest with a new test device and a new patient sample

- Result must be interpreted together with STANDARD F *C. difficile* GDH FIA.
- False positive result can be caused by intake of antibiotics or other medicine.
- Retest should be performed with newly collected fecal specimen.

EXPLANATION AND SUMMARY

■ Introduction

Clostridium difficile infection (CDI) has significant clinical impact especially on the elderly and/or immunocompromised patients. The pathogenicity of *Clostridium difficile* is mainly mediated by two exotoxins: toxin A (TcdA) and toxin B (TcdB). Toxin A is mildly cytotoxic but induces large fluid shifts and mucosal inflammation. Toxin B is intensely cytotoxic but its role in the disease process is not clearly understood. Variant strains which are Toxin A-negative, Toxin B-positive are known to exist, are fully pathogenic, and capable of producing the full spectrum of disease. The prevalence of these variant strains varies widely by institution and geographic location. The clinical presentation of *C. difficile* infection includes, in increasing order of severity, range from mild, self-limited watery diarrhea to fulminant pseudomembranous colitis, toxic megacolon, and death. This pathogen is capable of causing disease that could be severe or fatal if not diagnosed on time and treated.

■ Intended use

STANDARD F *C. difficile* Toxin A/B FIA is an *in vitro* diagnostic use to qualitative measure the *C. difficile* Toxin A/B. This is only for professional use. It provides only an initial screening test result. STANDARD F *C. difficile* Toxin A/B FIA should be used with the appropriate analyzer, STANDARD F Analyzers, manufactured by SD BIOSENSOR.

■ Test principle

STANDARD F *C. difficile* Toxin A/B FIA is based on the immunofluorescence technology with STANDARD F Analyzer manufactured by SD BIOSENSOR to measure the *C. difficile* Toxin A/B in human sample. The specimen from human should be processed for the preparation using the components of STANDARD F *C. difficile* Toxin A/B FIA. After applying the sample mixture to the test device, the complex will be formed on the membrane as the result of the antigen-antibody reaction. The intensity of fluorescence light produced on the membrane. STANDARD F Analyzers can analyze the *C. difficile* Toxin A/B of the clinical specimen based on a pre-programmed algorithms and display the test result on the screen.

KIT STORAGE AND STABILITY

Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

WARNINGS AND PRECAUTIONS

- STANDARD F *C. difficile* Toxin A/B FIA is for *in vitro* diagnostics use only.
- Carefully follow the instructions and procedures described in this instructions before testing.
- STANDARD F *C. difficile* Toxin A/B FIA should be with STANDARD F Analyzer.
- STANDARD F *C. difficile* Toxin A/B FIA should remain in its original sealed pouch until ready to use. Do not use if the pouch is damaged or the seal is broken.
- STANDARD F *C. difficile* Toxin A/B FIA is single use only. Do not re-use it.
- Do not use frozen sample or artificial materials.
- Place the analyzer on a flat surface when in use.
- Wash your hands in warm, soapy water. Rinse well and dry completely before testing.
- Discard the used test kit according to the proper method.
- Check the expiration date printed on the pouch or package.
- Use STANDARD F *C. difficile* Toxin A/B FIA at 15-32°C / 59-90°F.
- All kit components are must be at room temperature 30 minutes before running the assay.
- Do not write on the barcode or damage the barcode of the test device.

QUALITY CONTROL

■ Calibration

The calibration set test of STANDARD F Analyzers should be conducted according to the analyzers' manual.

[When to use calibration set]

- Before using the analyzer for the first time
- When you drop the analyzer
- Whenever you do not agree with your result
- When you want to check the performance of an analyzer and test device

[How to use calibration set]

Calibration set test is a mandatory function that ensures optimal performance by checking the internal analyzer optics and functions.

- Select the 'Calibration' menu.
 - The specific calibration set is included with the analyzer.
 - Insert the CAL-1 first, and then insert the CAL-2 for UV-LED testing and the CAL-3 for RGB-LED testing in order.
- STANDARD F analyzer automatically calibrate and identify the optical performance through measuring the membrane of the test device whenever the test is conducted in 'Standard Test' mode. If 'EEE' message displays on the screen, it means that the analyzer has a problem, so check with test device. Contact the SD BIOSENSOR local distributor if the 'EEE' message still appears.

■ External quality control

Quality control testing should be run to check the performance of STANDARD F *C. difficile* Toxin A/B FIA and STANDARD F *C. difficile* Toxin A/B Control manufactured by SD BIOSENSOR should be used for quality control testing. Control test should be conducted in accordance with the instructions of STANDARD F *C. difficile* Toxin A/B Control.

Control test should be run:

- once for each new lot.
- once for each untrained operator.
- as required by test procedures in instructions for use of STANDARD F *C. difficile* Toxin A/B Control and in accordance with local, state and federal regulations or accreditation requirements.

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■ Clinical evaluation

1. Limit of Detection

The limit of detection of STANDARD F *C. difficile* Toxin A/B FIA is as follows;

Antigen	Limit of Detection
Toxin A (Sigma)	0.49 ng/ml
Toxin B (Bionote)	0.54 ng/ml

2. Precision

STANDARD F *C. difficile* Toxin A/B FIA totally meet the acceptance criteria of repeatability and reproducibility.

1) Toxin A

Operator		Operator 1			Operator 2			Operator 3		
Site		Site 1			Site 2			Site 3		
Lot		Lot 1	Lot 2	Lot 3	Lot 1	Lot 2	Lot 3	Lot 1	Lot 2	Lot 3
Panel No.	CIFF-TOXIN A/B-1801-P02	100%	100%	100%	100%	100%	100%	100%	100%	100%
	CIFF-TOXIN A/B-1801-P01	98%	97%	98%	98%	98%	98%	98%	97%	97%
	CIFF-TOXIN A/B-1801-N00	0%	0%	0%	0%	0%	0%	0%	0%	0%

2) Toxin B

Operator		Operator 1			Operator 2			Operator 3		
Site		Site 1			Site 2			Site 3		
Lot		Lot 1	Lot 2	Lot 3	Lot 1	Lot 2	Lot 3	Lot 1	Lot 2	Lot 3
Panel No.	CIFF-TOXIN A/B-1801-P02	100%	100%	100%	100%	100%	100%	100%	100%	100%
	CIFF-TOXIN A/B-1801-P01	98%	98%	97%	97%	98%	100%	97%	98%	97%
	CIFF-TOXIN A/B-1801-N00	0%	0%	0%	0%	0%	0%	0%	0%	0%

3. Relative sensitivity & specificity (LFIA product compared)

Relative sensitivity : 95% (64/67)

Relative specificity : 100% (70/70)

4. Cross Reactivity Study

STANDARD F *C. difficile* Toxin A/B FIA did not show any sign of cross reactivity with the concentrations list of each cross-reacting microorganisms.

#	Cross-reactivity substance	Concentration
1	Adenovirus type 1	2.0x10 ⁶ cfu/mL
2	Adenovirus type 2	2.0x10 ⁶ cfu/mL
3	Adenovirus type 3	2.0x10 ⁶ cfu/mL
4	Adenovirus type 4	2.0x10 ⁶ cfu/mL
5	Adenovirus type 5	2.0x10 ⁶ cfu/mL
6	Adenovirus type 7	2.0x10 ⁶ cfu/mL
7	Adenovirus type 11	2.0x10 ⁶ cfu/mL
8	Adenovirus type 18	2.0x10 ⁶ cfu/mL
9	Adenovirus type 23	2.0x10 ⁶ cfu/mL
10	Adenovirus type 40	2.0x10 ⁶ cfu/mL
11	Adenovirus type 55	2.0x10 ⁶ cfu/mL
12	Bacteriodes fragilis	2.0x10 ⁶ cfu/mL
13	Candida albicans	2.0x10 ⁶ cfu/mL
14	Escherichia coli	2.0x10 ⁶ cfu/mL
15	Klebsiella pneumoniae	2.0x10 ⁶ cfu/mL
16	Peptostreptococcus anaerobius	2.0x10 ⁶ cfu/mL
17	Pseudomonas aeruginosa	2.0x10 ⁶ cfu/mL
18	Staphylococcus aureus	2.0x10 ⁶ cfu/mL
19	Staphylococcus epidermidis	2.0x10 ⁶ cfu/mL

5. Interference Substances

STANDARD F *C. difficile* Toxin A/B FIA is not affected by below substances.

#	Interfering substance	Concentration
-	No addition	No addition
1	Mucin	3.5% w/v
2	Barium sulfate	5% w/v
3	Human whole blood	40% v/v
4	Metronidazole	0.25% w/v
5	Imodium	5% w/v
6	Pepto-Bismol	5% w/v
7	Vancomycin	0.25% w/v

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SYMBOL

REF	Reference number	Caution	Use by	LOT	Batch code	Consult Instructions for Use	Do not re-use
IVD	In vitro Diagnostics	Note	Manufacturer		Date of manufacture	Contains Sufficient for <= 10 tests	Keep away from sunlight
	Indicate that you should keep the product dry	To indicate the temperature limitations in which the transport package has to be kept and handled.		CE	Fulfill the requirements of Directive 98/79/EC on in vitro diagnostic medical devices	Do not use if packaging is damaged	