

LIQUIXX

CREATININE (CRE)

Jaffe's Method, Initial Rate



INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Creatinine in human serum, plasma and urine by Jaffé method.

CLINICAL SIGNIFICANCE

Creatinine is a waste product formed in muscle from the high energy storage compound, creatine phosphate. The amount of creatinine produced is fairly constant (unlike Urea) and is primarily a function of muscle mass. It is not greatly affected by diet, age, sex or exercise. Creatinine is removed from plasma by glomerular filtration and then excreted in urine without any appreciable resorption by the tubules.

Creatinine is used to assess renal function, however, serum creatinine levels do not start to rise until renal function has decreased by at least 50%.

METHODOLOGY

Modified Jaffe's reaction.

PRINCIPLE

Creatinine reacts with alkaline picrate to produce an orange-yellow colour (the Jaffe's reaction). Specificity of the assay has been improved by the introduction of an initial rate method. However , Cephalosporin antibiotics are still major interferants.

The absorbance of the orange-yellow colour formed is directly proportional to creatinine concentration and is measured photometrically at 500-520 nm.

REAGENT COMPOSITION

REAGENT 1: Picric Acid Reagent

Picric acid	25.8 mmol/L
-------------	-------------

REAGENT 2: Sodium Hydroxide Reagent

Sodium Hydroxide	95 mmol/L
------------------	-----------

Creatinine Standard

Creatinine Standard	2 mg/dl (0.166 mmol/L)
---------------------	------------------------

REAGENT PREPARATION

Reagent (R1 & R2) are ready to use. Mix equal volume of Reagent (R1 & R2) to prepare working reagent.

STORAGE AND STABILITY

Unopened Reagents 1, 2 and standard are stable till the expiry date stated on the label.

***Keep the Standard vial plugged after use, in order to avoid deterioration.**

The Working Reagent is stable for 21 days at 2-8°C. The absorbance of the reagent blank should be approximately < 0.3 at 505 nm when read against distilled water.

SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA) or urine. It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum / plasma:	7 days	at 4-25°C
	at least 3 months	at -20°C
in urine:	2 days	at 20-25°C
	6 days	at 4-8°C
	6 months	at -20°C

For the determination in urine use 24 hours specimen. It is important to exactly measure the volume of collected urine. Dilute urine samples in 1+19 ratio with distilled water and multiply results by 20. Discard contaminated specimens.

ASSAY PARAMETERS

Mode	Fixed Time e
Wavelength (nm)	505
Sample Volume (µl)	50/100
Reagent Volume (µl)	500/1000
Lag Time (Sec.)	20
Kinetic Interval (Sec.)	60
No. of readings	1
Reaction Temperature (°C)	37
Reaction Direction	Increasing
Normal Low (mg/dl)	0.6
Normal High (mg/dl)	1.4
Linearity Low (mg/dl)	0
Linearity High (mg/dl)	25
Absorbance Limit (Max.)	0.3
Blank with	DI Water
Standard Concentration (mg/dl)	2
Units	mg/dl

Programme parameters for specific clinical analysers are available on request.

ASSAY PROCEDURE

Pipette	Standard	Test
Working Reagent	1000 µl	1000 µl
Standard	100 µl	-
Test	-	100 µl

Mix well and read initial absorbance (A₁) 20 seconds after mixing and final absorbance (A₂) 80 seconds after mixing.

COMPARISON

A comparison between LIQUIXX CREATININE (y) and a commercially available test (x) using 40 samples gave following results:

$$r = 0.997$$

$$y = 0.998x + 2.564 \mu\text{mol/l}$$

CALCULATION

Calculate the results as follows:

$$\Delta A = A_2 - A_1$$

$$\text{Creatinine} = \frac{\Delta A \text{ of Test}}{\Delta A \text{ of Standard}} \times \text{Concentration of Standard (mg/dl)}$$

LINEARITY

This assay is linear upto 25 mg/dl. For higher values it is recommended to dilute the samples with normal saline and repeat the assay. Multiply the results by dilution factor.

NOTES

1. The reagent and sample volumes may be altered proportionally to accommodate various analyzer requirements.
2. The temperature of the kinetic assay must be maintained constant as the rate of colour development is strongly temperature sensitive.
3. Gross haemolysis causes falsely elevated results.
4. **Bilirubin, acetone, ascorbic acid, pyruvic acid, barbiturates, and protein interferes with Jaffe's reaction. For more comprehensive review of factors affecting creatinine assays refer to the publications of Young et al.**

QUALITY CONTROL

For quality control ERBA NORM and ERBA PATH are recommended.

EXPECTED VALUES

Male	0.7 – 1.4 mg/dl
Female	0.6 – 1.2 mg/dl
Urine	0.06 – 0.2 mg/dl

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

INTERFERENCES

Following substances do not interfere:
bilirubin up to 15 mg/dl, hemoglobin up to 10 g/l, triglycerides up to 1000 mg/dl.

Revision No.: 7- CRE

Date of Issue : 09/11/2016

ISO 9001, ISO 13485
QUALITY SYSTEM CERTIFIED

PERFORMANCE DATA

Precision :

	Within run	
	LEVEL I	LEVEL II
Number of samples (n)	20	20
Mean (mg/dl)	1.7	4.84
S. D.	0.06	0.08
C. V. %	3.53	1.65

	Between run	
	20	20
Number of samples (n)	20	20
Mean (mg/dl)	1	4.71
S. D.	0.023	0.14
C. V. %	2.3	2.97

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.
The reagent R2 contains irritating 1.0% sodium hydroxide.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.
Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

REFERENCES

1. Bowers, L.D. (1980) Clin Chem. 26; 551.
2. Bartel, H. (1972) Clin. Chem. Acta 37, 193.
3. Slot, C. (1965) Scand. J. Clin Lab. Invest 17, 381.
4. Young D.S. (1975) Clin. Chem. 21; 266D.

PACK PRESENTATION

Product Code	Pack size	Reagent-1 Picric Acid Reagent	Reagent-2 Sodium hydroxide Reagent	Creatinine Standard
120246	4 x 60ml	2 x 60ml	2 x 60ml	1 x 5ml

SYMBOLS:

The following symbols are used in the labelling of ERBA Mannheim kits:

 REF	Catalogue No	 CE	CE Mark - Device comply with the Directive 98/79/EC
 LOT	Batch Code	 IVD	In Vitro Diagnostics
	Expiry Date (Last day of the month)		Consult Instruction for Use
	Manufactured by		Storage temperature
 CRE	Product Name	 CONT	Content