



Abbott

YOUR FIRST POINT OF CARE

NYCOCARD[™] READER II

NycoCard[™] READER II is a multi-assay instrument for point-of-care testing. With NycoCard READER II, healthcare professionals in doctors' offices, clinics and hospitals will obtain accurate, quantitative results within minutes.



NycoCard™ READER II

NYCOCARD – AN ESTABLISHED TECHNOLOGY

The NycoCard™ READER II leverages robust technology with proven quality. The instrument is factory calibrated to ensure accurate results. It is a portable, point-of-care instrument designed for rapid and reliable measurements of the NycoCard test. Lot-specific calibration by the operator is not necessary. Only a simple black/white calibration at start-up and the instrument is ready for use.

BENEFITS AT A GLANCE

- Two tests on one platform: HbA1c and CRP
- Easy to use
- Rapid, quantitative results
- Ideal for measurement of small series
- Local language selection
- Easy maintenance
- Portable, battery-powered (optional)

SIMPLE TEST PROCEDURE

1 Adjusting/ black calibration



2 White calibration



3 Measure NycoCard™ test



TECHNICAL DATA

PRODUCT INFORMATION

Product Code	1116086
Weight	540 g (Including Reader Pen and Batteries)
Display Type	LCD, 2 × 16 Characters
Communication Interface	RS-232

DIMENSIONS

Instrument Box	200 mm × 170 mm × 70 mm
Reader Pen	L: 144 mm, D: 29.5 mm

POWER SUPPLY

Batteries	Rechargeable NiMH Batteries Only, Size AA 1.2 V, Recharge 8–18 hours
Battery Charger	Mascot Type 9725 or Equivalent

OPERATING CONDITIONS

Temperature	15°C–35°C (The recommended range is 18°C–25°C.)
Relative Humidity	0%–90%. Avoid direct sunlight/too dim light.

TESTS AND CONTROLS

NYCOCARD™ HbA1c
NYCOCARD™ CRP

NYCOCARD™ HbA1c CONTROL
AFINION™ CRP CONTROL