


TRUSTwell® Scrub Typhus IgM ELISA Kit



CODE#
AE0620

PACK SIZE
96 Test


SPECIMEN
Serum / Plasma

Scrub Typhus infection is caused by *Orientia tsutsugamushi*, transmitted to humans via the bites of infected chiggers. Among the lab test methods, IgM ELISA detects early infection with high sensitivity and throughput. TRUSTwell Scrub Typhus IgM ELISA utilizes a specific Recombinant Protein to improve sensitivity and specificity for better performance .

TRUSTwell® Scrub Typhus IgM ELISA Kit

- Highly specific to Gilliam, Kato & Karp stains
- Evaluated with total of 1,227 positive & negative specimens across 4 geographical zones of India
- 2 wash steps, resulting in 70 minutes
- 18 months longer shelf life

S. typhus Ab detection window



FEATURES



Highly specific to IgM detection by using targeted sensitizers



Ready to use conjugate & substrate



Increased sensitivity with high S/Co ratio