



ISO 13485:2016
ICMED 13485

HEPACARD

One Step, Rapid, Visual & Qualitative Immunoassay for HBsAg detection



Sensitivity: 0.5ng/ml


In
Collaboration with
path U.S.A.
Program for Appropriate
Technology in Health

HEPACARD


Salient Feature

- One step, rapid, visual and qualitative immunoassay
- See-through device introduced for the first time in India
- No sample preparation required
- Ready to use
- Shelf Life 30 months at 2-30°C
- 100% user safe as there is no direct contact of user with sample & strip. Thus limits the infection risk.

Simple to Perform



Test Interpretation



Performance Characteristics

- n It can detect all the **eleven** sub-types of HBsAg
- n Sensitivity: 0.5 ng/ml
- n **Evaluation Reports:**
 - **PATH (USA) Evaluation:** Sensitivity: 100%* Specificity: 100%*
 - **CMC, Vellore, India:** Sensitivity: 100%* Specificity: 100%*

*This information is provided for the Scientific Community Enquiring for an independent evaluation other than company's in house evaluation. It is not for commercial or promotional purpose.

Kit Presentation

- 20 Test Pack
- 100 Test Pack

