

WONDFO BIOTECH
WeAreWorkingForYourHealth

Wondfo

CRP

(C-reactive protein)

Finecare™

CRP RAPID QUANTITATIVE TEST


(Infection Markers)

Combination of CRP and blood routine test
increases diagnosis efficiency.

Finecare™

CRP RAPID QUANTITATIVE TEST

- ✓ A sensitive marker helps infection differentiating and antibiotic prescribing
- ✓ Results ready in 3 minutes
- ✓ Simple operation with just a drop of blood
- ✓ hsCRP and CRP dual results display
- ✓ CE Certified



CLINICAL APPLICATION

Bacterial / Viral Infection Differentiation
Help Antibiotic Prescribing
Disease Treatment Monitoring
Observation of Antibiotic Effect
Evaluation of CVD Risk




ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Loading Volume	Reaction Time	Detection Range	Storage
CRP	W201	Fingerstick Blood/ Venous Whole Blood/ Serum/Plasma	8.5 µL of whole blood or 5 µL serum or plasma	75 µL	3 minutes	0.5~200 mg/L	4~30°C for 24 months

*Wondfo CRP is the whole range CRP(hs-CRP+CRP).

SPECIAL DESIGNED CAPILLARY TUBE FOR FINGERSTIC BLOOD




ADAPTATION INSTRUMENTS

- 1. Finecare™ FIA Meter Plus (FS-113)
- 2. Finecare™ II Plus SE(FS-114)
- 3. Finecare™ III Plus (FS-205)






Point of Care
(GP, Clinics, Pharmacy)




Laboratory



In-patient
Care



Respiratory
Department



Pediatrics



Emergency
Department