

WONDFO BIOTECH
WeAreWorkingForYour **Health**



IgE

Immunoglobulin E


Finecare™

Total IgE Rapid Quantitative Test

Helping to diagnose allergic diseases

Total IgE Rapid Quantitative Test

- Whole blood samples supported
- Instant results in 10 minutes
- Room temperature storage up to 24 months




CLINICAL APPLICATION

- ✓ Helping to diagnose allergic diseases
- ✓ Predicting the development of allergic disorders
- ✓ Auxiliary diagnosis of allergic rhinitis
- ✓ Auxiliary diagnosis of asthma

ITEMS SPECIFICATION

Item	Catalog No.	Sample Type	Sample Volume	Loading Volume	Reaction Time	Detection Range	Storage
Total IgE	W298 P0001	whole blood, serum and plasma	75 µL	75 µL	10 minutes	1~1000 IU/mL	4~30°C for 24 months

TEST PROCEDURE



ADAPTATION INSTRUMENTS

- 1. Finecare™ FIA Meter Plus (FS-113)
- 2. Finecare™ II Plus SE(FS-114)
- 3. Finecare™ III Plus (FS-205)

