



QuickVET


Rapid Test range
for Companion
Animals

Diagnostics for all



QuickVET

Building on over 60 years of experience in the *In Vitro* Diagnostics industry we are aiming to improve the quality of animal care by facilitating accurate diagnosis for veterinarians worldwide. We value high quality and aim to provide innovative diagnostic solutions for all veterinarians specialising in Companion and Livestock animals.



Features and Benefits

- ✓ Rapid – One Step Immunochromatographic assays
- ✓ Accurate – High sensitivity and good correlation with other test methods
- ✓ Simple to use – One step after sample collection
- ✓ Fast – Reliable result in 10 minutes
- ✓ Immediate – Treatment and Management (vaccination/quarantine) after identifying the cause
- ✓ Economical kit size – 1T | 10T | 25T
- ✓ Shelf Life – 24 months
- ✓ Storage – +2°C - +30°C
No refrigeration required



Canine Range

Code	Product Name	Specimen	Sensitivity	Specificity
Q002-01	Canine Giardia Ag	Faeces	100% vs commercial kit	100% vs commercial kit
Q001-01	Canine E. canis Ab	WB/P/S	92% vs IFA	97% vs IFA
Q003-01	Canine Distemper Ag (CDV)	Nasal/Ocular Secretions	100% vs RT-PCR	98.7% vs RT-PCR
Q004-01	Canine Parvovirus Ag (CPV)	Faeces	100% vs HA	99% vs HA
Q005-01	Canine Brucellosis Ab	WB/P/S	93% vs RSAT	99% vs RSAT
Q007-01	Canine Rabies Ag	Saliva/Brain Homogenate	94% vs RT-PCR	100% vs RT-PCR
Q010-01	Canine Corona Virus Ag (CCV)	Faeces	94% vs PCR	96.8% vs PCR
Q014-01	Canine Adeno Virus Ag	Nasal/Ocular Secretions/Faeces	93.6% vs PCR	94.5% vs PCR
Q011-01	Canine Leishmania Ab	WB/P/S	96% vs IFA	98% vs IFA
Q025-01	Canine Heartworm Ag (CHW)	WB	93% vs Necropsy	100% vs Necropsy
Q029-01	Canine Leptospira IgM	WB/P/S	92% vs Culture	94% vs Culture
Q030-01	Canine Rota Ag	Faeces	99% vs RT-PCR	98.6% vs RT-PCR
Q038-01	CPV/CCV Combo	Faeces	CPV : 100% vs HA	CPV : 99% vs HA
			CCV : 94% vs PCR	CCV : 96.8% vs PCR
Q065-01	CPV/CDV Combo	CPV : Faeces	CPV : 100% vs HA	CPV : 99% vs HA
		CDV : Nasal/Ocular Secretions	CDV : 100% vs RT-PCR	CDV : 98.7% vs RT-PCR
Q041-01	CPV/Giardia-CCV Combo	Faeces	CPV : 100% vs HA	CPV : 99% vs HA
			Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit
			CCV : 100% vs RT-PCR	CCV : 96.8% vs PCR
Q060-01	CPV/CCV/Giardia/Rota Combo	Faeces	CPV : 100% vs HA	CPV : 99% vs HA
			CCV : 100% vs RT-PCR	CCV : 96.8% vs PCR
			Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit
			Rota : 99% vs RT-PCR	Rota : 98.6% vs RT-PCR
Q054-01	Ehrlichia/CHW/Lyme/Anaplasma Combo	Ehrlichia : WB/P/S	Ehrlichia : 92% vs IFA	Ehrlichia : 97% vs IFA
		CHW : WB	CHW : 93% vs Necropsy	CHW : 100%
		Lyme : WB/P/S	Lyme : 93% vs ELISA	Lyme : 99.3% vs ELISA
		Anaplasma : WB/P/S	Anaplasma : 97.5%	Anaplasma : 93.6%

WB: Whole blood | P: Plasma | S: Serum

Ordering Information : 1T | 10T | 20T | Uncut Sheets



Code	Product Name	Specimen	Sensitivity	Specificity
Q067-01	Rota/Adeno/ Giardia/CCV Combo	Rota : Faeces	Rota : 99% vs RT- PCR	Rota : 98.6% vs RT- PCR
		Adeno : Nasal/Ocular secretions/Faeces	Adeno : 93.6% vs PCR	Adeno : 94.5% vs PCR
		Giardia : Faeces	Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit
		CCV : Faeces	CCV : 6.61 X 10 ³ TCID 50/ml	CCV : 96.8% vs PCR
Q013-01	Canine Lyme Ab	WB/P/S	93% vs ELISA	99.3% vs ELISA
QU001-01	Canine Urine Chemistry Strip	Urine	97.20%	95.02%
Q037-02	Canine Anaplasma Ab	WB/P/S	96.2% vs ELISA	97.3% vs ELISA
Q040-01	CDV/CAV Combo	CDV : Conjunctiva	CDV :100% vs RT-PCR	CDV : 98.7% vs RT-PCR
		CAV : Nasal Fluid	CAV :93.6% vs PCR	CAV : 94.5% vs PCR
Q050-01	CPV-Giardia Combo	CPV : Faeces	CPV : 100% vs HA	CPV : 99% vs HA
		Giardia : Faeces	Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit
Q052-01	CHW-E.Canis Combo	CHW : WB	CHW : 93% vs Necropsy	CHW : 100% vs Necropsy
		Ehrlichia : WB/P/S	Ehrlichia : 92% vs IFA	Ehrlichia : 97% vs IFA
Q053-01	Canine Rota-Adeno Antigen Combo	Rota : Faeces	Rota : 99% vs RT- PCR	Rota : 98.6% vs RT- PCR
		Adeno : Nasal/Ocular secretions/Faeces	Adeno : 93.6% vs PCR	Adeno : 94.5% vs PCR
Q055-01	Canine Leishmania/ Ehrlichia Combo	Leishmania : WB/P/S	Leishmania : 92% vs ELISA	Leishmania :97% vs ELISA
		Ehrlichia : WB/P/S	Ehrlichia : 96% vs IFA	Ehrlichia : 98% vs IFA
Q056-01	Giardia/CPV/CDV Combo	Giardia : Faeces	Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit
		CPV : Faeces	CPV : 100% vs HA	CPV : 99% vs HA
		Nasal/Ocular Secretions	CDV : 100% vs RT- PCR	CDV : 98.7% vs RT- PCR
Q057-01	CPV/CCV/Giardia/CDV Combo	CPV : Faeces	CPV : 100% vs HA	CPV : 99% vs HA
		CCV : Faeces	CCV : 100% vs RT- PCR	CCV : 96.8% vs PCR
		Giardia : Faeces	Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit
		CDV : Nasal/Ocular Secretions	CDV : 100% vs RT- PCR	CDV : 98.7% vs RT- PCR
Q058-01	Giardia/CPV Combo	Faeces	Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit
			CPV : 100% vs HA	CPV : 99% vs HA
Q043-01	Canine ABST	Nasal/Ocular/Saliva Secretions or Lesions	95.80%	97.02%
Q007-03	Animal Rabies Ag	Saliva/Brain Homogenate	94% vs RT-PCR	100% vs RT-PCR
Q063-01	Canine Babesia Gibsoni Ab	WB/P/S	92%	96.8%

WB: Whole blood | P:Plasma | S:Serum

Ordering Information : 1T | 10T | 20T | Uncut Sheets



Feline Range

Code	Product Name	Specimen	Sensitivity	Specificity
Q004-02	Feline Parvovirus Ag (FPV)	Faeces	99.5% vs HA	100% vs HA
Q016-01	Feline Leukaemia Virus Ag (FeLV)	WB/P/S	94.5% vs Virus Isolation	98.9% vs Virus Isolation
Q017-01	Feline Immunodeficiency Virus Ab (FIV)	WB/P/S	97.5% vs Western Blot	98.8% vs Western Blot
Q027-01	Feline Infectious Peritonitis Ab (FIP)	WB/P/S	95% vs. ELISA	95% vs ELISA
Q012-02	Feline Toxoplasma IgG Ab	WB/P/S	98.8% vs ELISA	99% vs ELISA
Q025-02	Feline Heartworm Ag (FHW)	WB/P/S	98% vs ELISA	100% vs ELISA
Q039-01	FeLV-FIV Combo	WB/P/S	94.5% vs Virus Isolation	98.9% vs Virus Isolation
Q049-01	Feline Panleukopenia Ag	Faeces	97.05% vs PCR	98.06% vs PCR
Q010-02	Feline Coronavirus Ag	Faeces	94% vs PCR	96.8% vs PCR
Q012-01	Feline Toxoplasma-IgG/IgM	WB/P/S	98% vs ELISA	100% vs ELISA
Q012-03	Feline Toxoplasma Ab	WB/P/S	98.5% vs ELISA	99.45% vs ELISA
Q072-01	FIP-Toxoplasma IgG/IgM Combo	WB/P/S	FIP : 95% vs ELISA	FIP : 95% vs ELISA
			Toxoplasma IgG/IgM : 98.8 % vs ELISA	Toxoplasma IgG/IgM : 99.45 % vs ELISA
Q069-02	FCV-FPV Combo	Faeces	FCV : 94% vs PCR	FCV : 96.8 % vs PCR
			FPV : 99.5% vs HA	FPV : 100% vs HA
Q070-02	FIV-FeLV-FCV Combo	FIV : WB/P/S	FIV : 97.5% vs Western Blot	FIV : 98.9% vs virus isolation
		FeLV : WB/P/S	FeLV : 94.5 % vs virus isolation	FeLV : 98.9% vs virus isolation
		FCV : Faeces	FCV : 94% vs PCR	FCV : 96.8 % vs PCR
Q041-02	FPV-Giardia-FCV Combo	Faeces	FPV : 99.5% vs HA	FPV : 99.5% vs HA
			Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit
			FCV : 94% vs PCR	FCV : 96.8 % vs PCR
Q058-02	FPV-Giardia Combo	Faeces	FPV : 99.5% vs HA	FPV : 99.5% vs HA
			Giardia : 100% vs commercial kit	Giardia : 100% vs commercial kit

WB: Whole blood | P:Plasma | S:Serum

Ordering Information : IT | 10T | 20T | Uncut Sheets

Health Checks...



For Top Dogs...

Veterinary ELISA

Code	Product Name	Specimen	Sensitivity	Specificity
QE012-27	Toxo IgG-CF EIA	Serum/Plasma	98.90%	99.20%
QE077-27	HSV 1/2 IgG-CF EIA	Serum/Plasma	97.80%	97.90%
QE072-27	Chlamydia IgG-CF EIA	Serum/Plasma	98.80%	98.90%
QE073-27	Mycoplasma IgG-CF EIA	Serum/Plasma	98.99%	97.90%
QE017-22	FIV IgG EIA	Serum/Plasma	99.80%	98.09%
QE005-21	Brucella IgG-C EIA	Serum/Plasma	98.18%	99.19%
QE074-21	Candida IgG-C EIA	Serum/Plasma	97.68%	97.59%
QE004-29	Parvo IgG EIA	Serum/Plasma	98.58%	96.60%
QE003-21	CDV IgG EIA	Serum/Plasma	99.18%	96.69%
QE007-11	Rabies Virus- C EIA	Serum/Plasma	97.19%	98.02%
QE075-27	Nematoda IgG-CF EIA	Serum/Plasma	99.18%	97.90%
QE076-27	Cestoda IgG- CF EIA	Serum/Plasma	96.92%	96.92%
QE078-11	Canine CRP EIA	Serum/Plasma	97.80%	97.90%
QE101-11	cTSH EIA	Serum/Plasma	95.39%	95.62%
QE102-17	T3-CF EIA	Serum/Plasma	97.80%	96.29%
QE104-17	T4-CF EIA	Serum/Plasma	94.78%	95.22%

Ordering Information: 1 Plate | 5 Plate Kit Format: 12*8 well strip

*For research use only



Code	Product Name	Specimen	Sensitivity	Specificity
QE103-17	ft3-CF EIA	Serum/Plasma	98.82%	94.90%
QE105-17	ft4-CF EIA	Serum/Plasma	98.35%	98.62%
QE106-11	Canine LH EIA	Serum/Plasma	96.18%	94.90%
QE107-17	Progesterone- CF EIA	Serum/Plasma	97.34%	96.52%
QE108-17	Estradiol- CF EIA	Serum/Plasma	98.34%	99.09%
QE109-17	Testosterone- CF EIA	Serum/Plasma	98.22%	96.25%
QE110-7	Cortisol- CF EIA	Serum/Plasma	96.28%	95.19%
QE000-41	Canine IgE EIA	Serum/Plasma	98.18%	94.90%

*For research use only

Ordering Information: 1 Plate | 5 Plate Kit Format: 12*8 well strip



QuickVET

**Affordable
Accessible
Accurate**

