

STANDARD F**HBsAg FIA**

STANDARD F HCV Ab qualitatively analyzes antibodies specific to hepatitis C virus in serum, plasma and whole blood using Fluorescent Immunoassay (FIA). With a simple test method, the test result can be obtained within 5 minutes with its high sensitivity and specificity.

**Specification**

Specimen Type	Whole Blood, Serum, Plasma
Specimen Volume	100 µl
Testing Time	20 minutes
Storage Condition	2 - 30°C / 36 - 86°F

TEST PERFORMANCE

Reference	Sensitivity	Specificity
CLIA	99.77% (439/440)	100% (1,210/1,210)

ORDERING INFORMATION

Category	Product	Pack Size
Hepatitis	STANDARD F HCV Ab FIA	25 Test

STANDARD F**HCV ab FIA**

STANDARD F HCV Ab qualitatively analyzes antibodies specific to hepatitis C virus in serum, plasma and whole blood using Fluorescent Immunoassay (FIA). With a simple test method, the test result can be obtained within 5 minutes with its high sensitivity and specificity.

**Specification**

Specimen Type	Whole Blood, Serum, Plasma
Specimen Volume	10 µl
Testing Time	15 minutes (5 mins. early detection for strong positive samples)
Storage Condition	2 - 30°C / 36 - 86°F

TEST PERFORMANCE

Reference	Sensitivity	Specificity
CLIA	99.77% (439/440)	100% (1,210/1,210)

ORDERING INFORMATION

Category	Product	Pack Size
Hepatitis	STANDARD F HCV Ab FIA	25 Test