

STANDARD F

S.pneumoniae Ag FIA

SD BIOSENSOR

STANDARD F S.pneumoniae Ag FIA test system (Analyzer+Test device) detects S.pneumoniae antigen in urine sample for the patients with pneumonia and in cerebral spinal fluid sample for the patients with meningitis.




CE

Specification

Intended use	Qualitative detection of Streptococcus pneumoniae antigen in urine specimens.
Specimen volume	100 µl
Specimen type	Urine, CSF
Measurement range	Quantitative result with Cut off index(COI)
Testing time	10 mins (Early detection available)
Storage condition	2 – 30°C / 36 - 86°F

STANDARD TEST MODE



READ ONLY MODE



TEST PERFORMANCE

Specification	Sensitivity	Specificity
Flourescent Immunoassay	96.2% (25/26)	100% (40/40)

ORDERING INFORMATION

Category	Product	Pack Size
Respiratory	STANDARD F Legionella Ag FIA	25 Test
	STANDARD F Legionella Ag Control	Pos x 10 / Neg x 10