

STANDARD F iFOB FIA



STANDARD F iFOB FIA is a fluorescent immunoassay for the quantitative measurement of hemoglobin in fecal sample. This test is offered as a screening test for the early detection of bowel cancer in patients without symptoms.


- iFOB is a tumor marker for the diagnosis of bowel cancer
- Automatic strip recognition by reading the information stored in 2D barcode




Specification

Specimen Type	Feces
Specimen Volume	3 drops
Measurement Range	25 - 1,000 ng/ml (5 - 200µg Hb/g feces)
Reference Range	<100 bg/mL (Hb/g feces)
Testing Time	5 mins.
Storage Conditions	2-30°C / 36-86°F

SAMPLE PREPARATION



TEST PROCEDURE



ORDERING INFORMATION

Category	Product	Pack Size
Tumor Markers	STANDARD F iFOB FIA	20 Test
	STANDARD F iFOB Control	Lv1 x 10 / Lv2 x 10