

STANDARD Q

Dengue/Chikungunya Trio

 SD BIOSENSOR

STANDARD Q Arbo Panel I (Z/D/C/Y) test is an immunochromatographic assay for the detection of IgM to Zika or Dengue or Chikungunya or Yellow fever virus and Dengue NS1 antigen in human serum, plasma or whole blood. Screening test results can be obtained within 15 minutes with high sensitivity and specificity.

GENERAL INFORMATION:

- Detects Dengue IgM/IgG or Dengue NS1 antigen and IgM/IgG Chikungunya IgM/IgG specific to Dengue and Chikungunya
- Easy to use, Easy to read and interpret
- Room temperature storage



CE

Specification

Specimen	Serum/Plasma/Whole Blood
Test Time	15-20 mins. (Do not read after 20 mins.)
Specimen Volume	NS1:100µl, IgM:10µl
Storage Temperature	2-40°C / 36-104°F

ORDERING INFORMATION

Category	Product	Pack Size
Vector Borne Disease	STANDARD Q Dengue/Chikungunya Trio Test	10 Test

STANDARD Q

Zika / Dengue Fast Trio

 SD BIOSENSOR

STANDARD Q Zika / Dengue Trio test is an immunochromatographic assay for the detection of IgM and IgG antibodies specific for Dengue virus and NS1 antigen as well as Zika virus specific IgM and IgG antibodies in human serum, plasma and whole blood. Screening test results can be obtained within 15 minutes with high sensitivity and specificity.

GENERAL INFORMATION:

- Detects Zika and Dengue simultaneously at once
- Easy to use, Easy to read and interpret
- Room temperature storage



CE

Specification

Specimen	Serum/Plasma/Whole Blood
Test Time	15-20 mins. (Do not read after 20 mins.)
Specimen Volume	NS1:100µl, IgM:10µl
Storage Temperature	2-40°C / 36-104°F

ORDERING INFORMATION

Category	Product	Pack Size
Vector Borne Disease	STANDARD Q Zika /Dengue Fast Trio Test	10 Test