

STANDARD Q**Dengue IgM/IgG****SD BIOSENSOR**

STANDARD Q Dengue IgM/IgG Test is an immunochromatographic assay for the detection of IgM/IgG antibodies against Dengue virus in human serum, plasma or whole blood samples. this test kit is for in vitro use only. This is intended for professional use only for an initial screening test

GENERAL INFORMATION:



- Detection of specific IgM and IgG to Dengue virus
- Easy to use, Easy to read and interpret
- Room temperature storage

**Specification**

Specimen	Serum/Plasma/Whole Blood
Test Time	15 mins. (Do not read after 20 mins.)
Specimen Volume	10µl
Storage Temperature	2-40°C / 36-104°F

TEST PROCEDURE

Dengue IgM/IgG

**INTERPRETATION****TEST PERFORMANCE**

Description	Reference	Sensitivity	Specificity
Dengue IgM	ELISA	97.5% (77/79)	96.6% (346/358)
Dengue IgG	ELISA	97.2% (140/144)	96.2% (282/293)

ORDERING INFORMATION

Category	Product	Pack Size
Vector Borne Disease	STANDARD Q Dengue IgM/IgG	25 Test