


STANDARD Q**Dengue NS1 Ag** **SD BIOSENSOR**


STANDARD Q Dengue NS1 Ag Test is a immunochromatographic assay for the detection of Dengue NS1 antigen in human serum, plasma or whole blood samples. This test kit is for in vitro use only. This is intended for professional use, only for initial test.

GENERAL INFORMATION:


- Detects NS1 antigens
- Easy to use, Easy to read and interpret
- Room temperature storage

**Specification**

Specimen	Serum/Plasma/Whole Blood
Test Time	15 mins. (Do not read after 20 mins.)
Specimen Volume	100µl
Storage Temperature	2-40°C / 36-104°F

TEST PROCEDURE

Read at
15 mins

INTERPRETATION**TEST PERFORMANCE**

Reference	Sensitivity	Specificity
RT-PCR	92.9% (184/198)	98.7% (222/225)

ORDERING INFORMATION

Category	Product	Pack Size
Vector Borne Disease	STANDARD Q Dengue Duo Test	25 Test