

STANDARD Q HAV IgM



STANDARD Q HAV IgM Test is an immunochromatographic assay for the detection of hepatitis A virus IgM antibody in human serum, plasma or whole blood sample. This test kit is for in vitro use only. This is intended for professional use, only for an initial screening test.

GENERAL INFORMATION:

- Easy to use, Easy to read and interpret
- Room temperature storage




Specification

Specimen	Serum/Plasma/Whole Blood
Test Time	15 mins. (Do not read after 20 mins.)
Specimen Volume	10µl
Storage Temperature	2-40°C / 36-104°F

TEST PROCEDURE

HAV IgM



INTERPRETATION



TEST PERFORMANCE

Description	Reference	Sensitivity	Specificity
HAV IgM	CLIA	100% (26/26)	98.04% (450/459)

ORDERING INFORMATION

Category	Product	Pack Size
Vector Borne Disease	STANDARD Q HAV IgM Test	25 Test