

## STANDARD Q

# H. pylori Ab



STANDARD Q H. pylori Ab performs qualitative analysis by detecting antibodies specific to Helicobacter pylori in serum, plasma and whole blood using immunochromatography. With a simple test method, with its high sensitivity and specificity, the screening test result can be obtained within 10 minutes.

### GENERAL INFORMATION:

- Easy to use, Easy to read and interpret
- Room temperature storage




### Specification


Specimen	Serum/Plasma/Whole Blood
Test Time	10 mins. (Do not read after 20 mins.)
Specimen Volume	Whole Blood: 20µl, Serum/Plasma: 10µl
Storage Temperature	2-40°C / 36-104°F

### TEST PROCEDURE

#### Whole blood



#### Serum / Plasma



Read at 10 mins

### INTERPRETATION



### TEST PERFORMANCE

Description	Reference	Sensitivity	Specificity
H. pylori Ab	ELISA	98.5% (64/65)	100% (35/35)

### ORDERING INFORMATION

Category	Product	Pack Size
Gastrointestinal Disease	STANDARD Q H. pylori Ab Test	25 Test