

STANDARD Q H.pylori Ag SD BIOSENSOR

STANDARD Q H.pylori Ag Test is an immunochromatographic assay for the detection of Helicobacter pylori antigen in the human stool sample. Users can easily and rapidly get the accurate result of the test within 10 mins.


GENERAL INFORMATION:

- Easy to use, Easy to read and interpret
- Room temperature storage




Specification	
Specimen	Feces
Test Time	10 mins. (Do not read after 15 mins.)
Specimen Volume	40-70 mg
Storage Temperature	2-30°C / 36-104°F

TEST PROCEDURE



INTERPRETATION



TEST PERFORMANCE

Description	Reference	Sensitivity	Specificity
H. pylori Ag	ELISA	98.5% (64/65)	100% (35/35)

ORDERING INFORMATION

Category	Product	Pack Size
Vector Borne Disease	STANDARD Q H.pylori Ag Test	25 Test