

**STANDARD Q****Strep A Ag****SD BIOSENSOR**



STANDARD Q Strep A Ag performs qualitative analysis by detecting Group A streptococcal antigens using immunochromatography. With a simple test method, presence/absence of infection can be discriminated, and with its high sensitivity and specificity, the screening test result can be obtained within 5 minutes.

**GENERAL INFORMATION:**

- Detection of Group A Streptococcal antigens
- Easy to use, Easy to read and interpret
- Room temperature storage

**Specification**

Specimen	Throat Swab
Test Time	5 mins (Do not read after 15 minutes)
Specimen Volume	200µl
Storage Temperature	2-40°C / 36-104°F

**TEST PROCEDURE****INTERPRETATION****TEST PERFORMANCE**

Reference	Sensitivity	Specificity
Flourescent Immunoassay	98.2% (56/57)	99.26% (135/130)

**ORDERING INFORMATION**

Category	Product	Pack Size
Respiratory Disease	STANDARD Q Strep A Ag Test	25 Test
	STANDARD Q Strep A Ag Control	Pos x 10 / Neg x 10