

STANDARD Q**Zika IgM****SD BIOSENSOR**

STANDARD Zika IgM Test is a rapid chromatographic Immunoassay for the detection of Zika IgM present in serum, plasma or whole blood. This test is for in vitro professional diagnostic use and intended as an aid to early diagnosis of Zika virus infection in patient with clinical symptoms. It provides only an initial screening test result.

GENERAL INFORMATION:


- Detection of specific IgM antibody to Zika virus
- Easy to use, Easy to read and interpret
- Room temperature storage

**Specification**

Specimen	Serum/Plasma/Whole Blood
Test Time	15 mins. (Do not read after 20 mins.)
Specimen Volume	10µl
Storage Temperature	2-40°C / 36-104°F

TEST PROCEDURE

Zika IgM

**INTERPRETATION****TEST PERFORMANCE**

Reference	Sensitivity	Specificity
MAC-ELISA / PCR	92.9% (184/198)	98.7% (222/225)

ORDERING INFORMATION

Category	Product	Pack Size
Vector Borne Disease	STANDARD Q Zika IgM Test	10 Test